

## GeneProof Mycoplasma pneumoniae PCR Kit



### 1. LIST OF PRODUCT VARIANTS

Product	Package	REF
GeneProof Mycoplasma pneumoniae PCR Kit	25 reactions	MP/ISEX/025
GeneProof Mycoplasma pneumoniae PCR Kit	100 reactions	MP/ISEX/100

### 2. INTENDED PURPOSE AND USE

<b>Indication</b>	In vitro diagnostic medical device
<b>Regulatory Status</b>	CE IVD / EC Directive 98/79/EC
<b>Function</b>	Diagnostics and aid to diagnosis
<b>What is Detected / Target</b>	<i>Mycoplasma pneumoniae</i>
<b>Automated / Manual detection</b>	Manual
<b>Type of Analysis</b>	Qualitative
<b>Validated Specimen</b>	BAL, sputum, swab
<b>Intended User</b>	For professional use in laboratories with trained staff
<b>Test Principle</b>	Real-time polymerase chain reaction (PCR) – amplification of the specific Target Sequence and detection using TaqMan probes with fluorophore-based detection

### 3. TECHNICAL SPECIFICATION

<b>Target Sequence</b>	M181 gene encoding the CARDS toxin		
<b>Specificity</b>	<i>Mycoplasma pneumoniae</i> (MP), 100 %		
<b>Sensitivity (LoD with 95 % probability)</b>	0.46 cp/µl		
<b>Extraction/Inhibition Control</b>	PCR inhibition and DNA extraction efficiency control by Internal Control (IC)		
<b>Validated Extraction Methods</b>	croBEE 201A Nucleic Acid Extraction Kit GeneProof PathogenFree DNA Isolation Kit		
<b>Applied Instruments</b>	<b>Instrument Name</b>	<b>MP</b>	<b>Internal Control (IC)</b>
	croBEE Real-Time PCR System	FAM	HEX
	AMPLilab Real-Time PCR System	FAM	HEX
	Applied Biosystems 7300 / 7500 Real-Time PCR System	FAM	JOE
	AriaMx Real-Time PCR System	FAM	HEX
	BioQuant-96 Real-Time PCR System	FAM	HEX
	CFX Connect™ / CFX96™/ Dx Real-Time PCR Detection System	FAM	HEX
	Gentier 96E/96R Real-Time PCR System	FAM	HEX
	LightCycler® 2.0 / 480	FAM	HEX
	LineGene 9600 / 9600 Plus	FAM	HEX
	Mic qPCR Cycler	FAM	HEX
	QuantStudio™ 3 / 5 Real-Time PCR System	FAM	VIC
	Rotor-Gene 3000 / 6000 / Q	FAM	JOE
SLAN® Real-Time PCR System	FAM	HEX	
StepOne™/StepOne Plus™ Real-Time PCR System	FAM	JOE	
<b>Detection Channels</b>	FAM (MP), HEX/JOE/VIC (IC)		
<b>External Quality Assessment</b>	Regularly tested in QCMD and INSTAND e.V. External Quality Assessment Panels - results at <a href="http://www.geneproof.com">www.geneproof.com</a>		

### 4. KIT CONTENT

Reagent	Content	Vial Title	Cap Colour	Guaranteed Volume [µl]	Number of Vials	
					MP/ISEX/025 – 25 rxn	MP/ISEX/100 – 100 rxn
Master Mix	Mixture of PCR enzymes, target specific primers and TaqMan probes in buffer	MasterMix <i>Mycoplasma pneumoniae</i>	Blue	750	1	4
Positive Control	DNA oligonucleotide in buffer	Positive Control <i>Mycoplasma pneumoniae</i> 10 <sup>2</sup> cp/µl	White	200	1	2
Internal Control	Plasmid DNA in buffer	Internal Control <i>Bordetella</i> <i>Chlamydia pneumoniae</i> <i>Mycobacterium tuberculosis</i> <i>Mycoplasma pneumoniae</i>	Red	1000	1	2

## 5. CALIBRATOR INFORMATION

No calibrators - qualitative detection only.

## 6. TRANSPORT AND STORAGE

<b>Storage Conditions</b>	(-20 ± 5) °C
<b>Transport Conditions</b>	-20 °C and below
<b>In-use Stability</b>	Thaw a maximum of 5 times or use within 30 days after the first use of a particular vial, whichever comes first

## 7. ASSAY PROCEDURE

### SPECIMEN COLLECTION, TRANSPORTATION AND HANDLING

1. Samples for DNA extraction must be collected and transported following professional guidelines, preferably at the temperature (5 ± 3) °C.
2. Samples for DNA extraction must be transported and treated by the laboratory in the shortest possible time (preferably within 24 hour).

### NUCLEIC ACID PURIFICATION

3. Prepare specimens for the assay according to the corresponding extraction kit manual.
4. Thaw required amount of Internal Control (IC) vials, mix and briefly centrifuge.
5. Add Internal Control (IC) directly into the sample at the beginning of the extraction process so that 1 µl of the resulting elution volume contains 0.1 µl of the IC:

<b>Elution Volume</b>	25 µl	50 µl	100 µl	200 µl
<b>Internal Control (IC)</b>	2.5 µl	5 µl	10 µl	20 µl

6. Continue extraction according to the appropriate protocol.

### PCR SETUP PROTOCOL

7. Thaw required vials and reagents completely.
8. Gently vortex and briefly centrifuge all vials before setting up the PCR run.

**NOTE:**

Keep the reagents at (5 ± 3) °C for the shortest time possible until the PCR reaction is set up.

9. Add 30 µl of Master Mix into PCR tubes.
10. Add 10 µl of the extracted nucleic acid sample or 10 µl of Positive Control into the individual PCR tubes and mix by pipetting. The total reaction mix volume is 40 µl.
11. Close the tubes, centrifuge shortly, insert them into the real-time PCR device and amplify according to the following PCR profile.

**NOTE:** It is recommended to perform at least 1 negative control and at least 1 positive control (for a qualitative kit) for each individual PCR run. Use your own negative control (not provided) in the form of nuclease-free water. For more information see chapter 9. Run Validity.

### AMPLIFICATION PROFILE:

Follow the thermocycler manufacturer's manual for setting the instrument and for analysis.

#### DNA PCR Profile

Step	Process	Temperature [°C]	Time	Cycles	Fluorescence Acquisition
1	UNG decontamination	37	2 min	1 cycle	
2	Initial denaturation	95	10 min	1 cycle	
3	Denaturation	95	5 s	45 cycles	FAM, HEX/JOE/VIC
	Annealing	60	40 s		
	Extension	72	20 s		

## 8. INTERPRETATION OF RESULTS

Channel FAM (MP)	Channel HEX/JOE/VIC (IC)	Result	Interpretation
+	+	Valid	<i>M. pneumoniae</i> positive
+	-	Valid	<i>M. pneumoniae</i> positive
-	+ (Ct<38)	Valid	<i>M. pneumoniae</i> negative
-	+ (Ct>38]	Invalid	-
-	-	Invalid	-

**NOTE:** For interpretation of PCR run see chapter 9. Run Validity.

## 9. RUN VALIDITY

### OVERALL VALIDITY OF DETECTION

	Signal	Channel	Run Validity	Recommendation
Positive Control	+	FAM	Valid	-
Positive Control	-	FAM	Invalid	Repeat PCR run
Negative control	-	FAM	Valid	-
Negative control	+	FAM	Invalid	Repeat PCR run

**NOTE:** If the issue persists, please contact Customer Support, see Contact information.

## 10. QUANTITATIVE DETECTION EVALUATION

Qualitative detection only.

## 11. ADDITIONAL PRODUCTS

No additional products.

## 12. MATERIAL AND DEVICES REQUIRED BUT NOT PROVIDED

### CONSUMABLE MATERIAL

96-well reaction plates or PCR tubes (0.2 ml volume), pipette tips with filters, powder-free gloves, biohazard waste bin, nuclease-free water.

### DEVICES

Real-time PCR instrument (see chapter 3. Technical Specification), nucleic acid extraction system or kit (see chapter 3. Technical Specification), desktop centrifuge (for 0.2 ml PCR tubes or 96-well plates), vortex mixer, freezer (-20 ± 5) °C, refrigerator (5 ± 3) °C, cooling rack

## 13. WARNINGS, PRECAUTIONS AND PROCEDURE LIMITATIONS

- Patient management decisions should never be made based solely on the results from this test. Other laboratory and clinical factors must also be considered in making clinical decisions.
- Any serious incident occurred in relation to the using of GeneProof PCR Kit shall be reported to the manufacturer and to the competent local authority.
- Use all necessary protective equipment (protective disposable gloves, a laboratory coat and eye protection) when handling specimens and kit reagents.
- Avoid microbial and ribonuclease contamination of the reagents when removing aliquots from reagent vials.
- Use RNase- and DNase-free filter pipette tips only.
- Use new tips for each pipetting step.
- Use separate working places for sample preparation / nucleic acid extraction and amplification reactions. Never introduce an amplified product in reagent and/or nucleic acid extraction (sample preparation) area.
- Close the kit components vials immediately after use and never interchange lids.
- Do not pool reagents from different lots or from different vials within the same lot.
- Do not substitute the reagents between different lots.
- Do not use reagents from damaged or leaking vials.
- Dispose of unused reagents and waste in accordance with national, federal, state or local regulations.
- Be very careful when handling the Positive Control or the clinical material; incorrect handling could result in contamination and the consequent impairment of the kit components! The manufacturer is not responsible for the kit impairment due to incorrect handling.










### Procedure Limitations:

- Read the whole Instruction for Use properly before starting the manipulation. Not following these instructions can lead to an erroneous result which can cause misdiagnosis or inappropriate treatment!
- Appropriate specimen collection, transport, storage and processing procedures are required for the optimal performance of this test.
- Do not use kit after the expiry date.
- The presence of UDG decontamination step reduces the risk of amplicon contamination. However, contamination from MP-positive controls and clinical specimens can be avoided only by good laboratory practices and careful adherence to the procedures specified in this Instruction for use. All reagents should be closely monitored for impurity and contamination. Any suspicious reagents should be discarded.
- This product is designed for use with the applied PCR instruments and validated extraction methods mentioned in chapter 3. Technical Specification

### Clinical Limitations:

- Use only with validated specimens (see chapter 2. Intended Purpose and Use) otherwise incorrect results could occur.

## 14. EXPLANATION OF SYMBOLS

Symbol	Explanation	Symbol	Explanation
	This product complies with the relevant EU requirements		Lot number
	<i>in vitro</i> diagnostic medical device		Contains sufficient amount for n-tests
	Catalogue number		Temperature limitation
	Manufacturer		Expiry date
	Read Instruction for Use		

### Customer Support

Tel.: +420 730 176 222  
E-mail: support@geneproof.com

### Orders

Tel.: +420 543 211 679  
E-mail: sales@geneproof.com



GeneProof a.s.  
Václavská 101/119 /Dolní Heršpice/ CZ-619 00 Brno / +420 543 211 679  
info@geneproof.com

