HLA-B5701 RealFast™ Assay

ViennaLab Diagnostics GmbH
Gaudenzdorfer Guertel 43-45
A-1120 Vienna, Austria
Phone: (+43-1) 8120156-0
info@viennalab.com
www.viennalab.com

REF 7-610 / 7-613 Σ 100 / 32 reactions
-30°C / -15°C CE IVD

1. Intended Use

The HLA-B5701 RealFast™ Assay is a fast and accurate real-time PCR test for detection of the HLA-B5701 allele, a specific variant of the *human leukocyte antigen B (HLA-B)* gene that is strongly associated with abacavir hypersensitivity. The kit is designed for genetic risk stratification of HIV-infected patients prior to initiation of abacavir therapy. HLA-B5701 positive patients must be excluded from abacavir treatment. The qualitative assay discriminates the presence or absence of HLA-B5701 in a human genomic DNA extract. The kit does not interfere with the closely related HLA-B5702 and HLA-B5703 variants.

Reference sequence: HGVS: NG_023187.1

2. Introduction

Abacavir is a nucleoside analogue reverse-transcriptase inhibitor used to treat HIV-1 infection. Approximately 5% to 8% of Caucasians with HIV-1 infection who undergo a combination antiretroviral therapy containing abacavir develop a hypersensitivity reaction. Manifestations usually appear within the first 6 weeks after start of the treatment and include multisystem involvement resulting in skin rash, fever, constitutional, gastrointestinal tract or respiratory symptoms that become more severe with continued dosing. Immediate and permanent discontinuation of treatment leads to a rapid reversal of side effects. Conversely, re-challenge with abacavir after a hypersensitivity reaction can provoke potentially life-threatening conditions. Prospective employment of routine screening for HLA-B5701 before starting therapy significantly reduces the incidence of hypersensitivity reactions due to abacavir.

3. Kit Contents

100 / 32 Rxn

| | | | |
|-----------------------------|--------|--------------|---------------|
| RealFast™ 2x Genotyping Mix | 1 vial | □ white cap | 1000 / 320 µl |
| HLA-B5701 Assay Mix | 1 vial | ■ purple cap | 550 / 550 µl |
| HLA-B5701 Positive Control | 1 vial | ■ green cap | 75 / 75 µl |
| HLA-B5701 Negative Control | 1 vial | ■ red cap | 75 / 75 µl |

The RealFast™ 2x Genotyping Mix comprises HotStart Taq DNA polymerase and dNTPs in an optimized buffer system.

The HLA-B5701 Assay Mix consists of gene-specific primers and dual-labeled hydrolysis probes for *HLA-B5701* and a control gene.

The kit contains reagents for 100 / 32 reactions in a final volume of 20 µl each.

A positive and a negative control for HLA-B5701 are supplied with the kit.

4. Storage and Stability

HLA-B5701 RealFast™ Assay is shipped on cooling blocks. On arrival, store the kit at -30 to -15°C. Alternatively, store at 2 to 8°C for short-term use within one month. The kit withstands up to 20 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to intense light. If stored correctly, the kit will retain full activity until the expiration date indicated on the label.

5. Product Description

5.1. Principle of the Test

The test is based on the fluorogenic 5' nuclease assay, also known as TaqMan® assay. Each reaction contains gene-specific primer pairs which amplify a 97 bp fragment of the *HLA-B5701* gene and a 147 bp fragment of a control gene, the latter serving as PCR control. Further components are two dual-labeled, gene-specific hydrolysis probes which hybridize to the target sequence of the corresponding fragment. The proximity of the 5'-fluorescent reporter and 3'-quencher dye on intact probes prevents the reporter from fluorescing. During the extension phase of PCR the 5' – 3' exonuclease activity of the Taq DNA polymerase cleaves the 5'-fluorescent reporter from the hybridized probe. The physical separation of the fluorophore from the quencher dye generates a fluorescent signal in real-time, which is proportional to the accumulated PCR product.

In samples positive for HLA-B5701 both, the **FAM-labeled HLA-B5701** probe as well as the **HEX-labeled PCR control** probe bind to the appropriate gene fragment. A strong fluorescence signal is detected in the FAM channel (520 nm) and in the HEX channel (556 nm). In samples negative for HLA-B5701 only the HEX-labeled PCR control probe hybridizes to the complementary strand of the control gene fragment. A strong fluorescence signal is detected in the HEX channel and no or only a baseline signal in the FAM channel.

5.2. Real-time PCR Instrument Compatibility

The HLA-B5701 RealFast™ Assay is validated for use with the AB 7500 Fast instrument.

The kit is compatible with various common real-time PCR instruments capable of recording FAM and HEX fluorescence:

- ✓ AB 7500 Fast (Applied Biosystems®)
- ✓ AB StepOne™ (Applied Biosystems®)
- ✓ CFX96™ (Bio-Rad)
- ✓ LightCycle® 480 (Roche)
- ✓ MIC qPCR Cycloer (bms)
- ✓ Mx3005P (Agilent Technologies)
- ✓ Rotor-Gene® 6000 (Qiagen)

» **Note:** RealFast™ Variant Detection QuickGuides for setting up and analyzing experiments on different types of instruments can be downloaded from www.viennalab.com.

When using AB StepOne™, set passive reference dye to "ROX"! «

The kit is supplied **without ROX**. For use with real-time PCR instruments requiring high ROX for normalization of data (e.g. Applied Biosystems® instruments: StepOne™, 7300, 7900/7900HT), add ROX at a final concentration of 1 µM to the 2x Genotyping Mix.

5.3. Assay Performance Specifications

Determination of **sensitivity** was performed on 64 alleles testing positive for the HLA-B5701 allele with a reverse SSO based reference kit. The HLA-B5701 RealFast™ Assay determined all 64 alleles as positive, which equaled a true positive rate of 100%.

Determination of **specificity** was performed on 66 alleles testing negative for the HLA-B5701 allele with a reverse SSO based reference kit. The HLA-B5701 RealFast™ Assay determined all 66 alleles as negative, which equaled a true negative rate of 100%.

Limit of detection: 0.2 ng genomic DNA (per reaction). Recommended DNA concentration: 2 to 20 ng/µl genomic DNA.

Detectable HLA-B5701 alleles are listed in the Productnote 03: HLA-B5701 (www.viennalab.com/products/pharmacogenetics/hla-b5701).

6. Materials Required but not Supplied

Real-time PCR instrument with FAM (520 nm) and HEX (556 nm) filters, instrument-compatible reaction vessels, disposable powder-free gloves, vortexer, mini-centrifuge for 2.0 ml tubes, tube racks, set of calibrated micropipettes (0.5 – 1000 µl), sterile tips with aerosol-barrier filter, molecular grade water, DNA extraction system, freezer, biohazard waste container.

7. Experimental Protocol

7.1. DNA Extraction

DNA extraction reagents are **not supplied** with the kit.

DNA isolated from various specimens (e.g. whole peripheral blood, dried blood spots, buccal swabs or saliva) can be used. Ensure extracted DNA is suitable for amplification in terms of concentration, purity and integrity.

For accurate genotype calling, the DNA amount per reaction should be within the range of 10 to 100 ng for all samples.

7.2. PCR Controls

Always include a **No Template Control (NTC)** in each experiment to confirm absence of potential contamination. It is advisable to run the NTC (use PCR-grade water instead of DNA) in duplicate.

Always include the HLA-B5701 **Positive Control** as positive reference signal.

Always include the HLA-B5701 **Negative Control** as negative reference signal for threshold setting in the FAM channel.

» **Note:** The Controls are potential sources of contamination. Make sure to handle them carefully. «

7.3. Preparation of HLA-B5701 RealFast™ Master Mix

Gently vortex and briefly centrifuge all solutions after thawing. Set up PCR at room temperature. Prepare sufficient **Master Mix** for all your reactions (N samples + positive control + negative controls) plus at least one additional reaction to compensate for pipetting inaccuracies:

| Component | per reaction | e.g. 24+1 reactions |
|-----------------------------|--------------|---------------------|
| RealFast™ 2x Genotyping Mix | 10 µl | 250 µl |
| HLA-B5701 Assay Mix | 5 µl | 125 µl |
| Master Mix | 15 µl | 375 µl |

Dispense **15 µl Master Mix** into each well. Add **5 µl purified DNA or Control** template to reach a final reaction volume of 20 µl.

To minimize risk of contamination, always pipette templates in the following order: first NTC, then samples, last controls. Immediately close reaction vessels.

» **Note:** Avoid creating bubbles in the final reaction mix and avoid touching the optical surface of the cap or sealing film without gloves. Both may interfere with fluorescence measurements. Centrifuge briefly if needed. «

7.4. PCR Program

Program the real-time PCR instrument according to the manufacturer's instructions for quantitation experiments with two targets / reporter dyes. Place the samples into the thermal cycler and run the following program:

AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480, Mx3005P and other Peltier heating block-based instruments:

| Cycles | Temp | Time | Steps |
|--------|------|--------|---|
| 1 | 95°C | 3 min | Initial denaturation |
| 40 | 95°C | 15 sec | Denaturation |
| | 60°C | 1 min | Annealing/Extension – Data acquisition on FAM and HEX channel |

MIC qPCR Cycler, Rotor-Gene® 6000*):

| Cycles | Temp | Time | Steps |
|--------|-------------------------------------|--------|---|
| 1 | 95°C | 3 min | Initial denaturation |
| 40 | 95°C | 15 sec | Denaturation |
| | 60°C *for 36-well rotor: 56°C | 1 min | Annealing/Extension – Data acquisition on Green and Yellow channel |

8. Data Analysis / Interpretation of Results

The presence or absence of the HLA-B5701 allele is defined by whether there is a signal in the **FAM channel** or not. Successful PCR can be verified by amplification of the control gene detected in the **HEX channel** (PCR control). Thus, genomic DNA samples positive for HLA-B5701 as well as the HLA-B5701 Positive Control show amplification in both, the HEX and the FAM channel. HLA-B5701 negative samples and the HLA-B5701 Negative Control show amplification in the HEX channel only. Fluorescent levels and corresponding amplification curves are automatically displayed in amplification plots in the real-time PCR software.

| Sample Type | Amplification in FAM channel (520 nm) | Amplification in HEX channel (556 nm) |
|----------------------------|---------------------------------------|---------------------------------------|
| HLA-B5701 positive | YES | YES |
| HLA-B5701 negative | NO | YES |
| HLA-B5701 Positive Control | YES | YES |
| HLA-B5701 Negative Control | NO | YES |
| NTC | NO | NO |

Some instrument software may need manual threshold settings for accurate analysis.

Recommendations for Threshold Settings (C_q):


Set threshold value for the FAM channel just above the background fluorescent signal generated by the HLA-B5701 Negative Control.

Samples crossing the threshold line beyond C_q 37 give invalid results and must be repeated.

To analyze acquired data, please follow your instrument software instructions.

9. Warnings and Precautions

- For *in vitro* diagnostics use only.
- Always use disposable powder-free gloves and wear suitable lab coat when handling specimens and reagents.
- Perform reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- Use pipettes dedicated for PCR setup only, use aerosol-guarded pipette tips.
- Use instrument-compatible reaction vessels with optically clear caps or sealers.
- Do not mix reagents from different lots.
- Do not use expired kits or kit components.

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AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480, Mx3005P e outros instrumentos de bloco de calor Peltier: