

Protocol 1, part A: Fluid Sample Preparation

Kit 1: Buffers & Consumables (+18 to +25°C)

Kit 3: Consumables (+18 to +25°C)

- Pipette 10ml of the fluid specimen in the *50ml Centrifuge tube* (Kit 3).
- In case of less sample volume available, pipette the fluid into the *50ml Centrifuge tube* and fill up to 10ml with buffer *SU* (Kit 1) (use the measure line of the tube).

Continue with the instructions of the scheme Ultra-Deep Microbiome Prep10 / Protocol 1, part B: DNA Isolation (page 3, short manual).



***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.

Protocol 2, part A: Tissue Sample Preparation

Kit 1: Buffers & Consumables (+18 to +25°C)

Kit 2: Enzymes & Reagents (-15 to -25°C)

- Transport the specimen under conditions avoiding contamination to the laboratory for analysis.
- Per specimen, place a *Sample tube* (ST, Kit 1,) in a rack and mark. Pipette 180µl of buffer *PKB* (Kit 1) into the *Sample tube* (ST).
- Transfer the specimen to a sterile support (e.g., Petri dish). Cut the specimen (~0.5x0.5cm) into small pieces by using a sterile scalpel.
- Transfer the cut specimen to the *Sample tube* (ST) filled with buffer *PKB*.
Add 20µl of *Proteinase K* (Kit 2), vortex for 15s.
Incubate at 56°C, 10min, 1,000rpm (thermomixer).
- Fill up to 1ml with the transport solution, if available, or with buffer *TSB* (use the measure line of the tube).



Continue with the instructions of the scheme Ultra-Deep Microbiome Prep10 / Protocol 2, part B: DNA Isolation (page 5, short manual).

***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.

Kit 1: Buffers & Consumables (+18 to +25°C)

Kit 2: Enzymes & Reagents (-15 to -25°C)

Arrange bottles according to the sequence of steps as below:

CM – DB1 – RS – RL – RP – CS – AB – WB – 70% Ethanol – Deionized Water

Continued from Ultra-Deep Microbiome Prep10 / Fluid Sample Preparation,
part A (page 1, short manual).

Per sample:

1. Add **4ml buffer CM**, vortex for 15s.
Let stand at room temperature (+18 to +25°C) for 5min.
2. Add **4ml buffer DB1**.
Add **10µl MolDNase B** (Kit 2), vortex for 15s.
Incubate at room temperature (+18 to +25°C) for 15min.
3. Centrifuge at 9,500xg, 10min.
Remove supernatant by pipetting and discard.
4. Resuspend pellet in **1ml buffer RS** by pipetting.
5. Transfer the suspension by pipetting into a *Sample tube* (ST, Kit 1) and centrifuge at $\geq 12,000xg$, 5min.
Remove supernatant by pipetting.
(Optional: freeze pellet at -15 to -25°C for storage).
6. Resuspend pellet in **80µl buffer RL**, briefly centrifuge tube.
7. Add **20µl BugLysis** (Kit 2).
Add **1.4µl β -mercaptoethanol** (Kit 2), vortex for 15s.
Take care not to inhale.
Incubate at 37°C, 30min, 1,000rpm (thermomixer).
8. Briefly centrifuge.
Add **150µl buffer RP**.
Add **20µl Proteinase K** (Kit 2), vortex for 15s.
Incubate at 56°C, 10 min, 1,000 rpm (thermomixer).

Continue on page 4

Removal of Human DNA

Lysis of Pathogens



***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.

During 10 min incubation:

Kit 1: Buffers & Consumables

Unpack *Spin columns (SC)*, 2 ml *Collection tubes (CT)* and 1.5 ml *Elution tubes (ET)*, label; heat **Deionized Water** (100µl each sample) vial to **70°C** (thermomixer).

9. Briefly centrifuge.

Add **250µl buffer CS**, vortex for 15s.

10. Briefly centrifuge.

Add **250µl buffer AB**, vortex for 15s.

11. Briefly centrifuge to clear lid.

Pipette lysate into a *Spin column*.

Centrifuge: $\geq 12,000 \times g$, 30 to 60s.

12. Remove column and place in a new 2 ml *Collection tube*.

Add **400µl buffer WB**.

Centrifuge: $\geq 12,000 \times g$, 30 to 60s.

13. Remove column and place in a new 2 ml *Collection tube*.

Add **400µl 70% Ethanol**.

Centrifuge: $\geq 12,000 \times g$, 3min.

14. Carefully remove column and place in a 1.5 ml *Elution tube*.

15. Add **100µl Deionized Water** heated to 70°C.

Incubate at room temperature (+18 to +25°C) for 1min.

Centrifuge: $\geq 12,000 \times g$, 1min.

Discard column, close lid of *Elution tube*.

16. Store eluted DNA (1.5 ml *Elution tube*) at -15 to -25°C.

Lysis of
Pathogens

DNA Purification

DNA Elution

***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.

Kit 1: Buffers & Consumables (+18 to +25°C)

Kit 2: Enzymes & Reagents (-15 to -25°C)

Arrange bottles according to the sequence of steps as below:

CM – DB1 – RS – RL – RP – CS – AB – WB – 70% Ethanol – Deionized Water

Continued from Ultra-Deep Microbiome Prep10 / Tissue Sample Preparation, part A (page 2, short manual).

Per sample:

1. Add **250µl buffer CM**, vortex for 15s.
Let stand at room temperature (+18 to +25°C) for 5min.
2. Briefly centrifuge.
Add **250µl buffer DB1**.
Add **10µl MolDNase B** (Kit 2), vortex for 15s.
Incubate at room temperature (+18 to +25°C) for 15min.
3. Centrifuge at $\geq 12,000 \times g$, 10min.
Remove supernatant by pipetting and discard.
4. Resuspend pellet in **1ml buffer RS**.
5. Centrifuge at $\geq 12,000 \times g$, 5min.
Remove supernatant by pipetting.
(Optional: freeze pellet at -15 to -25°C for storage).
6. Resuspend pellet in **80µl buffer RL**, briefly centrifuge tube.
Add **20µl BugLysis** (Kit 2).
Add **1.4µl β -mercaptoethanol** (Kit 2), vortex for 15s.
Take care not to inhale.
Incubate at 37°C, 30min, 1,000rpm (thermomixer).
7. Briefly centrifuge.
Add **150µl buffer RP**.
Add **20µl Proteinase K** (Kit 2), vortex for 15s.
Incubate at 56°C, 10 min, 1,000 rpm (thermomixer).

Removal of Human DNA

Lysis of Pathogens

Continue on page 6

***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.

During 10 min incubation:

Kit 1: Buffers & Consumables

Unpack *Spin columns (SC)*, 2 ml *Collection tubes (CT)* and 1.5 ml *Elution tubes (ET)*, label; heat **Deionized Water** (100µl each sample) vial to **70°C** (thermomixer).

8. Briefly centrifuge.

Add **250µl buffer CS**, vortex for 15s.

9. Briefly centrifuge.

Add **250µl buffer AB**, vortex for 15s.

10. Briefly centrifuge to clear lid.

Pipette lysate into a *Spin column*.

Pipette the fluid phase in the column.

Avoid transfer of any unresolved particles!

Centrifuge: $\geq 12,000 \times g$, 30 to 60s.

11. Remove column and place in a new 2 ml *Collection tube*.

Add **400µl buffer WB**.

Centrifuge: $\geq 12,000 \times g$, 30 to 60s.

12. Remove column and place in a new 2 ml *Collection tube*.

Add **400µl 70% Ethanol**.

Centrifuge: $\geq 12,000 \times g$, 3min.

13. Carefully remove column and place in a 1.5 ml *Elution tube*.

14. Add **100µl Deionized Water** heated to 70°C.

Incubate at room temperature (+18 to +25°C) for 1min .

Centrifuge: $\geq 12,000 \times g$, 1min.

Discard column, close lid of *Elution tube*.

15. Store eluted DNA (1.5 ml *Elution tube*) at -15 to -25°C.

Lysis of
Pathogens

DNA Purification

DNA Elution

***Please note: Before using this scheme inform yourself of the details of the procedure. Please consult the manual.** Special care is required for working under DNA-free conditions and secure working conditions, please consult the manual for more information.