


GeneProof HIV type 1 (HIV-1) PCR Kit



In vitro diagnostic medical device

The kit has been manufactured according to the EC Directive 98/79/EC as an *in vitro* diagnostic medical device and it has been designed for professional use in specialized clinical and research laboratories.

KIT CONTENT

REF:	HIV1/ISEX/025 25 rxn	HIV1/ISEX/100 100 rxn
MasterMix HIV-1	1x 750 µl	4x 750 µl
Calibrator A HIV-1 10 ⁵ IU/µl	1x 400 µl	1x 400 µl
Calibrator B HIV-1 10 ⁴ IU/µl	1x 400 µl	1x 400 µl
Calibrator C HIV-1 10 ³ IU/µl	1x 400 µl	1x 400 µl
Calibrator D HIV-1 10 ² IU/µl	1x 400 µl	1x 400 µl

STORAGE AND TRANSPORTATION CONDITIONS

The kit must be transported at temperature below -20 °C. The kit will remain stable at least until the expiry date printed on the package, if the storage temperature is kept (-20 ± 5) °C. The components are stable for a maximum of 4 repeated freezing/thawing cycles after the first use of a particular vial. The component must be used before the expiry date or 25 days after the first use of a particular vial (whichever comes first).

TECHNICAL SPECIFICATION

Target Sequence	LTR sequence and GaG gene
Analytical Specificity	HIV genotypes A - D, AE, F, AG-GH, Group N, Group O, BF, H, K, CRF03_AB, 100 %
Analytical Sensitivity	273.971 IU/ml i.e. 153.424 cp/ml (on HIV-1 NIBSC 16/194 using GeneProof PathogenFree RNA Isolation Kit)
(LoD with the probability of 95%)	548.121 IU/ml i.e. 306.948 cp/ml (on HIV-1 NIBSC 16/194 using croBEE 201A Nucleic Acid Extraction Kit) 98.59 IU/ml i.e. 55.21 cp/ml (on Acrometrix HIV-1 Panel UI/ml using SpinStar Viral Nucleic Acid Kit 1.0 with SpinStar Pretreatment Solution)
Diagnostic Specificity	100 % (CI _{95%} : 99.10 % - 100 %)
Diagnostic Sensitivity	93.66 % (CI _{95%} : 87.96 % - 96.88 %)
Linear Range	10 ⁹ - 10 ^{2.5} IU/ml with precision of ± 0.5 log (GeneProof PathogenFree RNA Isolation Kit) 10 ⁹ - 10 ³ IU/ml with precision of ± 0.5 log (croBEE 201A Nucleic Acid Extraction Kit)
Dynamic Range	10 ⁹ - 273.971 IU/ml (GeneProof PathogenFree RNA Isolation Kit) 10 ⁹ - 548.121 IU/ml (croBEE 201A Nucleic Acid Extraction Kit)
Reporting units	IU/ml (1 IU = 0.56 cp)
Validated Specimen	plasma
Quality Control	regularly tested by QCMD and Instand e.V. External Quality Assessment Panels
Regulatory Status	CE ₁₀₂₃ IVD

Quality management system is certified in compliance with the requirements of the standard ČSN EN ISO 13485 ed.2:2016

INTERFERENCES

The interferences testing was performed using negative plasma with set level of biochemical markers, which can be potential endogenous interferences. The negative plasma was spiked with HIV positive material at 3x LoD. Elevated levels of bilirubin (342 µmol/l), albumin (60 g/l), haemoglobin (2 g/l), urea (42.9 mmol/l), uric acid (1.4 mmol/l) and D-glucose (55 mmol/l) in samples have been tested in the presence and absence of HIV DNA. The evaluation and settings of pathological values for interference testing were performed according to CLSI guidelines EP7-A2.

PLASMA

Tested substance	Tested level(s)	Observed interference	Tested substance	Tested level(s)	Observed interference
Albumin	60 g/l	None	Haemoglobin	2 g/l	None
Bilirubin	342 µmol/l	None	Urea	42.9 mmol/l	None
Glucose	55 mmol/l	None	Uric acid	1.4 mmol/l	None

The tested endogenous interferences were shown not to interfere with the GeneProof HIV type 1 (HIV-1) PCR Kit in plasma samples.

METHOD PRINCIPLE

The HIV-1 detection is based on the amplification of a specific region of the LTR sequence and a region of the *GaG* gene sequence and the measurement of fluorescence increase. The kit provides maximum sensitivity and specificity for all variants of the HIV-1 virus from the M group (including groups N and O) and for the virus CRF variants. The presence of HIV-1 is indicated by the increased FAM fluorophore fluorescence. The human GAPDH mRNA is used as an Internal control (IC), it serves as a control for the whole diagnostic process, i.e. RNA extraction efficiency, reverse-transcription step efficiency (transcription of RNA into cDNA) and PCR amplification efficiency (PCR inhibition). The IC positive amplification is detected in the HEX fluorophore fluorescence channel. The PCR kit is designed as Ready-to-Use Master Mix for *in vitro* diagnostics for both qualitative and quantitative detection and it utilizes the "hot start" technology minimizing non-specific reactions and ensuring maximum sensitivity.

USER MANUAL

SAMPLING AND SAMPLE STORAGE

Whole blood samples should be centrifuged and the plasma removed to secondary tube within 4 hours from blood collection. If there is not possibility to process the sample within 4 hours, it is recommended to centrifuge whole blood and freeze the plasma sample at -20 °C or -70 °C.

NUCLEIC ACID PURIFICATION

Nucleic acid extraction should be performed by extraction kits available at the market according to protocols for the particular clinical material extraction. The manufacturer recommends the following products:

- GeneProof PathogenFree RNA Isolation Kit
- croBEE 201A Nucleic Acid Extraction Kit

PCR SETUP

1. Gently vortex and briefly centrifuge the Master Mix and Calibrators' tubes.
2. Add 30 µl of Master Mix into PCR tubes.
3. Add 20 µl of the isolated nucleic acid sample or 20 µl of Calibrator into the individual PCR tubes and mix by pipetting. The total reaction mix volume will be 50 µl. *It is necessary to keep all components at (2 – 8) °C during the PCR preparation. The customer must use his own negative control in the form of isolate of negative clinical material in each test. All 4 Calibrators must be used for setting the standard curve for quantitative detection.*
4. Close the tubes, centrifuge shortly, insert them into the device and let them amplify according to the following PCR profile.
Be very careful when handling the Calibrators or the clinical material; incorrect handling could result in contamination and the consequent impairment of the kit components! The manufacturer is not responsible for the kit impairment due to incorrect handling.

AMPLIFICATION PROFILE

Step	Temperature	Time	Data Collection	Cycles
Hold	37 °C	2 min		1
Hold	50 °C	5 min		1
Hold	95 °C	20 s		1
PCR	55 °C	20 s		3
	72 °C	5 s		
	95 °C	5 s		
PCR	95 °C	5 s		45
	58 °C	35 s	FAM + HEX	

INSTRUMENTS

GeneProof HIV type 1 (HIV-1) PCR Kit is designed for use with real-time devices from various manufacturers:

croBEE Real-Time PCR System

- AMPLilab Real-Time PCR System
- Applied Biosystems 7300 / 7500 Real-Time PCR System
- AriaMx Real-Time PCR System
- CFX Connect™ / CFX96™/ Dx Real-Time PCR Detection System
- LightCycler® 480
- LineGene 9600 Plus

- Mic qPCR Cycler
- QuantStudio™ 3 / 5 Real-Time PCR System
- Rotor-Gene 6000 / Q*
- SLAN® Real-Time PCR System

Required channels: FAM, HEX

GeneProof diagnostic kits are continually verified with various types of devices. Current list is available at www.geneproof.com or can be requested at support@geneproof.com

*NOTE: For Rotor-Gene instruments, the third channel must also be set to switch on.

CLINICAL SAMPLE ANALYSIS EVALUATION

Channel FAM (HIV-1)	Channel HEX (IC)	Result	Interpretation
+	+	Valid	HIV-1 positive
+	-	Valid	HIV-1 positive
-	+	Valid	HIV-1 negative
-	-	Invalid	-

QUANTITATIVE DETECTION EVALUATION

Use the following formula to calculate the virus concentration in IU/ml for **manual** extraction (using GeneProof PathogenFree RNA Isolation Kit):

$$\text{IU/ml} = \frac{\text{SC} \times \text{EV}}{\text{IV}}$$

SC - Sample concentration (IU/μl)

EV - Elution volume (μl)

IV - Isolation volume (ml)

To easily calculate pathogen concentrations using manual or automated extraction, you can use the calculator at www.geneproof.com

WARNING

A single valid Instruction for Use for a specific kit is included in the package or to be requested for the particular lot from the manufacturer. Use only combination of components from the specific PCR kit. The kit should be disposed of after use according to the current legal regulations considering the fact, that the kit doesn't contain any dangerous, infectious or toxic components that would be subject to special safety regulations, and the packaging materials are made of paper and polypropylene. If you have any questions, please contact our Customer care.

Customer care and Technical support

Tel.: +420 730 176 222

Email: support@geneproof.com

Orders

Tel.: +420 543 211 679

Email: sales@geneproof.com



GeneProof a.s.

Vídeňská 101/119 / Dolní Heršpice / CZ-619 00 Brno / +420 543 211 679 / info@geneproof.com

Version: IFU_0125_A01_1.0, Effective date: 9. 5. 2022

