

# Hot MolTaq 16S/18S

## DNA-free Taq DNA polymerase

From *Thermus aquaticus* BM, recombinant (*E.coli*)

Deoxynucleoside-triphosphate:  
DNA deoxynucleotidyltransferase, *EC 2.7.7.7*

**- For research use only -**



© 2017 - 2020 Molzym, all rights reserved

**Version 02**

Date of first release: 01/2017

Last update: 10/2020

**Contents**

<b>Content</b> .....	<b>Page</b>
Product Information .....	2
PCR Protocol .....	6
Troubleshooting .....	10
Related Products .....	11
Contact .....	12

**Product Information**

**Contents**

<b>Hot MolTaq 16S/18S (Taq DNA polymerase), microbial DNA-free</b>	<b>100 units</b>	<b>500 units</b>
<i>Hot MolTaq 16S/18S (Taq DNA polymerase, 5U/μl)</i>	1x 100U	1x 500U
<i>10x PCR buffer (1.5mM Mg<sup>2+</sup> final concentration)</i>	1x 1.5ml	2x 1.5ml
<i>10x KCl buffer basic (without Mg<sup>2+</sup>)</i>	1x 1.5ml	2x 1.5ml
<i>100mM MgCl<sub>2</sub></i>	1x 1.5ml	1x 1.5ml
<i>DNA-free water (PCR grade)</i>	2x 1.25ml	10x 1.25ml

<b>Catalog Numbers (Cat. No.):</b>	<b>P-080-0100</b>	<b>100 units</b>
	<b>P-080-0500</b>	<b>500 units</b>

**Buffers and Reagents**

- PCR buffer:** 10 x concentrated; contains ammonium sulfate (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; contains 1.5 mM Mg<sup>2+</sup> final concentration, microbial DNA-free
- KCl buffer basic:** 10 x concentrated, contains potassium chloride (KCl); without Mg<sup>2+</sup> ; microbial DNA-free
- MgCl<sub>2</sub>:** 100 mM, microbial DNA-free
- DNA-free water:** PCR grade, microbial DNA-free

**Application of the Product:**

- 16S, 18S and 23S rRNA gene amplification (without background)
- Routine PCR
- Real-Time PCR
- GC-Rich PCR
- Nested PCR
- MultiPlex PCR
- Colony PCR
- DNA labeling reactions

**For research use only**

## **Additional Equipment and Reagents Required:**

- Template DNA, gene-specific PCR primer pair
- dNTPs, PCR grade
- PCR tubes, DNA-free
- DNA-free polypropylene tubes for preparation of mastermixes and dilutions (e.g., Biosphere®, Sarstedt, Germany)
- PCR cycler (e.g., Mastercycler Gradient, Eppendorf, Germany)
- UV irradiated workstation (e.g., Airclean 600, StarLab, Germany)

## **Product Description**

**Hot MolTaq 16S/18S** is a mixture of a genetically engineered, thermostable, DNA-dependent DNA polymerase (EC 2.7.7.7.) originally isolated from *Thermus* sp. and an aptamer-based inhibitor that reversibly binds to the enzyme. *Hot MolTaq 16S/18S* is DNA-free and thus especially useful for 16S, 18S and 23S rDNA gene amplification and cloning, because the enzyme preparation is free of genomic and plasmid DNA from the production strain. The polymerase activity is inhibited at ambient temperatures. This allows reaction set-up at room temperature without loss of specificity or sensitivity.

### **The enzyme has the following activities:**

5'-3' polymerase (60 to 150 nucleotides/s; approx. 1 kb/min), 5'-3' exonuclease (strand displacement) and 3' terminal deoxynucleotidyl transferase (addition of single dATP to the duplex DNA). *Hot MolTaq 16S/18S* incorporates modified nucleotides (dNTP $\alpha$ S, c7GTP, biotin-11-dUTP, digoxigenin-11-dUTP and fluorescein-12-dUTP but not biotin-16-dUTP) at high rates. The error frequencies for misincorporation and frameshift mutations are 10<sup>-4</sup> and 2 x 10<sup>-5</sup>, respectively. Only the last three bases adjacent to the 3' end of the primer need to be correctly base-paired in order to initiate polymerization. The 5' region of the primer is less sensitive to base mismatches. Thus, new restriction sites can be easily introduced into an amplification product. Up to 9 kb can be amplified from lambda DNA and up to 5 kb from genomic DNA.

### **Unit definition:**

One unit of enzyme incorporates 10 nmoles of dNTPs into an acid-insoluble form in 30min at 72°C under standard assay conditions, using a DNA template. Specific activity is 50.000U/mg protein.

### **Heat activation:**

The enzyme is activated during normal cycling conditions. An extra activation step is not required.

## Stability

Stable for 20 months from the date of manufacturing under proper storage condition. Guarantee for full performance of reagents and buffers is given through the expiration date printed on the label at the outer box, if the packed material is undamaged upon arrival and the reagents are unopened.

## Enzyme Properties

---

Volume activity	5U/μl
Error rate*	approx. $1.3 \times 10^{-5}$
Optimal enzyme concentration	1 to 5U per 25μl reaction
Optimal pH	approx. 9 (adjusted at 20°C)
Optimal polymerization temperature	around 72°C
Optimal Mg <sup>2+</sup> concentration	varies from 1.5 to 5mM
Standard Mg <sup>2+</sup> concentration	1.5mM when using 200μM dNTP each
PCR product size	optimized for up to 5kb
PCR cloning	TA-cloning; addition of A-overhangs ("extendase activity")
Incorporation of modified nucleotides	accepts modified nucleotides like radiolabeled nucleotides, DIG-dUTP, biotin-dUTP
Thermostability	Over 80% activity after 30 cycles (1min 95°C, 1min 37°C, 3min 72°C)

---

\* According to the *lacI* assay (Frey, B. and Suppmann, B. (1995) *Biochemica* 2, 8-9.)

## Packaging, Storage and Handling

The purification of the *Hot MolTaq 16S/18S* and its confectioning are done under standard precautions for the avoidance of air-borne and handling-based DNA contaminations. The *Hot MolTaq 16S/18S* is supplied in a DNA-free screw cap vial. Store all vials in the kit at -15 to -25°C upon receipt.

**Important: Replace *Hot MolTaq 16S/18S* to the freezer (-15 to -25°C) after handling.**

For usage, the *10x PCR buffer*, *10x KCl buffer basic*, *100mM MgCl<sub>2</sub>* and the *DNA-free water* of the kit are thawed at room temperature (+18 to +25°C). After use, these components can be stored in the refrigerator (+4 to +12°C) for further use at the same day, but should be replaced to -15 to -25°C for longer storage.

## Avoidance of DNA Contamination:

PCR analysis demands special care with respect to avoidance of contamination from exogenous sources. Take care to separate places of DNA preparation from places where PCR reagents are handled for runs, in particular preparation of mastermixes, pipetting into PCR tubes and performance of PCR runs. Wear a protective lab coat, sterile protective gloves and sleeve covers, protective goggles and a disposable mask at any handling step, also during DNA preparation. Use only DNA-free disposables. In particular, use only PCR strips or tubes and pipette filter tips from manufacturers guaranteeing DNA-free products. Please contact Molzym for further information regarding our products and other suppliers of DNA-free plastic consumables.

Controls and DNA sample template:

Generally, for each analysis, run positive and negative controls to check for proper performance of the reaction and sterility of reagents and buffers used.

PCR negative controls using *DNA-free water* instead of template DNA are used for analysis of contamination of microbial DNA in the purified final mastermix. Ensure that negative controls are prepared first and the tube is sealed before positive controls and samples are pipetted. Take care to avoid cross contamination by aerosols.

For DNA sample template preparation, make sure to use kits guaranteeing microbial DNA-free reagents, buffers and matrices like Molzym's products (e.g. MolYsis™ kits) for the sensitive detection of pathogens in blood, other body fluids and other matrices.

PCR positive controls are run using known amounts of genomic DNA. This can be extracted and purified by Molzym's *PrestoSpin D Bug* or *PrestoSpin D Fungi* (see page 11) from a microorganism for the specific assay. Alternatively, use Molzym's Positive Control DNA (see page 11) for gam-positive bacteria and yeasts.

### **Quality Control and Specifications**

Each lot of *Hot MolTaq 16S/18S* Taq DNA polymerase is produced and quality-controlled according to Molzym's recorded quality management system.

The *Hot MolTaq 16S/18S* DNA polymerase is tested for performance in PCR runs, using different primer pairs for amplification of sequences between 0.2 and 1 kb.

Guarantee of the DNA-free product is given for the absence of signals in PCR negative controls at a rate of  $\geq 97\%$  for up to 40 PCR cycles (provided the avoidance of contamination). DNA-free *Hot MolTaq 16S/18S* is defined as giving no bacterial and fungal DNA-specific signal. In negative control runs, the absence of banding in gel electrophoretic analysis is demonstrated.

Positive controls are done using known amounts of genomic DNA from *Staphylococcus aureus* or other bacteria. Alternatively, Molzym's Positive Control DNA (P1) is used as PCR positive control (cat. no. S-200-050).

### **Material Safety Data Sheet (MSDS)**

Please contact Molzym for an actual MSDS (according to Regulation (EC) No. 1272/2008) of this product:

**Tel.:** +49(0)421 69 61 62 0 • **E-Mail:** info@molzym.com

### **Patents/Disclaimer**

Some applications in which this product can be used may be covered by patents issued and applicable in certain countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a patent license depending upon the particular application and country in which the product is used.

## PCR Protocol

### Standard PCR Procedure

#### General Considerations

The optimal conditions (incubation times and temperatures, concentrations of enzyme, template DNA,  $Mg^{2+}$ ) depend on the system used and have to be determined individually. In particular, the  $Mg^{2+}$  concentration and the amount of enzyme used per assay should be titrated for optimal efficiency of DNA synthesis. As a starting point, use the following guidelines:

- **Optimal enzyme concentration:** 1 to 5U/25 $\mu$ l
- **Optimal  $Mg^{2+}$  concentration** can vary from 1.5 to 5mM; in most cases  $Mg^{2+}$  concentration of 1.5mM will produce satisfactory results when using 200 $\mu$ M dNTPs (each).  
The 10x PCR buffer supplied contains 1.5mM  $Mg^{2+}$ . The *10x KCl buffer basic* supplied contains no  $Mg^{2+}$ .
- **dNTP concentration:** Always use balanced solutions of all four dNTPs. The final concentration of each dNTP should be between 50 and 500 $\mu$ M; the most commonly used concentration is 200 $\mu$ M. Increase concentrations of  $Mg^{2+}$  when increasing the concentration of dNTPs.
- **Template concentration:**
  - Genomic DNA template: 10 to 250ng (e.g. human genomic DNA template)
  - Plasmid DNA template: 0.1 to 15ng,
  - Bacterial DNA template: 10pg to 10ng (*10x KCl buffer basic*, 3mM  $Mg^{2+}$  is recommended)
- **Dilution buffer for template DNA:** The optimal dilution buffer is either simply sterile double distilled water or 5-10 mM Tris (pH 8-9).  
Storage of DNA should be in 10mM Tris pH 8.0 to 8.5 or TE buffer (10mM Tris pH 8.0; 1mM EDTA).

### PCR and Real-Time PCR Amplification of the 16S/18S rDNA

Follow the instructions for avoidance of DNA contamination as indicated on page 4.

#### Fluorescent dye:

When using SYBR Green I as fluorescent dye the *10x KCl PCR buffer basic* and 100mM  $MgCl_2$  (3mM  $Mg^{2+}$  final concentration) are recommended for better performance.

For detection of bacteria in whole blood and other body liquids using Molzym's kits MolYsis™ kits, we have made good experience following the protocol B described on page 8.

## Preparation of the Reaction Mixtures

Take care that all handling is done in a DNA-free environment. **Work under a UV irradiated workstation.** Make sure that plastic consumables (including PCR tubes, pipette tips, polypropylene tubes) are free of contaminating microbial DNA when used in combination with the amplification reaction mixture. With each series of PCR, run a negative control and a positive control. Work according to the sequence of steps of protocol A or protocol B below (25µl final reaction volume per assay):

### Protocol A:

#### **Hot MolTaq 16S/18S with 10x PCR buffer (1.5mM Mg<sup>2+</sup> final concentration)**

##### 1. Reagents Pretreatment

Thaw *10x PCR buffer* and *DNA-free water* at room temperature (+18 to +25°C). Vortex for a few seconds to mix and briefly centrifuge vials to clear the lids. Place the *Hot MolTaq 16S/18S* in the rack of the other components. After handling, place the *Hot MolTaq 16S/18S* vial and other reagents back to the freezer (-15 to -25°C).

##### 2. Preparation of the Mastermix

Pipette the following components (table below) in a DNA-free 1.5ml polypropylene tube. In case of more reactions, prepare a mastermix containing all PCR reagents in proportions as indicated in the table below. Please note: The volume of the *DNA-free water* added is corrected for the variable volumes of the components.

Vortex the prepared mastermix for few seconds to mix and briefly centrifuge to clear the lid of the tube.

Components	Volume for 1 reaction	Final concentration
<i>DNA-free water</i> (PCR-grade; supplied)	add up to 25µl	
<i>10x PCR buffer</i> (supplied)	2.5µl	1x (1.5mM Mg <sup>2+</sup> )
dNTP mix (10mM)	0.5µl	200µM (of each dNTP)
Forward primer	variable	0.1 to 0.6µM
Reverse primer	variable	0.1 to 0.6µM
<i>Hot MolTaq 16S/18S</i> (supplied)	0.2 to 1.0µl	1 to 5U/reaction
Template (sample DNA, DNA-free water or positive control DNA)	x* µl (variable)	variable
<b>Final volume</b>	25µl	

\* x µl = volume of template DNA added

3. Pipette aliquots of the mastermix to each PCR tube and add the supplied *DNA-free water* for a PCR negative control.
4. Place the PCR tubes in another rack designated for transport to the place where DNA is handled.  
Add x µl of the DNA template (sample DNA or positive control DNA) into PCR tubes containing the mastermix.
5. Seal tubes and place in a PCR machine immediately.
6. Start the program of the specific assay ('PCR Thermocycling Conditions', page 9).

**Protocol B:**

**Hot MolTaq 16S/18S with 10x KCl buffer basic (without Mg<sup>2+</sup>)**

**1. Reagents Pretreatment**

Thaw 10x KCl buffer basic, 100mM MgCl<sub>2</sub> and DNA-free water at room temperature (+18 to +25°C). Vortex for a few seconds to mix and briefly centrifuge vials to clear the lids. Place the Hot MolTaq 16S/18S in the rack of the other components. After handling, place the Hot MolTaq 16S/18S vial and other reagents back to the freezer (-15 to -25°C).

**2. Preparation of the Mastermix**

Pipette the following components (table below) in a DNA-free 1.5ml polypropylene tube. In case of more reactions, prepare a mastermix containing all PCR reagents in proportions as indicated in the table below. Please note: The volume of the DNA-free water added is corrected for the variable volumes of the components.

Vortex the prepared mastermix for few seconds to mix and briefly centrifuge to clear the lid of the tube.

Components	Volume for 1 reaction	Final concentration
DNA-free water (PCR-grade; supplied)	add up to 25µl	
10x KCl buffer basic (supplied, without Mg <sup>2+</sup> )	2.5µl	
100mM MgCl <sub>2</sub>	0.375 to 1.25µl	1.5 to 5mM
dNTP mix (10mM)	0.5µl	200µM (of each dNTP)
Forward primer	variable	0.1 to 0.6µM
Reverse primer	variable	0.1 to 0.6µM
Hot MolTaq 16S/18S (supplied)	0.2 to 1.0µl	1 to 5U/reaction
Template** (sample DNA <sup>#</sup> , DNA-free water or positive control DNA)	x* µl (variable)	variable
<b>Final volume</b>	<b>25µl</b>	

\* x µl = volume of template DNA added

The following table shows an overview of the final Mg<sup>2+</sup> concentration and the corresponding volume of the 100mM MgCl<sub>2</sub> reagents:

Final Mg <sup>2+</sup> concentration (25µl PCR volume)								
Final concentration (mM)	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0
Required volume (µl) of 100mM MgCl <sub>2</sub> per reaction	0.375	0.5	0.625	0.75	0.875	1.0	1.125	1.25

- Pipette aliquots of the mastermix to each PCR tube and add the supplied DNA-free water for a PCR negative control.
- Place the PCR tubes in another rack designated for transport to the place where DNA is handled.  
Add x µl of the DNA template (sample DNA or positive control DNA) into PCR tubes containing the mastermix.
- Seal tubes and place in a PCR machine immediately.
- Start the program of the specific assay ('PCR Thermocycling Conditions', page 9).



## PCR Thermocycling Conditions

Place the samples in the thermal cycler, and start cycling. For single extension time use the protocol below:

Method	Target Temperature	Incubation time	Cycles
Initial Denaturation*	95°C	1 to 3min	1
Cycling	Denaturation*	15 to 30s	25 to 40
	Annealing#	30 to 60s	
	Elongation <sup>+</sup>	30s to 1min	
Final elongation	72°C	10min	1
Cooling	4°C	unlimited time	∞

\* The denaturation temperature can vary from 92 to 95°C.

# Optimal annealing temperature depends on the melting temperature of the primers and on the system used.

\* For PCR products longer than 1kb, use an extension time of approximately 1min per kb DNA.

After amplification, the PCR products can be stored overnight at +4 to +12°C, or at -15 to -25°C for longer storage. Avoid frequent freeze-thaw cycles, because this may result in a decay of the DNA (in particular at low DNA concentrations).

### Possible further procedures:

Check the PCR product on an agarose gel for size and specificity using an appropriate size marker.

---

Please address any questions relating this product to the support hotline:

**E-Mail:** support@molzym.com • **Tel.:** +49(0)421 69 61 62 0

## Troubleshooting

Possible cause	Recommendation
<b><i>Little or no PCR product</i></b>	
DNA template problems	Check quality and concentration of template: <ul style="list-style-type: none"> <li>• Analyze an aliquot with photometric methods or on an agarose gel to check for possible degradation.</li> <li>• Make a control reaction on template with an established primer pair or PCR system.</li> <li>• Check or repeat purification of template.</li> </ul>
Enzyme concentration too low	<ul style="list-style-type: none"> <li>• Increase the amount of polymerase in 0.5U steps.</li> </ul>
Mg <sup>2+</sup> concentration too low	Increase the Mg <sup>2+</sup> concentration in steps of 0.25mM. (1.5mM Mg <sup>2+</sup> is the final concentration in the mastermix.)
Cycle conditions not optimal	<ul style="list-style-type: none"> <li>• Decrease annealing temperature.</li> <li>• Increase cycle number.</li> <li>• Make sure that the final elongation step was carried out.</li> </ul>
Primer design not optimal	Design alternative primers
Primer concentration not optimal	<ul style="list-style-type: none"> <li>• Both primers must have the same concentration.</li> <li>• Titrate primer concentration (0.1 – 0.6 μM)</li> </ul>
Primer quality or storage problems	<ul style="list-style-type: none"> <li>• If you use an established primer pair, check performance on an established PCR system (control template).</li> <li>• Make sure that the primers are not degraded.</li> <li>• Always store primers at -15 to -25°C.</li> </ul>
<b><i>Multiple bands or background smear</i></b>	
Annealing temperature too low	Increase annealing temperature according to the primer length. <ul style="list-style-type: none"> <li>• Review primer design.</li> </ul>
Primer design or concentration not optimal	<ul style="list-style-type: none"> <li>• Titrate primer concentration (0.1 – 0.6μM).</li> <li>• Both primers must have the same concentration.</li> <li>• Perform nested PCR with nested primers.</li> </ul>
Difficult template e.g. GC-rich template	Perform PCR with PCR enhancer.
DNA template problems	Use serial dilution of template.
<b><i>PCR products in negative control experiments</i></b>	
Carry-over contamination	<ul style="list-style-type: none"> <li>• Set up PCR reactions in an area separate from that used for PCR product analysis (UV workstation e.g., Airclean 600, StarLab, Germany) and use sterile disposables (gloves, sleeves and bouffant covers).</li> <li>• Exchange all reagents, use DNA-free reagents (e.g., Molzym´s DNA-free PCR reagents).</li> <li>• Use aerosol-resistant and microbial DNA-free pipette tips.</li> <li>• Use DNA-free consumables (PCR tubes, polypropylene tubes (e.g., 1.5ml micro tubes, Biosphere®, Sarstedt, Germany)</li> <li>• Use dUTP (600μM) instead of dTTP (200μM) in combination with thermolabile UNG (1U/50μl reaction) and increase Mg<sup>2+</sup> concentration up to 4mM at most.</li> </ul> <p><b>Further information see part Product Information – Avoidance of DNA contamination, page 4</b></p>
<b><i>Specific problems in RT-PCR application</i></b>	
No product, additional bands, background	<ul style="list-style-type: none"> <li>• The volume of cDNA template (RT-reaction) should not exceed 10% of the smear final volume of the PCR reaction.</li> <li>• Follow trouble shooting above.</li> <li>• Increase Mg<sup>2+</sup> by titration in steps of 0.25mM.</li> </ul>

**Related Products**

<b>Product</b>	<b>Contents</b>	<b>Cat. No.</b>
<b>Isolation of microbial DNA from whole blood and other body liquids (DNA-free reagents)</b>		
<b>MolYsis™ Complete5</b>	50 microbial DNA isolations from samples	D-321-050
≤1ml and 5ml sample volumes	100 microbial DNA isolations from samples	D-321-100
<b>PCR Mastermix, DNA-free (incl. MolTaq 16S/18S)</b>		
<b>Mastermix 16S Complete</b>	100 reactions	S-020-0100
Universal 16S rDNA PCR and Real-Time PCR assay for detection of bacteria	250 reactions	S-020-0250
	1000 reactions	S-020-1000
<b>Mastermix 16S Primer</b>	100 reactions	S-021-0100
PCR assay for universal PCR detection of bacteria	250 reactions	S-021-0250
	1000 reactions	S-021-1000
<b>Mastermix 18S Complete</b>	100 reactions	S-070-0100
Universal 18S rDNA PCR and Real-Time PCR assay for detection of fungi.	250 reactions	S-070-0250
	1000 reactions	S-070-1000
<b>Mastermix 16S/18S Dye</b>	100 reactions	S-030-0100
Premixed reagents and fluorescent dye for Real-Time PCR with custom primers	250 reactions	S-030-0250
	1000 reactions	S-030-1000
<b>Mastermix 16S/18S Basic</b>	100 reactions	S-040-0100
Premixed reagents for PCR analysis with custom primers	250 reactions	S-040-0250
	1000 reactions	S-040-1000
<b>PCR Reagents</b>		
<b>DNA-free water (PCR grade)</b>	Microbial DNA-free PCR grade water for molecular analysis 10x 1.7ml	P-020-0003
<b>Positive Control DNA (P1)</b>	Positive run control in universal PCRs and Real-Time PCRs. 50 reactions	S-200-050
<b>Panfungal Sequencing Primer</b>	For sequencing analysis of PCR amplicons produced with <i>Mastermix 18S Complete</i> and <i>Mastermix 18S Primer</i>	S-785-100
<b>Set of Eubacterial Sequencing Primers</b>	For sequencing analysis of PCR amplicons produced with <i>Mastermix 16S Complete</i> and <i>Mastermix 16S Primer</i> . 100 reactions, each primer	S-775-100

**Order Hotline:**

**Tel.:** +49(0)421 69 61 62 0 • **Fax:** +49(0)421 69 61 62 11 • **E-Mail:** order@molzym.com

**Contact**



**Molzym**

**Molzym GmbH & Co. KG**

Mary-Astell-Str. 10

28359 Bremen, Germany

**Tel.:** +49(0)421 69 61 62 0 • **Fax:** +49(0)421 69 61 62 11

**E-Mail:** [info@molzym.com](mailto:info@molzym.com); [order@molzym.com](mailto:order@molzym.com) • **Web:** [www.molzym.com](http://www.molzym.com)