

Product Name

Coxsackievirus A6/A10/A16/Enterovirus 71 RNA Diagnostic Kit(PCR-Fluorescence)

Specification

12 reactions/kit

Intended Use

The kit uses real-time multiplex fluorescent PCR technology, combined with an Automated Fully Enclosed qPCR Instrument, for the in vitro qualitative detection of Coxsackievirus A6 (CA6), Coxsackievirus A10(CA10), Coxsackievirus A16(CA16) and Human Enterovirus 71(EV71)and other serotype nucleic acids. The test results can be used for the auxiliary diagnosis of enterovirus.

The hand, foot and mouth disease is a common childhood disease caused by enteroviruses. The incidence is mainly in children aged 5 and below, especially children aged 3 and below has the highest incidence. Most patients will have fever, skin rashes and ulcers in the hands, feet, mouth, etc., A few patients may have fatal complications such as meningitis, encephalitis, myocarditis, pneumonia, and severe cases can lead to death. The enteroviruses that cause hand, foot and mouth disease are more common with CA16 and EV71 . The methods of Laboratory diagnostic is mainly include virus isolation and nucleic acid detection.

The kit uses molecular biology methods and requires experimental operators to have received professional training in gene amplification or molecular biology method detection, who also have relevant experimental operating qualifications and the laboratory should have reasonable biosafety facilities and protective procedures.

Test Principle

The Automated Fully Enclosed qPCR Instrument uses iCassette technology to automatically perform nucleic acid extraction, amplification, data reading and result analysis. The instrument includes two parts: nucleic acid extraction and multiple fluorescent PCR. At the same time, the built-in QR code scanner can automatically identify the execution process corresponding to this kit. The uniquely designed software is used to execute the extraction and PCR process of the kit, display the results and analyze the results. The kit contains 12 disposable CA6/CA10/CA16/EV71 iCassettes. The reagents stored in the iCassette include nucleic acid extraction reagents and PCR reaction reagents, which are suitable for Automated Fully Enclosed qPCR Instrument. Since the Automated Fully Enclosed qPCR Instrument is fully enclosed and automatically performs nucleic acid extraction and PCR reactions, it reduces direct cross-contamination of samples. For a complete description of the instrument, please refer to the user manual of the corresponding instrument.

The kit includes nucleic acid extraction reagents and PCR reaction reagents. The nucleic acid extraction reagents contain internal standards which are used to monitor the entire process of nucleic acid extraction, confirm whether the sample is adequately processed and monitor whether there are PCR inhibitors to avoid PCR false negatives. At the same time, before the start of the PCR reaction, the Automated Fully Enclosed qPCR Instrument will first test the fluorescent signal of the probe to monitor whether the PCR reagents are fully reconstituted, whether the samples are added to the PCR reaction tube and the integrity and stability of the probe.

The kit uses specific region ofCA6, CA10, CA16 and EV71 as the detection target area(CA6 is labeled with ROX, CA10 is labeled with CY5, CA16 is labeled with HEX, EV71 is labeled with FAM, and the internal standard is labeled with CY5.5).

Reagents are designed with specific primers and fluorescent probes to perform one-step RT-PCR which can simultaneously detect four serotypes of CA6 , CA10, CA16 and EV71 in one PCR reaction tube.(CA6, ROX labeled; CA10, Cy5 labeled; CA16, HEX labeled; EV71, FAM labeled; Internal control, Cy5.5 labeled.)

COMPONENTS

Kit Components		Content	Quantity	
CA6/CA10/CA16/EV71 iCassette (12 pcs)	Nucleic acid extraction reagent	Magnetic beads	Magnetic microspheres	11μL/iCassette
		CA6/CA10/CA16/EV71 lyophilization A	Proteinase K	1pc/iCassette
		CA6/CA10/CA16/EV71 lyophilization B	Pseudovirus lyophilization contains internal control fragments	1pc/iCassette
		Virus lysate solution	Guanidine isothiocyanate	450μL/iCassette
		Virus binding solution	Guanidine isothiocyanate	175μL/iCassette

	Virus washing solution	Sodium chloride	900μL/iCassette
	Virus elution	Tris-HCl	100μL/iCassette
	CA6/CA10/CA16/EV71 cosolvent solution	PCR Buffer, MgCl ₂	35μL/iCassette
CA6/CA10/CA16/EV71 PCR reagent	CA6/CA10/CA16/EV71 lyophilization	Specific primer probe, dNTP, enzyme	1 pc/iCassette
Control	CA6/CA10/CA16/EV71 positive control	pseudovirus containing target fragment	1 tube (200μL)
	CA6/CA10/CA16/EV71 negative control	Sterilizing purified water	1 tube (200μL)

Note:Components in kits with different batch numbers are not interchangeable

Storage condition and Shelf life

- The kit can be stored at 2-8℃ and the shelf life is 9 months.
- The transportation temperature range of the kit should be kept at 2-8℃.
- Please do not open the iCassette lid before adding the sample. If you open the iCassette lid, it should be used within 30 minutes.

Applicable instruments

Automated Fully Enclosed qPCR Instrument: Galaxy Nano, Galaxy Lite and Galaxy Pro.

Materials Required but Not Provided

- Automated Fully Enclosed qPCR Instrument : Galaxy Nano, Galaxy Lite or Galaxy Pro.
- Leak-proof, sterile, screw-capped specimen collection containers.
- Disposable gloves, eye protection, laboratory coats, and labels or permanent marking pen.
- Timer.
- Pipettes.
- Sterile pipette tips.

Sample Requirements

- Applicable sample types: oropharyngeal swab
- Sample collection

2.1 Sampling equipment requirements

2.1.1.Swab: The head used to sample should be the synthetic fiber (eg, polyester fiber), and use aluminum or plastic as a rod.

2.1.2. Sampling tube: Contains 3mL physiological saline or sampling liquid (containing protein stabilizers, antibiotics, buffer to prevent bacteria and fungal growth).

2.2. Sample collection: Use a swab to wipe the posterior pharyngeal wall and tonsils on both sides with moderate force, avoid touching the tongue; quickly put the swab into a sampling tube for storage.

3. Sample storage and transportation

All collected respiratory samples should be airtight and sent for inspection immediately, and stored at 2–8℃ for testing within 12 hours, or placed at -20±5℃ as soon as possible, but not longer than 1 month.

4. Principles of Biosafety Protection

All operations should comply with local relevant laws and regulations.

Test Method and Operation
1. Prepare CA6/CA10/CA16/EV71 iCassette

1.1Processing samples or control materials in the samples preparation room,first vortex the collection tube or controls for 10-15 seconds.

1.2 Open the package of CA6/CA10/CA16/EV71 lyophilization, observe whether the CA6/CA10/CA16/EV71 lyophilization is intact, replace the empty PCR tube on the iCassette with the PCR tube containing the CA6/CA10/CA16/EV71 lyophilization, and make sure the PCR tube is screwed up.

1.3 Open the lid of the iCassette, and pipette 200 μl of sample or controls to the CA6/CA10/CA16/EV71 iCassette sample compartment as shown in Figure 1. slowly, close the lid tightly.(Note: After the sample is added to the iCassette, the test needs to start running within 30 minutes)

1.4 Place the iCassette into the instrument.




Figure 1 CA6/CA10/CA16/EV71 iCassette (Lateral View)

2. Test Operation

2.1 Turn on the power, press the switch button on front side of the instrument, and the blue light is on, i.e. "the instrument is on". After the instrument is turned on, click the control software icon on the desktop to enter the login interface.

2.2 Log in to the software for the first time with the administrator account (Admin/123456), click "OK" to complete the login.

2.3 Click the "Open" button in the initial interface to open the compartment door of the instrument.

2.4 Put the iCassette that has been added with sample into the instrument and click the "Close" button in the initial interface to close the compartment door. (For multi-channel instrument, iCassette tray holder is required to load the iCassette. Input the sample information at the corresponding channel position of the control software while loading the iCassette)

2.5 After the compartment door is closed, the instrument scans the corresponding QR code on the iCassette automatically.

2.6 After scanning, click the "Run" button, the software will automatically load the program corresponding to the QR code of the iCassette, and click "OK" to start the program.

2.7 After the program starts , the progress of the instrument running will be displayed in the main interface.

2.8 After the amplification is completed, the compartment door will open automatically.

For the detailed steps of test operation, please refer to user manual of the instrument.

3. Result Analysis

After the experiment is finished, the instrument will automatically save the results, and output the Ct value and amplification curves with results interpretation in the interface.

Quality Control

The results of positive control are positive in all five channels (FAM, HEX, ROX, CY5 and CY5.5) .

The results of negative control are negative in four channels (FAM, HEX, ROX and CY5), CY5.5 channel is positive.

Internal control: The CY5.5 channel is positive (the internal control is specifically competing with the sample. When the sample is positive, the internal control may be tested negative. When the sample is negative, the internal control must be positive. Otherwise it may be a sampling failure or test failure).

The above requirements must be met at the same time in the same experiment, otherwise, this experiment is invalid and needs to be repeated.

Reference Interval

Based on the analysis of clinical sample test results, using the ROC curve method, the final determined Ct positive threshold value for this reagent kit is 36.

Interpretation of Test Results

The sample to be tested will be determined according to the standards in the table below:

FAM Channel	HEX Channel	ROX Channel	CY5 Channel	CY5.5 Channel	Results
√	×	×	×	√	EV71 positive(Figure 2)
×	√	×	×	√	CA16 positive(Figure 3)
×	×	√	×	√	CA6 positive(Figure 4)
×	×	×	√	√	CA10 positive(Figure 5)
×	×	×	×	√	Negative(Figure 6)

√	√	√	√	√	Positive control.
×	×	×	×	√	Negative control.

Note: "√" indicates the result "has an obvious logarithmic amplification curve and the value of Ct is less than or equal to 36"; "×" indicates the result "No logarithmic amplification curve and the value of Ct is more than 36".




Figure 2 EV71 positive




Figure 3 CA16 positive




Figure 4 CA6 positive




Figure 5 CA10 positive




Figure 6 Negative

Limitations of test methods

- The test results of this reagent should be combined with the patient's clinical symptoms and other relevant medical examination results for comprehensive analysis, and should not be used as a basis for patient management alone.
- Unreasonable sample collection, transportation and processing, as well as improper experimental operation and experimental environment may lead to false negative or false positive results.
- Variations in the target sequence of the enterovirus or sequence changes caused by other reasons may cause false negative results.
- Because the positive rates of samples from different stages of the disease are inconsistent, the most suitable sample type for detection and the best sampling time after infection may not be confirmed. Therefore, collecting samples from the same patient multiple times will reduce the possibility of false negative results.
- Variations in the target sequence of enteroviruses or sequence changes caused by other reasons may lead to false negative results.
- Cross-contamination between samples, which may lead to false positive results.

Product performance index

- Analysis sensitivity

S1-S4 are derived from virus strain samples, which are serially diluted with RNA virus diluent and calibrated (see the table below for details).

Reference number	Type	Cont:copies/mL
S1	CA6	1.0×10 ³
S2	CA10	1.0×10 ³
S3	CA16	1.0×10 ³
S4	EV71	1.0×10 ³

The test results of the 4 sensitivity reference products are all positive and the coincidence rate was 100%.

2. Analysis accuracy

P1-P4 are derived from virus strain samples, which are serially diluted and calibrated with RNA virus diluent (see the table below for details).

Reference number	Type	Cont:copies/mL
P1	CA6	1.0×10 ⁴
P2	CA10	1.0×10 ⁴
P3	CA16	1.0×10 ⁴
P4	EV71	1.0×10 ⁴

The test results of the four positive reference products were all positive, and the compliance rate was 100%.

3. Analysis specificity

The following conclusions are drawn after the reagent specificity test: Pathogenic microorganisms with the same site of infection or similar symptoms of infection (Coxsackievirus A2, Coxsackievirus B3, Coxsackievirus B5, EB virus, Rubella virus, Measles virus, Influenza A virus, Parainfluenza virus, Influenza B virus, Cytomegalovirus, Respiratory syncytial virus, Mumps virus, *Klebsiella pneumoniae*, *Escherichia coli*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Salmonella*, *Shigella*, *Haemophilus influenzae*, *Pseudomonas aeruginosa* and human spore virus) have no cross reaction.

4. Precision

The coefficient of variation of intra-assay precision CV≤5%.

Precautions

- The kit must be used within the shelf life.
- Do not open the iCassette lid before use. If you open it, you must use it within 30 minutes.
- If you find that the iCassette is oscillated and leaked after adding the sample, do not continue to use it.
- Prevent sample contamination; experiment personnel should be well protected and wear disposable gloves and masks.
- Each iCassette is single-test, please do not reuse it.
- In order to avoid any potential biological hazards in the sample, the test sample should be regarded as an infectious substance and avoid contact with the skin and mucous membranes; sample manipulation and handling must comply with relevant local regulatory requirements
- This kit is only used for in vitro diagnosis.

References

- Epidemiological and etiological characteristics of hand, foot, and mouth disease in Ningbo, China, 2008–2011
- Ni H, Yi B, Yin J, et al. Epidemiological and etiological characteristics of hand, foot, and mouth disease in Ningbo, China, 2008–2011[J]. Journal of Clinical Virology the Official Publication of the Pan American Society for Clinical Virology, 2012, 54(4):342.
- He S J, Han J F, Ding X X, et al. Characterization of enterovirus 71 and coxsackievirus A16 isolated in hand, foot, and mouth disease patients in Guangdong, 2010[J]. International Journal of Infectious Diseases Ijidd Official Publication of the International Society for Infectious Diseases, 2013, 17(11):e1025.
- Chen Q, Hu Z, Zhang Q, et al. Development and evaluation of a real-time method of simultaneous amplification and testing of enterovirus 71 incorporating a RNA internal control system[J]. Journal of Virological Methods, 2014, 196(1):139-144.
- Reid S M, Mioulet V, Knowles N J, et al. Development of tailored real-time RT-PCR assays for the detection and differentiation of serotype O, A and Asia-1 foot-and-mouth disease virus lineages circulating in the Middle East[J]. Journal of Virological Methods, 2014, 207(207):146-153.
- Cheng H Y, Huang Y C, Yen T Y, et al. The correlation between the presence of viremia and clinical severity in patients with enterovirus 71 infection: a multi-center cohort study.[J]. BMC Infectious

Diseases, 2014, 14(1):417.

7.Zhang S1, Wang J2, Yan Q2, et al. A one-step, triplex, real-time RT-PCR assay for the simultaneous detection of enterovirus 71, coxsackie A16 and pan-enterovirus in a single tube.[J]. PloS One, 2014, 9(7):e102724.

8.Liu N, Xie J, Qiu X, et al. An atypical winter outbreak of hand, foot, and mouth disease associated with human- enterovirus 71, 2010. BMC Infect Dis, 2014, 4(14):123.


















Instruction Version

Version: A/1

Date of Issue: May,2019

Symbols

The following symbols may appear on the product labeling:

	In vitro diagnostic medical device		Do not re-use
	Use-by date		Consult instructions for use or consult electronic instructions for use
	Caution		Manufacturer
	Temperature limit		Batch code
	Contains sufficient for <n>		Keep dry
	Keep away from sunlight		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Biological risks
	Catalogue number		CE marking of conformity
	Authorized representative in the European Community		



IGENESIS(SHANGHAI)CO. , LTD.

Address: Floor 3, building 1, Lane 500, Furonghua Road, Pudong New Area, 201318 Shanghai, P.R. China.
Tel: +86-21-38016598



Lotus NL B.V.

Address: Koningin Julianaplein 10, 1e Verd,2595AA, The Hague, Netherlands.
Email: peter@lotusnl.com