**BIOSYNEX AMPLIQUICK Malaria**

**REF** 3120019\_SEC01/ 3120019\_SEC02

FR (FRANÇAIS/FRENCH).....1  
EN (ENGLISH) .....15

**CONTENTS**

## BIOSYNEX AMPLIQUICK Malaria

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## BIOSYNEX AMPLIQUICK MALARIA

PCR FOR THE QUALITATIVE DETECTION OF *PLASMODIUM SPP* AND IDENTIFICATION OF *PLASMODIUM FALCIPARUM* FROM WHOLE BLOOD SAMPLES.

*In vitro* diagnostic medical device for professional use.

**REF** 3120019\_SEC01/ 3120019\_SEC02

### 1. INTENDED PURPOSE

BIOSYNEX AMPLIQUICK Malaria is a molecular *in vitro* diagnostic test for the qualitative detection by qPCR of 5 Plasmodium species - *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale* and *Plasmodium knowlesi* (PAN) and identification of *Plasmodium falciparum* (FALCI). The test is an aid for the diagnosis of Malaria performed using a DNA extract obtained from whole blood sample. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

### 2. CLINICAL SUMMARY

Malaria is a life-threatening disease mostly spreads to people through the bites of some infected female *Anopheles* mosquitoes of the genus *Plasmodium*. It represents a global public health concern: being one of the first parasitic endemic species with a high mortality rate.<sup>1</sup> Five species have been identified as human pathogens: *P. falciparum*, *P. vivax*, *P. ovale*, *P. knowlesi* and *P. malariae*. Among these species, *P. falciparum* has the large highest prevalence (>80%).<sup>2</sup> In addition, an infection by this species requires treatment within 24 hours after the first symptoms appear, to prevent the parasitosis moving towards a serious, sometimes fatal form. This explains the importance of a specific and sensitive identification test for this species. The BIOSYNEX AMPLIQUICK Malaria test was developed to highlight a *Plasmodium spp* core-genome (PAN) sequence and a specific sequence of *P. falciparum* (FALCI).

### 3. TEST PRINCIPLE

The BIOSYNEX AMPLIQUICK Malaria kit is an *in vitro* diagnostic test, based on real-time Polymerase Chain Reaction (PCR) technology by fluorescent probe hydrolysis. The assay consists of a multiplex qPCR allowing the specific amplification and simultaneous detection of molecular targets carrying the sequences of interest: the 18S rRNA gene sequences from *Plasmodium spp* and *Plasmodium falciparum*, as well as the internal control, CIEZ, which allows to identify any PCR inhibition, thus excluding false negative results.

The increase in fluorescence signal is only detected if the target complementary sequence to the fluorescent probe is present in the sample and amplified. The fluorescent signal is therefore directly proportional to the amplification of the target during the amplification phase. The Cq value (quantification cycle) corresponds to the number of cycles at which the fluorescence starts to increase exponentially in contrast to the background noise.

The BIOSYNEX AMPLIQUICK Malaria kit must be used with whole blood samples. This kit can be used with purified DNA extracts or samples treated with the BIOSYNEX Blood Lysis kit (ref 3150077).

#### 4. REQUIRED MATERIAL

**Material provided:**

- 1 divisible 96-wells microplate prefilled with the Master mix
- 1 tube of positive control (CONTROL+, red cap)
- 1 tube of negative control (CONTROL-, green cap)
- 1 bag of optical cap strips.

**Material required but not provided:**

- DNA extraction kit
- Powder-free disposable gloves
- Pipettes and filtered tips
- PCR microplate or microtubes centrifuge
- Vortex mixer
- qPCR Thermal Cycler

The PCR thermal cycler used for the test must have the following main characteristics:

- Open system
- Real-time quantitative PCR assays
- Programmable Thermal Cycler block

Reference	Thermal cycling block
3120019_SEC01	0.1mL low profile
3120019_SEC02	0.2mL high profile

- Excitation source: LEDs, lamp or laser
- Filter sets (excitation/emission wavelengths) suitable for the detection of "reporter" fluorophores of the FAM, HEX and Cy5 probes.
- Connection with a computer using specific analysis software that allow the recovery of fluorescence data, absolute quantification assays, and the interpretation of results.

The kit has been developed and validated with the following real-time PCR thermal cyclers:

- CFX96 Touch™ Real-Time PCR Detection System (BioRad)
- CFX96 Opus™ Real-Time PCR Detection System (BioRad)
- Quant Gene 9600 (Bioer)
- LightCycler480 (Roche)
- QuantStudio5 System (Applied Biosystems)
- ViiA™ 7 Real-Time PCR System (Applied Biosystems).

If another thermal cycler is used, please perform your own validation of BIOSYNEX AMPLIQUICK Malaria before using the test. We recommend using qualified samples and the negative and positive controls provided in the kit.

#### 5. PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of many samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination. Pre-handling PCR, post-PCR and DNA extraction should therefore be performed in separate rooms. Wear disposable gloves

in each zone and change them before moving from one zone to another.

- The test and cap are for single use only. Do not reuse. Do not open the PCR tubes at the end of the test.
- Do not use the test if aluminum foil is opened or damaged. Once the aluminum foil is removed, use the test immediately.
- Do not use the kit if it arrived unfrozen.
- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.
- Protect the kit from light.
- Centrifuge the tubes before opening, opening them one after the other by closing them well between each one to avoid any contamination.
- The positive control (CONTROL+) contains significant amounts of DNA sequences. It can therefore potentially contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store this component outside the kit at the first opening of the kit.
- The negative and positive controls included in the kit mimic results obtained with negative or positive samples respectively. They must be used with each new run.
- To ensure that there is no contamination during handling, it is left to the user to include as part of their internal quality approach, at the frequency chosen, an additional negative control well in which 8µL of molecular biology grade water are added to the Master mix.
- When using negative and positive controls with a series of patients, it is recommended to first deposit the negative control, then deposit the patient samples and finish with the deposit of the positive control.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste. Consider your local regulation on biological waste disposal.
- The device contains material of bacterial or animal origin and should be handled with caution.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Malaria, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed.

## 6. KIT STORAGE

The kit is shipped frozen. Kit components must arrive frozen. Store the kit at a temperature of -20°C or below. Under these conditions, the reagents are stable until the expiry date indicated on the kit label. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light.

The product can undergo 15 freeze/thaw cycles without the shelf-life and qualities of the product are altered.

## 7. SAMPLE COLLECTION AND STORAGE

BIOSYNEX AMPLIQUICK Malaria allows the amplification of *Plasmodium* DNA from whole blood samples. The whole blood sample should be collected in EDTA or heparinized tube.

The test can be performed using purified DNA extract or sample treated with the BIOSYNEX Blood Lysis (ref 3150077).

The transport of clinical samples must comply with local regulations for the transport of infectious agents.

Blood samples can be stored up to 7h between +21°C and +30°C, and up to 5 days at +4°C before being processed. Samples can be stored for 3 years at -20°C. Repeated freeze/thaw cycles of the sample should be avoided (a maximum of 1 cycle is allowed).

Blood samples, after treatment with the BIOSYNEX Blood Lysis kit, can be stored up to 7 days between +4°C and +30°C or up to 3 months between -20°C and -80°C before being tested.

## 8. TEST PROCEDURE

### Extraction of nucleic acids

The nucleic acids extraction must be performed before the amplification protocol using a suitable extraction system (not provided within the kit). Please follow the manufacturer's instructions for use.

The BIOSYNEX AMPLIQUICK Malaria kit has been validated with the following extraction kits:

- BIOSYNEX Blood Lysis - Reference 3150077 (Biosynex)
- QIAamp DNA Blood Mini kit - Reference 51106 (Qiagen)
- QIAamp DNA mini kit - Reference 51306 (Qiagen)
- NucleoMag Dx Pathogen - Reference 744215.4 (Macherey Nagel)
- Whole blood genomic DNA/RNA extraction reagent kit - Reference MTQB036 (Singuway)
- MagNA Pure LC DNA Isolation Kit I - Reference 03 003 990 001 (Roche).

If another DNA extraction method or kit is used, please perform your own validation method of BIOSYNEX AMPLIQUICK Malaria before using the test (do not use the controls provided).

It is not necessary to extract the positive and negative controls with the nucleic acid extraction kit.

### Amplification protocol

Thawed samples should be gently homogenized before performing the test. Remove and defrost only the required number of Master mix strips.

1. Take one breakable plate and take the number of strips required. If less than a whole strip of 8 wells is needed, the strip can be cut horizontally using a pair of scissors.
2. Centrifuge the strips for a few seconds to collect any droplets on the well edges or on the seal.
3. Carefully remove and discard the aluminum foil. Once the aluminum foil is removed, use the strip within 2 hours at room temperature (+21°C ±2°C) or within 16 hours at +4°C.
4. Add 8 µL of sample or control.
5. Close the wells with the transparent caps provided.
6. Centrifuge the strips for a few seconds.
7. Place the strips in the thermal cycler and run the following amplification program:

Step	Repetitions	Temperature	Duration	Acquisition
Taq Activation	1x	95°C	3 min	-
Denaturation	50x	94°C	10 sec	-
Hybridation / Elongation		60°C	20 sec	yes

Enter 20 µL of reaction volume into the Thermal Cycler program. Please refer to the operating instructions of the Thermal Cycler used for programming information.

Setting the detection channels:

Target	Fluorochrome
<i>Plasmodium spp</i> (PAN)	HEX
<i>Plasmodium falciparum</i> (FALCI)	FAM
Internal control (IC)	Cy5

## 9. DATA ANALYSIS AND INTERPRETATION OF RESULTS

### A. Trial Validation Criteria

#### **Negative control:**

The fluorescence emitted must be below the threshold except for the Cy5 channel. This is an indicator of non-specific amplification. If the fluorescence is above the threshold, check for an atypical curve. In the case of an amplification curve, contamination, or a distribution error in the microtubes is to be considered. Only the internal control CIEZ should be amplified.

#### **Positive control:**

The value of the positive control should preferably be detected before 30 cycles ( $Cq \leq 30$ ). In the absence of amplification of the positive control, the existence of an amplification or fluorescence detection problem (defective thermal cycler or instrument non-adapted to the method) should be considered.

#### **Internal amplification control:**

The exogenous internal control CIEZ ensures that the enzymes in the Master mix are functional. Indeed, the amplification curve of the CIEZ internal control should be observed in the Cy5 channel. The internal control value should ideally be detected within 35 cycles ( $Cq \leq 35$ ). Nevertheless, two situations of lack of amplification of the internal control can be observed:

- If the target genes are initially present in the sample with a high number of copies, the internal control may not be amplified. This result is consistent and does not invalidate the test. It should be interpreted as a positive result despite the lack of signal from the internal control. This phenomenon is the result of amplification competition between the Internal Control and targets present at high copy numbers.
- If the target genes in the HEX and FAM channels are not amplified, as well as the internal control in the Cy5 channel, then no result can be rendered. This situation highlights the presence of PCR inhibitors. PCR should be repeated starting from the primary sample and preferably on DNA extract.

### B. Qualitative interpretation (positive or negative)

Signals above the threshold, **and visually consistent with a classical PCR amplification curve**, are considered positive results. Some samples may show atypical curves that are not characteristic of amplification curves. In this case, the result should not be considered as interpretable, and the analysis of the sample must be repeated with the controls.

Detection channels			Interpretation
HEX (PAN)	FAM (FALCI)	Cy5 (Internal Control)	
+	+	+	Patient with <i>Plasmodium falciparum</i> specific DNA and may also have co-infection status
+	+	-	
+	-	+	Patient with <i>Plasmodium spp</i> other than <i>P. falciparum</i> DNA
+	-	-	
-	-	+	Patient without specific DNA of any of the 5 detectable <i>Plasmodium spp</i> species
-	-	-	Invalid result: poor sample quality, inhibition of PCR reaction, or problems during run or test process - perform a new test or re-sampling

## 10. PERFORMANCES

- Analytical sensitivity**

The detection limit (LoD) of the BIOSYNEX AMPLIQUICK Malaria kit is defined as the concentration that can be detected at 95% at least on a specific *Plasmodium* DNA sample.

### Detection limits on extracted DNA samples:

The LoD<sub>95</sub> on the WHO Standard extracted DNA is expressed in IU/μL. This quantified DNA from NIBSC (National Institute for Biological Standards and Control) is the first WHO International Standard for *Plasmodium falciparum* DNA Nucleic Acid Amplification Techniques (NIBSC code: 04/176).

Quantified DNA of the other *Plasmodium* species (*P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi*) are not available on the market, therefore the LoD<sub>95</sub> could not be assessed for those species.

The LoD<sub>95</sub> of the BIOSYNEX AMPLIQUICK Malaria on extracted DNA have been statistically determined for each targeted gene sequence. Results are given in the table below.

Target	LoD <sub>95</sub>
FALCI (FAM)	1.55466801 IU/μL (1554.66801 IU/mL)
PAN (HEX)	0.48978936 IU/μL (489.78936 IU/mL)

### Detection limits for samples treated with the BIOSYNEX Blood Lysis kit:

The LoDs were determined by performing a series of dilutions of the WHO Standard with a known concentration IU/μL treated with the BIOSYNEX Blood Lysis kit.

Quantified DNA of the other *Plasmodium* species (*P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi*) are not available on the market, therefore the LoD<sub>95</sub> could not be assessed for those species.

The LoD<sub>95</sub> of the BIOSYNEX AMPLIQUICK Malaria on DNA treated with BIOSYNEX Blood Lysis has been statistically determined for each targeted gene sequence. Results are given in the table below.

Target	LoD <sub>95</sub>
FALCI (FAM)	71.744765 IU/μL (71 744.765 IU/mL)
PAN (HEX)	26.028150 IU/μL (26 028.150 IU/mL)

- Analytical specificity**

The design of oligonucleotides (primers and probes) was validated *in silico* by BLAST alignment. The comparison of the sequences obtained shows a specific detection of all species of *Plasmodium spp.* and *Plasmodium falciparum*. No sequence homology was observed with human genomic DNA or other pathogens.

The different species of *Plasmodium* have been tested to ensure the specificity of the BIOSYNEX AMPLIQUICK Malaria kit.

A panel of negative whole blood samples was tested. No signals were detected in the FAM and HEX channels. Internal control was valid in all tests.

The BIOSYNEX AMPLIQUICK Malaria kit showed no positive signals with all the following pathogens:

Tested species	DNA detection		
	FAM - FALCI	HEX - PAN	Cy5 -Internal Control
<i>Babesia microti</i>	-	-	+
<i>Babesia divergens</i>	-	-	+
<i>Toxoplasma gondii</i>	-	-	+
<i>Leishmania sp</i>	-	-	+
<i>Hepatitis C virus</i>	-	-	+
<i>Hepatitis B virus</i>	-	-	+
<i>HIV 1</i>	-	-	+
<i>HIV 2</i>	-	-	+
<i>Dengue virus</i>	-	-	+
<i>Borrelia sp</i>	-	-	+

A panel of 79 DNA samples from a biobank was tested with the BIOSYNEX AMPLIQUICK Malaria kit. For all these samples, no amplification of the target DNA sequences was observed.

DNA		
<i>Acanthamoeba castellanii</i>	<i>Enterococcus faecium (VanA)</i>	<i>Mycobacterium ulcerans</i>
<i>Adenovirus</i>	<i>Epstein-Barr Virus</i>	<i>Mycoplasma genitalium</i>
<i>Adenovirus 41</i>	<i>Escherichia coli (EAEC)</i>	<i>Mycoplasma hominis</i>
<i>Aspergillus fumigatus</i>	<i>Escherichia coli (EIEC)</i>	<i>Neisseria gonorrhoeae</i>
<i>Bacillus cereus</i>	<i>Escherichia coli (ETEC)</i>	<i>Neisseria meningitidis Sg A</i>
<i>Bartonella henselae</i>	<i>Escherichia coli (VTEC)</i>	<i>Neisseria meningitidis Sg B</i>
<i>Bartonella Quintana</i>	<i>Francisella tularensis</i>	<i>Neisseria meningitidis Sg C</i>
<i>Bk Virus</i>	<i>Gardnerella vaginalis</i>	<i>Papillomavirus type 16</i>
<i>Bordetella holmesii</i>	<i>Giardia intestinalis</i>	<i>Papillomavirus type 18</i>
<i>Bordetella parapertussis</i>	<i>Haemophilus ducreyi</i>	<i>Parvovirus B19 (Plasmid)</i>
<i>Bordetella pertussis</i>	<i>Haemophilus influenzae</i>	<i>Rickettsia conorii</i>
<i>Borrelia afzelii</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>
<i>Borrelia garinii</i>	<i>Herpes simplex 1</i>	<i>Salmonella typhi</i>
<i>Borrelia burgdorferi</i>	<i>Herpes simplex 2</i>	<i>Staphylococcus aureus (MecA-)</i>
<i>Brucella abortus</i>	<i>Hhv-6</i>	<i>Staphylococcus aureus (MecA+)</i>
<i>Campylobacter jejuni</i>	<i>Hhv-8</i>	<i>Streptococcus agalactiae</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae (NDM-1)</i>	<i>Toxoplasma gondii</i>
<i>Candida auris</i>	<i>Legionella pneumophila</i>	<i>Treponema pallidum</i>
<i>Chlamydia trachomatis</i>	<i>Leishmania chagasi</i>	<i>Trichomonas vaginalis</i>
<i>Chlamydophila pneumoniae</i>	<i>Leishmania infantum</i>	<i>Trypanosoma cruzi</i>
<i>Chlamydophila psittaci</i>	<i>Listeria monocytogenes</i>	<i>Ureaplasma parvum</i>
<i>Clostridium difficile</i>	<i>Moraxella catarrhalis</i>	<i>Ureaplasma urealyticum</i>
<i>Coccidioides immitis</i>	<i>Mycobacterium avium</i>	<i>Varicella-Zoster Virus</i>
<i>Coxiella burnetii</i>	<i>Mycobacterium intracellulare</i>	<i>Vibrio cholerae</i>
<i>Cryptosporidium parvum</i>	<i>Mycobacterium kansasii</i>	<i>Yersinia enterocolitica</i>
<i>Cytomegalovirus</i>	<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pneumoniae</i>
<i>Enterococcus faecalis (VanB)</i>		

#### • Interference study

Whole Blood samples collected in Heparin or EDTA tubes do not interfere with the PCR reaction when the BIOSYNEX AMPLIQUICK Malaria kit is used with the BIOSYNEX Blood Lysis kit.

Whole blood interference tests with abnormal haematological or biochemical parameters on the BIOSYNEX AMPLIQUICK Malaria kit were performed using negative blood samples spiked or not with quantified *P. falciparum* genomic DNA from the ATCC. No positive or negative interference was demonstrated for the following haematological and biochemical parameters:

- Bilirubin  $\leq 158,6$  mg/L
- Triglycerides  $\leq 8,29$  g/L
- Leukocytes  $\leq 13\ 880$  cells/m<sup>3</sup>
- Erythrocytes  $\leq 6\ 300\ 000$  cells/m<sup>3</sup>

To this date, our clinical studies on samples with blood have not shown any interference on the results.

#### • Precision

The precision data in the context of the BIOSYNEX AMPLIQUICK Malaria kit was determined based on following conditions:

- Intra-assay variability (within one test run)
- Inter-assays variability (between different test runs)

- Inter-operators variability (between different operators)
- Inter-laboratory variability (between different locations)
- Inter-lots variability (between different batches)

The variability data are expressed in terms of mean value, standard deviation and coefficient of variation, based on the threshold cycle of quantification (Cq) values of the *Plasmodium Falciparum* DNA.

#### Intra-assay variability data:

Three dilutions of the WHO International Standard for *Plasmodium falciparum* or from a commercial genomic DNA for *P. knowlesi* or from an extracted DNA from a positive whole blood samples for *P. vivax*, *P. ovale* and *P. malariae* were tested: one high (+++), one medium (++) and one low (+), as well as a negative sample (-).

Each sample was tested 30 times.

#### Extracted DNA

Cv= Coefficient of variation

Negative Sample									
Samples	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample -	N/A	-	-	N/A	-	-	27.1	0.2	0.8
Plasmodium falciparum									
Samples	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	23.0	0.1	0.5	22.9	0.1	0.3	28.9	0.7	2.5
Sample ++	30.0	0.1	0.3	29.8	0.1	0.2	27.0	0.2	0.7
Sample + (FALCI)	37.1	0.7	1.9	36.8	0.4	1.0	27.1	0.2	0.7
Sample + (PAN)	Below LoD			38.6	0.6	1.5	27.0	0.2	0.7
Plasmodium malariae									
Samples	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	N/A	-	-	26.7	0.1	0.3	27.2	0.4	1.5
Sample ++	N/A	-	-	31.1	0.1	0.3	27.1	0.2	0.7
Sample +	N/A	-	-	38.2	0.7	1.7	27.2	0.2	0.6
Plasmodium ovale									
Samples	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	N/A	-	-	24.2	0.1	0.3	27.9	0.9	3.2
Sample ++	N/A	-	-	31.0	0.1	0.3	26.8	0.2	0.9
Sample +	N/A	-	-	37.9	0.4	1.2	26.9	0.2	0.7
Plasmodium vivax									
Samples	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	N/A	-	-	28.9	0.1	0.3	27.2	0.4	1.3
Sample ++	N/A	-	-	32.4	0.1	0.3	27.4	0.2	0.8
Sample +	N/A	-	-	39.3	0.7	0.7	27.3	0.3	1.3

Samples	<i>Plasmodium knowlesi</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	N/A	-	-	28.5	0.2	0.6	26.7	0.1	0.4
Sample ++	N/A	-	-	31.5	0.1	0.3	27.3	0.1	0.5
Sample +	N/A	-	-	39.5	0.5	1.3	27.4	0.1	0.5

Pre-treatment with BIOSYNEX Blood Lysis kit

Samples	<i>Negative Sample</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample -	N/A	-	-	N/A	-	-	27.2	0.2	0.7

Samples	<i>Plasmodium falciparum</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	26.4	0.1	0.4	26.1	0.0	0.2	26.9	0.2	0.9
Sample ++	31.0	0.1	0.3	30.7	0.1	0.3	27.0	0.2	0.7
Sample + (FALCI)	37.6	0.6	1.5	37.3	0.4	1.0	26.9	0.2	0.6
Sample + (PAN)	<i>Below LoD</i>			38.9	0.6	1.6	27.0	0.2	0.7

Samples	<i>Plasmodium malariae</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample ++	N/A	-	-	32.1	0.1	0.2	27.5	0.4	1.4
Sample +	N/A	-	-	37.4	0.6	1.7	27.4	0.1	0.5

Samples	<i>Plasmodium ovale</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample +++	N/A	-	-	28.5	0.1	0.2	27.2	0.2	0.7
Sample ++	N/A	-	-	31.1	0.1	0.2	27.0	0.2	0.6
Sample +	N/A	-	-	37.8	0.5	1.4	27.2	0.2	0.7

Samples	<i>Plasmodium vivax</i>								
	FAM - FALCI			HEX - PAN			Internal control CIEZ		
	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %	Average Cq value	Standard deviation	CV %
Sample ++	N/A	-	-	30.7	0.1	0.3	27.4	0.5	1.9
Sample +	N/A	-	-	37.6	0.3	0.9	27.3	0.2	0.7

**Inter-assays variability data:**

Two sample dilutions of the *Plasmodium falciparum* WHO Standard were tested: one high (+++) corresponding to 2 500 copies/μL and one low (+) corresponding to 3x LoD<sub>95</sub>. Each sample is tested in duplicate each day, for 5 days. Negative sample give negative results.

Targets	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
FAM - FALCI	27.3	0.0	0.2	38.8	0.6	1.5
HEX - PAN	26.4	0.1	0.5	37.9	0.3	0.8
Cy5 – Internal Control	26.7	0.1	0.4	26.8	0.1	0.5

#### Inter-operators variability data:

Two sample dilutions of the *Plasmodium falciparum* WHO Standard were tested: one high (+++) corresponding to 2 500 copies/μL and one low (+) corresponding to 3x LoD<sub>95</sub>.

Each sample is tested in duplicate, by two different operators.

Negative sample give negative results.

Targets	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
FAM - FALCI	27.3	0.1	0.4	39.2	0.2	0.4
HEX - PAN	26.5	0.0	0.1	38.0	0.5	1.4
Cy5 – Internal Control	26.8	0.1	0.3	26.9	0.1	0.5

#### Inter-laboratory variability data:

Two sample dilutions of the *Plasmodium falciparum* WHO Standard were tested: one high (+++) corresponding to 2 500 copies/μL and one low (+) corresponding to 3x LoD<sub>95</sub>.

Each sample is tested in duplicate, by the same operator in two different laboratories.

Negative sample give negative results.

Targets	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
FAM - FALCI	27.4	0.2	0.7	39.0	0.0	0.1
HEX - PAN	26.4	0.1	0.2	37.8	0.2	0.5
Cy5 – Internal Control	26.8	0.0	0.0	26.7	0.1	0.5

#### Inter-lots variability data:

Two sample dilutions of the *Plasmodium falciparum* WHO Standard were tested: 1 high (+++) corresponding to 2 500 copies/μL and 1 low (+) corresponding to 3x LoD<sub>95</sub>.

Each sample is tested in duplicate, on three different lots.

Negative sample give negative results.

Targets	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
FAM - FALCI	27.4	0.3	1.1	38.8	0.3	0.8
HEX - PAN	26.5	0.1	0.3	38.2	0.2	0.5
Cy5 – Internal Control	27.0	0.2	0.9	27.0	0.2	0.8

- Clinical performances**

The clinical performances were determined on 416 samples for the DNA extraction/purification with a commercial kit and on 149 samples for samples treated with BIOSYNEX Blood Lysis kit.

Total samples from all studies for DNA extraction/purification with a commercial kit:

	Number of Samples	Malaria Qualification					
		Negative Samples	Positive Samples				
			<i>P. falciparum</i>	<i>P. vivax</i>	<i>P. ovale</i>	<i>P. malariae</i>	Co-Infections
Study 1	67	0	31	13	15	10	2
Study 2	93	0	93	0	0	0	0
Study 3	100	100	0	0	0	0	0
Study 4	48	10	34	0	2	2	0
Study 5	108	108	0	0	0	0	0
<b>Total</b>	<b>416</b>	<b>218</b>	<b>158</b>	<b>13</b>	<b>17</b>	<b>12</b>	<b>2</b>

Total samples from all studies for samples treated with BIOSYNEX Blood Lysis kit:

	Number of Samples	Malaria qualification					
		Negative Samples	Positive Samples				
			<i>P. falciparum</i>	<i>P. vivax</i>	<i>P. ovale</i>	<i>P. malariae</i>	Co-Infections
Study 4	41	10	27	0	2	2	0
Study 5	108	108	0	0	0	0	0
<b>Total</b>	<b>149</b>	<b>118</b>	<b>27</b>	<b>0</b>	<b>2</b>	<b>2</b>	<b>0</b>

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Malaria kit on samples after DNA extraction/purification with a commercial kit:

		Sample qualification <i>Plasmodium spp</i>	
		Positive	Negative
BIOSYNEX AMPLIQUICK Malaria	Positive	198	0
	Negative	0	218

**Sensitivity: 100.00%**      **Specificity: 100.00%**  
 (95%CI: 98.15% to 100.00%)      (95%CI: 98.32% to 100.00%)

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Malaria on samples treated with BIOSYNEX Blood Lysis kit:

		Sample qualification <i>Plasmodium spp</i>	
		Positive	Negative
BIOSYNEX AMPLIQUICK Malaria	Positive	31	0
	Negative	0	118

**Sensitivity: 100.00%**      **Specificity: 100.00%**  
 (95%CI: 88.78% to 100.00%)      (95%CI: 96.92% to 100.00%)

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Malaria kit on *Plasmodium falciparum* specific samples after DNA extraction/purification with a commercial kit:

		Sample qualification <i>Plasmodium falciparum</i>	
		Positive	Negative
BIOSYNEX AMPLIQUICK Malaria	Positive	158	0
	Negative	0	260

**Sensitivity: 100.00%** (95%CI: 97.69% to 100.00%)      **Specificity: 100.00%** (95%CI: 98.59% to 100.00%)

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Malaria on *Plasmodium falciparum* specific samples treated with BIOSYNEX Blood Lysis kit:

		Sample qualification <i>Plasmodium falciparum</i>	
		Positive	Negative
BIOSYNEX AMPLIQUICK Malaria	Positive	27	0
	Negative	0	122

**Sensitivity: 100.00%** (95%CI: 87.23% to 100.00%)      **Specificity: 100.00%** (95%CI: 97.02% to 100.00%)

For this study, we did not have access to clinical samples from patients positive for *P. knowlesi* (rare samples). The performances of the kit for that species were validated using synthetic DNA.





## 11. LIMITATIONS

1. Failure to follow any instruction of the IFU may adversely affect the test performance and/or invalidate the test result.
2. The possibility of obtaining false positive or false negative results cannot be ruled out (e.g. in case of accidental contamination, in case of poor quality of the sample or in case the test procedure is not correctly followed). As with all diagnostic tests, results must be considered with other clinical information available to the physician.

## 12. LITERATURE

1. World Health Organization (WHO) Malaria Fact Sheets 29 March 2023. <https://www.who.int/news-room/fact-sheets/detail/malaria>
2. ECDC Surveillance Report: Malaria, Annual Epidemiology Report for 2021. <https://www.ecdc.europa.eu/en/publications-data/malaria-annual-epidemiological-report-2021>

## 13. INTERNATIONAL SYMBOLS

	Consult instructions for use or electronic instructions for use		Contains sufficient for <n> tests
	Catalogue number		<i>In vitro</i> diagnostic medical device
	Positive control		Negative control
	Temperature limits		Use-by-date
	Manufacturer		Batch code
	Authorized representative in Switzerland		Do not use if package is damaged and consult the instructions for use
	Master Mix		Bag of strips of caps
	Keep away from sunlight		Importer
	Do not reuse		Microplate
	Unique Device Identifier		Qualitative detection by PCR

## 14. MANUFACTURER INFORMATION

### BIOSYNEX S.A.

22 boulevard Sébastien Brant  
67400 ILLKIRCH-GRAFFENSTADEN  
France

#### Standard:

Phone: +33 3 88 78 78 87

[www.biosynex.com](http://www.biosynex.com)

#### Contacts France:

Phone: +33 3 88 77 57 00

[service.clients@biosynex.com](mailto:service.clients@biosynex.com)

#### Contacts other countries:

Phone: +33 3 88 77 57 52

[sales@biosynex.com](mailto:sales@biosynex.com)

#### Customer service

Phone : +33 3 88 77 57 25

[tech.support@biosynex.com](mailto:tech.support@biosynex.com)

### BIOSYNEX SWISS S.A.

Route de Rossemaison 100  
2800 DELEMONT - Switzerland

## 15. HISTORY

Record version	Amended paragraph	Change details
V1	Entire IFU	Compliance with IVDR 2017/746 and multi-language format;
V2	Contents	Addition of the summary.
	6: Kit storage	Change « Store the kit at a temperature of -20°C or below.»
	7: Sample collection and storage	Change « Samples can be stored for 3 years at -20°C. »
	8: Test procedure	Removal « MagPurix Blood DNA Extraction kit - Reference ZP02001 (Zinexts LifeScience Corp)»
	10: Performances	Addition of results on blood interference tests with abnormal haematological or biochemical parameters.
	11: Limitations	Rewording of limitation No.2.
	15 : Modifications history	Creation of this paragraph.
V3	10 : Performances	Inter-lots variability data: Review of values