# BIO SYNEX

## Ampliquick Fecal Virology

**REF** 3150047\_SEC01 / 3150047\_SEC02 / 3150047\_BLK96

FR (FRANÇAIS / FRENCH) .....	2
EN (ENGLISH) .....	16

### CONTENTS

## BIOSYNEX AMPLIQUICK Fecal Virology

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## BIOSYNEX AMPLIQUICK Fecal Virology

QUALITATIVE DETECTION OF GASTROINTESTINAL VIRUSES FROM STOOL SAMPLES.

*In vitro* diagnostic medical device for professional use.

**REF** 3150047\_SEC01 / 3150047\_SEC02 / 3150047\_BLK96

### 1. INTENDED PURPOSE

BIOSYNEX AMPLIQUICK Fecal Virology is a molecular *in vitro* diagnostic test for the qualitative detection by RT-qPCR of the presence of gastrointestinal viruses: Rotavirus, Sapovirus, Astrovirus, Norovirus GI/GII and Adenovirus F40/41 in human stool samples. The test is an aid in the diagnosis of viral gastroenteritis performed using a purified DNA/RNA extract obtained from stool samples collected with BIOSYNEX AMPLIQUICK Fecal Pretreatment or with FecalSwab™ (Copan). For a syndromic diagnosis, the PCR can be performed simultaneously with the BIOSYNEX AMPLIQUICK Fecal Bacteriology test (ref. 3150046 and related references) using the same samples and the same amplification program. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

### 2. CLINICAL SUMMARY

Gastroenteritis is the inflammation of the gastrointestinal tract that causes diarrhea, abdominal pain, nausea, fever and sometimes vomiting. This inflammation affects nearly 1.7 billion children every year and is a leading cause of death in children under five years old mainly due to severe dehydration<sup>1</sup>. It can be provoked by many different infectious agents such as parasites (10-15%), bacteria (15-20%) and viruses (50-70%) that are spread by the fecal-oral route through the contamination of food and water.

Viruses commonly associated with gastroenteritis are Rotavirus, Sapovirus, Astrovirus, Norovirus and Adenovirus.

Rotaviruses are the most common cause of diarrhea in children worldwide<sup>2</sup>. These non-enveloped viruses consist of three concentric capsids that encompass a double-stranded RNA genome. There are 10 different species of Rotavirus (species A-J) with species A being the most prevalent in gastrointestinal infections.

Sapoviruses are associated with less than 5% of gastroenteritis cases worldwide. They are non-enveloped, small (30-38 nm), single-stranded RNA viruses classified in 5 different genogroups (GI-GV).

Astroviruses are small non-enveloped single-stranded RNA viruses classified into 8 different genotypes. The infection with Astroviruses is very common as an important proportion of the population carries antibodies against type1 Astrovirus.

Noroviruses are single-stranded RNA viruses with relatively small genome. These viruses are divided into 5 different genogroups (GI-GV): the GI, GII and GIV are known to infect humans. The GI and GII are the most associated genogroups with gastroenteritis. Around 685 million cases of Noroviruses are detected every year.

Adenoviruses are the third most common cause of gastrointestinal infections in children (9% of children with gastroenteritis). They are medium-sized 70-90 nm, non-enveloped viruses consisting of icosahedral nucleocapsid containing a double-stranded DNA genome (26-45 Kbp). There are 88 different Human Adenovirus types that are grouped into 7 species (A to G). F40 and F41 are the most common serotypes of Adenoviruses found in feces of patients.<sup>3,4</sup>

### 3. TEST PRINCIPLE

The AMPLIQUICK Fecal Virology test uses a single PCR (Polymerase Chain Reaction) step to simultaneously amplify and detect a specific region of the viral genome of Rotavirus, Sapovirus, Astrovirus, Norovirus, Adenovirus, an exogenous CIEZ internal control and an exogenous Rp49 control in process. The test is composed of 2 different master mixes: Mix A which targets the detection of Rotavirus, Sapovirus and Astrovirus and Mix B which targets the detection of Norovirus GI/GII and Adenovirus F40/41.


The test is based on real-time reverse transcription Polymerase Chain Reaction (PCR) technology by fluorescent probe hydrolysis. This method uses probes labelled with the fluorophores FAM, HEX and Texas Red for the

amplification of viral sequences and probes labelled with the fluorophore Cy5 for the detection of the sequences used as controls CIEZ and Rp49. It should be noted that detection of the CIEZ internal control validates the PCR amplification step and rules out a possible false negative result due to PCR inhibition. Rp49 is a known DNA sequence from *Drosophila* that is absent in the human and viral genomes and is therefore added to stool samples as an extraction control, also known as a control in process. The supplementation of the Rp49 exogenous sequence in the fecal samples and its detection allow to control the quality of the biological sample and to validate the genomic material purification step.

The increase in fluorescence signal is only detected if the target sequence is complementary to the amplified probe present in the sample. The fluorescent signal is therefore directly proportional to the amplification of the target during the amplification phase. The Cq value (quantification cycle) corresponds to the number of cycles during which the fluorescence increases exponentially above the background noise.

BIOSYNEX AMPLIQUICK Fecal Virology is available in two formats:

- In bulk tubes (ref. 3150047\_BLK96): each master mix is provided in a microtube for the user to dispense in PCR plates.
- In prefilled microplates (ref. 3150047\_SEC01 and 3150047\_SEC02): the master mixes are distributed in the wells of a PCR plate according to order shown below.

Regardless the packaging, the Master mixes are ready-to-use and contains dNTPs, MgCl<sub>2</sub>, primers and fluorescent probes, Taq polymerase enzyme, reverse transcriptase and reaction buffer.

BIOSYNEX AMPLIQUICK Fecal Virology must be used with stool sample collected with FecalSwab™ (Copan) or collected and treated with BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065). Regardless the collection system used, the nucleic acids must be extracted and purified before performing the amplification step.

#### 4. REQUIRED MATERIAL

3150047_SEC01	3150047_SEC02	3150047_BLK96
<b>Kit composition</b>		
<ul style="list-style-type: none"> <li>- 2 dividable microplates with 96 wells format for use with 0.1 mL low-profile (SEC01)/ 0.2 mL (SEC02) high-profile thermal cycler</li> <li>- 1 tube of positive control (CONTROL+, red cap)</li> <li>- 1 tube of negative control (CONTROL-, green cap)</li> <li>- 2 tubes of control in process (CONTROL IP, yellow cap)</li> <li>- 1 tube of control master mix (CONTROL Mmix, blue cap)</li> <li>- 2 bags with 12x8 wells strips of caps</li> </ul>		<ul style="list-style-type: none"> <li>- 4 tubes of Master MIX A (Mmix A, white cap marked V1)</li> <li>- 4 tubes of Master MIX B (Mmix B, white cap marked V2)</li> <li>- 1 tube of positive control (CONTROL+, red cap)</li> <li>- 1 tube of negative control (CONTROL-, green cap)</li> <li>- 2 tubes of control in process (CONTROL IP, yellow cap)</li> <li>- 1 tube of control Master Mix (CONTROL Mmix, blue cap)</li> </ul>

Material required but not provided	
<ul style="list-style-type: none"> <li>- BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or FecalSwab® (Copan) transport medium (ref. 470CE)</li> <li>- DNA/RNA extraction kit</li> <li>- Powder-free disposable gloves</li> <li>- Pipettes and filtered tips</li> <li>- qPCR Thermal Cycler</li> <li>- PCR microplate or microtubes centrifuge</li> </ul>	<ul style="list-style-type: none"> <li>- BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or FecalSwab® (Copan) transport medium (ref. 470CE)</li> <li>- DNA/RNA extraction kit</li> <li>- Powder-free disposable gloves</li> <li>- Pipettes and filtered tips</li> <li>- qPCR Thermal Cycler</li> <li>- PCR plates or tube strips + films or optical sealing caps</li> <li>- PCR microplate or microtubes centrifuge</li> </ul>

The PCR thermal cycler used for the test must have the following main characteristics:

- Open system
- Real-time quantitative PCR assays.
- Programmable Thermal Cycler block

Reference	Thermal cycling block
3150047_SEC01	0.1mL low profile
3150047_SEC02	0.2mL high profile

- Excitation source: Leds, lamp or laser
- Filter sets (excitation/emission wavelengths) suitable for the detection of "reporter" fluorophores of the FAM, HEX, TexasRed and Cy5 probes.
- Connection with a computer using specific analysis software that allow the recovery of fluorescence data, absolute quantification assays, and the interpretation of results.

The kit has been developed and validated with the following real-time PCR thermal cyclers:

- CFX96 Touch™ Real-Time PCR Detection System (BioRad)
- CFX96 OPUS™ Real-Time PCR Detection System (BioRad)
- QuantGene 9600 (Bioer)
- QuantStudio™ 5 System (Applied Biosystems™)

If another thermal cycler is used, please perform your own validation of the BIOSYNEX AMPLIQUICK Fecal Virology kit before using the test. We recommend using qualified samples and the negative and positive controls provided in the kit.

## 5. PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of many samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination. Pre-handling PCR, post-PCR and DNA/RNA extraction should therefore be performed in separate rooms. Wear disposable gloves in each room and change them before moving between different rooms.
- The provided strips of caps are for single use only. Do not reuse them. Do not open the PCR tubes or microplates at the end of the test.
- For the prefilled kit formats \_SEC01 and \_SEC02: Do not use the test if the aluminum foil is opened or damaged. (Once the aluminum foil is removed, use the test immediately.
- Do not use the kit if it is received defrosted.
- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.
- Protect the kit from light.
- Centrifuge the PCR strips before opening; open each PCR strip individually and rapidly close it after usage. Do not open several PCR strips simultaneously to avoid any cross-contamination
- The positive control (CONTROL +) contains significant amounts of DNA sequences. It can therefore potentially

contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store the CONTROL + separately from the other components of the kit, preferably outside the kit on the first opening.

- The negative and positive controls included in the kit mimic results obtained with negative or positive samples respectively. They must be used with each run.
- To ensure that there is no contamination during handling, it is left to the user to include as part of their internal quality approach, at the frequency chosen, an experiment negative control well in which 5µL of molecular biology grade water are added to the Master mix.
- When using negative and positive controls with a series of patients, it is recommended to first deposit the negative control, then deposit the patient samples, and to finish with the deposition of the positive control.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste.
- The device contains material of bacterial or animal origin and should be handled with caution.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Fecal Virology, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed.

## 6. KIT STORAGE

The kit is shipped frozen. Kit components must arrive frozen. Store the kit at a temperature below or equal to -20°C. Under these conditions, the reagents are stable until the expiry date indicated on the kit label. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light.

The Control Master mix can undergo 16 freeze-thaw cycles and the Master mixes in bulk format (\_BLK96) can undergo 10 freeze-thaw cycles without the shelf-life and qualities of the product being altered. The controls (including the control in process) can undergo 30 freeze-thaw cycles without the shelf-life and qualities of the product being altered.

*NB:* Only thaw the number of well strips or plates needed for the pre-filled kit format. As the plates and strips are ready to use, there is no need for repeated thaw/freeze cycles.

## 7. SAMPLE COLLECTION AND STORAGE

The stool sample should be collected with BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or with FecalSwab™ (Copan) from feces collected in a preservative-free container.

For sample collection and storage, please refer to the instructions for use of BIOSYNEX AMPLIQUICK Fecal Pretreatment or FecalSwab™ (Copan).

For samples in FecalSwab™: the Control in process can be added ahead of time. The fecal swab samples discharged in the transport medium and supplemented with the AMPLIQUICK Fecal Virology control in process can be stored up to 72 hours at room temperature or at 4°C and up to 1 week between -20°C and -80°C before proceeding to the nucleic acid extraction and testing with the BIOSYNEX AMPLIQUICK Fecal Virology PCR detection kit.

## 8. TEST PROCEDURE

### 1. Sample preparation:

Add 5µL of the control in process (CONTROL IP, yellow cap) in the sample preparation (i.e. Fecal Pretreatment tube or FecalSwab™ with the stool sample).

### 2. Extraction of nucleic acids:

The nucleic acids extraction must be performed before the amplification pcol using a suitable extraction system for viral DNA/RNA. Please follow the manufacturer's instructions for use.

The positive and negative controls provided in the kit should not undergo extraction and should be used as provided.

The BIOSYNEX AMPLIQUICK Fecal Virology kit has been validated with the following extraction kits:

- Macherey-Nagel NucleoMag® Dx Pathogen ref 744215 using the Magnetapure automation system provided by Dutscher
- QIAGEN QIAampFAST DNA Stool Mini Kit (ref. 51604)

- Singuway Nucleic Acid Extraction Reagents (ref. MTQM036) with the Singu20 Nucleic Acid Extractor Automated System (Singuway).


If another DNA/RNA extraction method or kit is used, please perform your own validation method of the BIOSYNEX AMPLIQUICK Fecal Virology kit using positive and negative qualified samples supplemented with the BIOSYNEX AMPLIQUICK Fecal Virology control in process.

- Amplification protocol:

It is recommended that the Master mixes be placed on a refrigerated rack or on ice while samples are being added.

**For microplates prefilled with ready-to-use master mixes (REF 3150047\_SEC01 & 3150047\_SEC02):**


1. Take from the freezer the number of strips or plates needed. One-8wells strip = 4 samples. Include wells for one negative control and one positive control per run.
2. Centrifuge for 10 seconds to remove any droplets from the walls of the wells or the aluminum sealing foil.
3. Remove the aluminum sealing foil.
4. Add 5  $\mu$ L of negative control to the intended wells.
5. Add 5  $\mu$ L of extraction eluate from the same sample into the wells containing the master mixes (A and B as shown on below picture).
6. Finally, add 5  $\mu$ L of positive control to the intended wells.
7. Close the plate or strips using the provided caps and centrifuge once more for 10 seconds to eliminate potential bubbles at the surface.

**With tubes of ready-to-use master mixes (REF 3150047\_BLK96):**

1. Remove the Master mixes tubes from the freezer and allow them to thaw on ice or refrigerated rack.
2. Centrifuge the tubes for 10 seconds to remove any droplets from the sides or caps.
3. Dispense 15  $\mu$ L of each Master mix into PCR tubes or plate. Each Master Mix must be deposited in a separate well. So, for each patient, 2 wells enable a diagnosis over the entire panel.
4. Dispense 15  $\mu$ L of control master mix in two empty wells for the controls.
5. Add 5  $\mu$ L of negative control to one of the tubes/wells containing the control master mix.
6. Add 5  $\mu$ L of extraction eluate from the same sample to the dedicated reaction wells (Master mixes A and B).
7. Finally, add 5  $\mu$ L of positive control in the last tube/well containing the control master mix.
8. Close the wells with optical caps or film.
9. Centrifuge for 10 seconds.
10. Return the master mixes tubes to the freezer.

*N.B.* To limit the loss of reagent, it is important to pipette the Master Mix without plunging the tip completely into the volume of Master Mix. Alternatively, you could use the “low-binding” or “low retention” tips to limit reagent loss.

Place the plate, strips or tubes in the thermal cycler and run the following amplification program:

Step	Repetitions	Temperature	Duration	Acquisition
Reverse transcription	1x	45°C	15 min	-
Taq Activation	1x	95°C	2 min	-
Denaturation	50x	95°C	5 sec	-
Hybridation / Elongation		60°C	20 sec	yes

If you use :

- CFX96 Touch™ Real-Time PCR Detection System (BioRad)
- CFX96 OPUS™ Real-Time PCR Detection System (BioRad)

To analyse, please, select **Apply fluorescent drift correction**.

Enter 20 µL as the reaction volume into the Thermal Cycler program.

Please refer to the operating instructions of the Thermal Cycler used for programming information.

Detection channels settings :

Master mixes	Channel	Targeted Virus
<b>A</b>	FAM	<i>Rotavirus</i>
	HEX	<i>Sapovirus</i>
	TxRed	<i>Astrovirus</i>
	Cy5	Internal amplification control
<b>B</b>	FAM	<i>Norovirus GI/GII</i>
	HEX	<i>Adenovirus F40/41</i>
	Cy5	Control in process

## 9. INTERPRETATION OF RESULTS

The threshold of the real-time PCR reaction corresponds to the level of signal that reflects a statistically significant increase in fluorescence emission compared to the baseline signal calculated; the value is set to distinguish a significant amplification signal from background noise. We recommend automatic setting of the threshold by the real-time PCR instrument software, rather than manual setting. If further settings should be done for the test, please contact the thermal cycler manufacturer's technical support team.

### Samples:

Signals above the threshold, **and visually consistent with a classical PCR amplification curve**, are considered positive results. Please take into consideration the amplification curves of the positive control from the kit as reference.

Some samples may show atypical curves that are not characteristic of amplification curves. In this case, the result should not be considered as interpretable, and the analysis of the sample should be repeated in a new PCR run

that includes positive and negative controls. Cq cut-off values should also be considered when interpreting the results.

Result interpretations and Cq cut-off values are given per target/channel in the table below.

Samples can be positive for one or more pathogens detected in the same well or in different wells.

## Quality control

### Internal control

The internal control CIEZ ensures that the enzymes in the Master mix are functional. Indeed, an amplification curve needs to be observed in Cy5 channel of Master mix A. The internal control value must preferentially be detected before 35 cycles (Cq ≤ 35). Nevertheless, two situations of lack of amplification of the internal control can be observed:

- If the target sequences are initially present in the sample with a high number of copies, the internal control may not be amplified. This result is consistent and does not invalidate the test. It should be interpreted as a positive result despite the lack of signal from the internal control. This phenomenon is the result of amplification competition between the Internal Control and targets that are present at high copy numbers.
- If the target sequences in the FAM, HEX and Texas Red channels are not amplified, as well as the internal control in the Cy5 channel, then no result can be rendered. This situation highlights the presence of PCR inhibitors or a technical problem occurring during test realization. The PCR should be repeated:
  - after diluting the obtained DNA/RNA extract at 1/10
  - or after realizing a new extraction / purification of the nucleic acids.

### Control in process (CONTROL IP, yellow cap)

The exogenous Control in Process Rp49 should be detected in the Cy5 channel of Master mix B. This procedural control validates the quality of the biological sample collection and validates the DNA/RNA extraction and purification steps.

In the absence of an amplification signal for the procedural control, the test is considered invalid when no target has been detected. In this case, the nucleic acids extraction or the sample treatment must be repeated.

### Negative control (CONTROL -, green cap)

The fluorescence emitted must be below the threshold except for the Cy5 channel. This is an indicator of non-specific amplification. If the fluorescence is above the threshold, check for an atypical curve. In the case of an amplification curve, a contamination or a distribution error in the microtubes should be considered.

### Positive control (CONTROL +, red cap)

The value of the positive control should preferably be detected before 30 cycles (Cq ≤ 30) in the four channels (FAM, HEX, TxRed and Cy5). In the absence of amplification of the positive control, the existence of an amplification problem can be suggested. Some samples show atypical curves that are not characteristic of sigmoid amplification curves. In this case, the result should not be considered interpretable, and the sample should be analyzed again.

Detection channels							Interpretation
Master mix A				Master mix B			
FAM (NSP3)	HEX (RdRp-VP1)	TxRed (ORF1b)	Cy5 (Internal control)	FAM (ORF1-ORF2)	HEX (hexon gene)	Cy5 (Control in process)	
+ Cq<41	-	-	+ Cq ≤ 35	-	-	Cq ≤ 35	Rotavirus
-	+ Cq<41	-	+ Cq ≤ 35	-	-	Cq ≤ 35	Sapovirus
-	-	+	+ Cq ≤ 35	-	-	+ Cq ≤ 35	Astrovirus
+ Cq<41	-	-	-	-	-	+ Cq ≤ 35	Rotavirus
-	+ Cq<41	-	-	-	-	+ Cq ≤ 35	Sapovirus
-	-	+	-	-	-	+ Cq ≤ 35	Astrovirus
-	-	-	-	-	-	-	Invalid
-	-	-	+ Cq ≤ 35	-	-	+ Cq ≤ 35	Negative for all targets
-	-	-	+ Cq ≤ 35	+	-	+ Cq ≤ 35	Norovirus GI/GII
-	-	-	+ Cq ≤ 35	+	-	-	Norovirus GI/GII

-	-	-	+ Cq ≤ 35	-	+ Cq<41	+ Cq ≤ 35	Adenovirus F40/41
-	-	-	+ Cq ≤ 35	-	+ Cq<41	-	Adenovirus F40/41

## 10. PERFORMANCES

### • Clinical performance

The clinical performances of the BIOSYNEX AMPLIQUICK Fecal Virology PCR kit were determined in two different studies:

- Study 1: 21 negative and 27 positive fecal swab samples from the University Hospital of Amiens, France
- Study 2: 50 negative and 395 positive fecal swab samples at the French National Reference Center of viral gastroenteritis laboratory

The determined clinical performances from Study 1:

		Sample Qualification*	
		Positive	Negative
BIOSYNEX AMPLIQUICK Fecal Virology	Positive	25**	0
	Negative	0	23

\*Original Sample Qualification by a CE marked commercially available rapid test. As the rapid test has a low specificity and a low sensitivity, the qualification results of a CE marked commercially available PCR diagnostic kit were considered in case of discrepancies.

\*\*12 out of 27 samples were found positive to others or additional viral pathogens detected by the AMPLIQUICK Fecal Virology PCR kit and confirmed by a CE marked commercially available PCR diagnostic kit but not by the CE marked rapid test.

The determined clinical performances from Study 2:

Master Mix	Target	True Positive	True Negative	False Positive	False Negative	Sensitivity (95% CI)	Specificity (95% CI)
A	Rotavirus	61	50	0	0	100% (94.13-100%)	100% (92.89-100%)
	Sapovirus	60	50	0	0	100% (94.04-100%)	100% (92.89-100%)
	Astrovirus	41	50	0	0	100% (91.40-100%)	100% (92.89-100%)
B	Norovirus GI/GII	173	50	0	1	99.45% (96.84-99%)	100% (92.89-100%)
	Adenovirus F40/F41	52	69 (21 non-F Adenovirus*)	0	1	98.25% (89.93-99.95%)	100% (94.79-100%)

\*21 samples were positive to non-F Adenoviruses but were not detected by the AMPLIQUICK Fecal Virology. The BIOSYNEX AMPLIQUICK Fecal Virology is specific to Adenovirus F40 and F41 serotypes responsible for gastroenteritis.

After combining the results from both clinical studies, the claimed clinical performances for the BIOSYNEX AMPLIQUICK Fecal Virology PCR kit are:

Master Mix	Target	True Positive	True Negative	False Positive	False Negative	Sensitivity (95% CI)	Specificity (95% CI)
A	Rotavirus	77	73	0	0	100% (95,32-100%)	100% (95,07-100%)
	Sapovirus	61	50	0	0	100% (94,13-100%)	100% (92,89-100%)
	Astrovirus	41	50	0	0	100% (91,40-100%)	100% (92,89-100%)
B	Norovirus GI/GII	181	73	0	1	99,45% (96,98-99,99%)	100% (95,07-100%)
	Adenovirus F40/F41	56	92	0	1	98,25% (90,61-99,96%)	100% (96,07-100%)

\*21 samples were positive to non-F Adenoviruses but were not detected by the AMPLIQUICK Fecal Virology. The BIOSYNEX AMPLIQUICK Fecal Virology is specific to Adenovirus F40 and F41 serotypes responsible for gastroenteritis.

- Analytical Sensitivity**

The detection limit of the BIOSYNEX AMPLIQUICK Fecal Virology Real-Time PCR Kit is defined as the concentration, in number of copies, that can be 95% detected in a DNA/RNA sample specific to each target.

**Detection limits on DNA/RNA samples per targeted virus:**

Quantified RNA samples were purchased from Vircell (Amplirun®) for Rotavirus and from ATCC for Sapovirus, Astrovirus, Norovirus G1 and Norovirus G2. Quantified Adenovirus F41 DNA sample was purchased from Vircell (Amplirun®). Quantified Adenovirus F40 DNA sample is not available on the market and therefore the sample used was a quantified synthetic DNA purchased from idtDNA (gBlock).

The LoD<sub>95</sub> of the BIOSYNEX AMPLIQUICK Fecal Virology on DNA/RNA samples was statistically determined for each pathogen and targeted nucleic sequence. Results are given in the table below:

Master mix	Target	LoD <sub>95</sub> (DNA/RNA copy/μL)	LoD <sub>95</sub> (DNA/RNA copy/rxn)
A	Rotavirus	4.39	21.95
	Sapovirus	5.07	25.35
	Astrovirus	5.32	26.6
B	Norovirus G1	43.77	218.85
	Norovirus G2 (ATCC strain)	2.34	11.7
	Norovirus G2 (Vircell strain)	8.34	41.7
	Adenovirus F40	1.21	6.05
	Adenovirus F41	0.43	2.15

- Analytical Specificity**

The design of oligonucleotides was validated *in silico* by BLAST alignment. The comparison of the primers and probes sequences shows that the BIOSYNEX AMPLIQUICK Fecal Virology PCR kit is specific to the detection of Rotaviruses, Sapovirus, Astrovirus, Norovirus G1, Norovirus G2, Adenovirus F40 and Adenovirus F41.

**Cross-reactivity**

The specificity of the BIOSYNEX AMPLIQUICK Fecal Virology PCR kit was confirmed by testing a panel of 42 RNA and 95 DNA of different microorganisms. No cross-reactivity was detected between any of the following pathogens:

RNA		
Coronavirus OC43	<i>Chikungunya virus</i>	Parainfluenza 1
Coronavirus	<i>Rubella virus</i>	Parainfluenza 2
Coronavirus SARS (2003)	SARS-CoV-2	Parainfluenza 3
Coxsackie A6	<i>Tick-borne-encephalitis virus</i>	Parainfluenza 4 A
Coxsackie B1	<i>Human parainfluenza 1</i>	Respiratory syncytial virus (subtype A)
Coxsackie B5	<i>Influenza A H1</i>	Respiratory syncytial virus (subtype B)
Dengue 1 virus	<i>Influenza A H3</i>	West Nile virus
Dengue 2 virus	<i>Influenza A H5</i>	Yellow fever virus
Dengue 3 virus	<i>Influenza B</i>	Zika virus (Asian lineage)
Dengue 4 virus	<i>Measles</i>	Zika virus
Echovirus 5	<i>MERS-CoV</i>	-
Enterovirus 68	<i>Mumps</i>	-
Rhinovirus	<i>Novel influenza A H1N1</i>	-

DNA		
<i>Adenovirus c-type</i>	<i>Salmonella typhi</i>	<i>Mycoplasma genitalium</i>
<i>Aspergillus fumigatus</i>	<i>Shigella flexneri</i>	<i>Mycoplasma hominis</i>
<i>Bacillus cereus</i>	<i>Staphylococcus aureus (mecA+)</i>	<i>Neisseria gonorrhoeae</i>
<i>Bartonella henselae</i>	<i>Staphylococcus aureus (mecA-)</i>	<i>Neisseria meningitidis sg A</i>
<i>Bartonella quintana</i>	<i>Streptococcus pneumoniae</i>	<i>Neisseria meningitidis sg B</i>
<i>BK virus</i>	<i>Toxoplasma gondii</i>	<i>Neisseria meningitidis sg C</i>
<i>Bordetella holmesii</i>	<i>Treponema pallidum</i>	<i>Papillomavirus type 16</i>
<i>Borrelia afzelii</i>	<i>Trichomonas vaginalis</i>	<i>Papillomavirus type 18</i>
<i>Borrelia burgdorferi</i>	<i>Trypanosoma cruzi</i>	<i>Parvovirus B19 (plasmid)</i>
<i>Borrelia garinii</i>	<i>Ureaplasma parvum</i>	<i>Vibrio cholerae</i>
<i>Brucella abortus</i>	<i>Ureaplasma urealyticum</i>	<i>Yersinia enterocolitica</i>
<i>Campylobacter jejuni</i>	<i>Varicella-zoster virus</i>	<i>Aeromonas caviae</i>
<i>Candida albicans</i>	<i>Francisella tularensis</i>	<i>Aeromonas dhakensis</i>
<i>Candida auris</i>	<i>Gardnerella vaginalis</i>	<i>Vibrio vulnificus</i>
<i>Chlamydia trachomatis</i>	<i>Haemophilus ducreyi</i>	<i>Aeromonas veronii</i>
<i>Chlamydomphila pneumoniae</i>	<i>Haemophilus influenzae</i>	<i>Aeromonas hydrophila</i>
<i>Chlamydomphila psittaci</i>	<i>Helicobacter pylori</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i>	<i>Herpes simplex 1</i>	<i>Shigella boydii</i>
<i>Coccidioides immitis</i>	<i>Herpes simplex 2</i>	<i>Shigella flexneri</i>
<i>Coxiella burnetii</i>	<i>HHV-6</i>	<i>Shigella dysenteriae</i>
<i>Cryptosporidium parvum</i>	<i>HHV-8</i>	<i>Campylobacter upsaliensis</i>
<i>Cytomegalovirus</i>	<i>Klebsiella pneumoniae (NDM-1)</i>	<i>Plesiomonas shigelloides</i>
<i>Enterococcus faecalis (vanB)</i>	<i>Legionella pneumophila</i>	<i>Campylobacter coli</i>
<i>Enterococcus faecium (vanA)</i>	<i>Leishmania chagasi</i>	<i>Escherichia coli O111/NM</i>
<i>Epstein-barr virus</i>	<i>Leishmania infantum</i>	<i>Plasmodium falciparum Strain 3D7</i>
<i>Escherichia coli (EAEC)</i>	<i>Listeria monocytogenes</i>	<i>Campylobacter fetus DNA</i>
<i>Escherichia coli (EIEC)</i>	<i>Moraxella catarrhalis</i>	<i>Yersinia pseudotuberculosis DNA</i>
<i>Escherichia coli (EPEC)</i>	<i>Mycobacterium avium</i>	<i>Adenovirus E04</i>
<i>Escherichia coli (ETEC)</i>	<i>Mycobacterium intracellulare</i>	<i>Adenovirus D20</i>
<i>Escherichia coli (VTEC)</i>	<i>Mycobacterium kansasii</i>	<i>Adenovirus B11</i>
<i>Rickettsia conorii</i>	<i>Mycobacterium tuberculosis</i>	<i>Adenovirus C02</i>
<i>Salmonella enteritidis</i>	<i>Mycobacterium ulcerans</i>	<i>Adenovirus A31</i>

- **Interference Study**

The presence of PCR inhibitors in the sample may induce a positive or negative interference on test results. The presence of various inhibitors in samples was tested.

The substances listed below were supplemented to fecal swab positive samples and tested on the BIOSYNEX AMPLIQUICK Fecal Virology:

Tested substances	Tested concentration	Interference observed
Betamethasone	2.5% v/v	None
Hemorrhoid cream (Titanoréine crème)	0.25% v/v	
Mucin	0.8% w/v	
Imodium (Lopéramide-Lyoc)	5% w/v	
Ampicillin sodium salt	152 µmol/L	
Fatty acid supplement	4.8%	
Pepto-Bismol	5%	
Starch	3% w/v	
Cellulose	4% w/v	
Pectin	0.01% w/v	
Lubricant	5% v/v	
Ethanol	0.2% v/v	
Urine	20% v/v	
Blood	40% v/v	

No positive or negative interference was found.

- **Precision**

The precision data in the context of the BIOSYNEX AMPLIQUICK Fecal Virology kit was determined based on 4 conditions:

- Intra-assay variability
- Inter-days variability
- Inter-operators variability
- Inter-laboratories variability
- Inter-batches variability

Intra-assay variability data for the BIOSYNEX AMPLIQUICK Fecal Virology PCR kit was determined on 30 replicate tests using negative and positive viral samples for each viral target. The Standard Deviation (SD) and Coefficient Variation (CV) determinations showed that the AMPLIQUICK Fecal Virology PCR kit is repeatable with CV<5%. Inter-days, inter-operators, inter-laboratories and inter-batches variabilities were determined on duplicate tests from medium and low concentrations of positive viral samples. The AMPLIQUICK Fecal Virology results are reproducible with CV<5%.

Rotavirus variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.80	0.16	0.56%	35.84	0.34	0.96%
Inter-days	28.44	0.41	1.46%	33.58	0.42	1.27%
Inter-operators	28.55	0.16	0.56%	34.67	0.82	2.37%
Inter-laboratories	28.45	0.15	0.53%	34.18	0.28	0.83%
Inter-batches	28.63	0.17	0.60%	36.25	0.32	0.89%

Sapovirus variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.17	0.37	1.37%	34.39	0.46	1.33%
Inter-days	28.72	0.48	1.68%	37.43	0.66	1.76%
Inter-operators	28.76	0.25	0.89%	37.23	0.23	0.61%
Inter-laboratories	28.65	0.12	0.42%	38.38	1.17	3.04%
Inter-batches	27.92	0.21	0.75%	35.44	0.46	1.29%

Astrovirus variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.83	0.30	1.05%	36.22	0.49	1.34%
Inter-days	28.05	0.53	1.90%	36.83	0.48	1.29%
Inter-operators	28.66	0.23	0.79%	37.33	0.27	0.71%
Inter-laboratories	28.1	0.15	0.53%	37.19	0.31	0.85%
Inter-batches	27.73	0.08	0.28%	35.00	0.21	0.60%

Norovirus GI variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.10	0.18	0.64%	36.29	0.78	2.15%
Inter-days	28.16	0.68	2.40%	31.80	0.48	1.52%
Inter-operators	28.31	0.91	3.21%	31.82	1.17	3.67%
Inter-laboratories	28.01	0.12	0.44%	32.17	0.14	0.43%
Inter-batches	29.19	0.74	2.52%	35.33	0.87	2.47%

Norovirus GII variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.26	0.40	1.40%	36.80	0.53	1.45%
Inter-days	27.33	0.49	1.79%	35.50	0.99	2.78%
Inter-operators	28.07	0.14	0.52%	35.79	0.76	2.11%
Inter-laboratories	27.37	0.02	0.06%	36.54	0.43	1.18%
Inter-batches	28.86	1.09	3.79%	37.37	1.16	3.11%

Adenovirus F40 variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	24.92	0.48	1.93%	34.63	0.56	1.62%
Inter-days	25.73	0.31	1.20%	35.77	0.68	1.90%
Inter-operators	26.05	0.15	0.57%	35.86	0.72	2.00%
Inter-laboratories	25.37	0.10	0.40%	34.93	0.02	0.05%
Inter-batches	26.80	0.59	2.21%	35.33	0.63	1.78%

Adenovirus F41 variability	"Medium" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	24.14	0.27	1.11%	35.01	0.68	1.95%
Inter-days	25.88	0.40	1.53%	35.54	0.09	0.24%
Inter-operators	26.14	0.59	2.27%	35.75	0.12	0.33%
Inter-laboratories	25.46	0.01	0.04%	36.11	0.89	2.46%
Inter-batches	26.93	0.44	1.64%	36.18	0.68	1.88%

## 11. LIMITATIONS

- Failure to follow any instruction of the IFU may adversely affect the test performance and/or invalidate the test result.
- This test cannot differentiate between Norovirus GI and GII as they are detected in the same fluorescent channel.
- This test cannot differentiate between Adenovirus F40 and F41 as they are detected in the same fluorescent channel.
- The possibility to obtain false positive or false negative results cannot be ruled out (e.g. in case of accidental contamination, in case of poor quality of the sample or in case the test procedure is not correctly followed). As with all diagnostic tests, results must be considered with other clinical information available to the physician.

## 12. LITERATURE

1. <https://www.who.int/fr/news-room/fact-sheets/detail/diarrhoeal-disease>
2. <https://www.who.int/fr/publications/m/item/vaccine-preventable-diseases-surveillance-standards-rotavirus>

3. Bányai, K., Estes, M.K., Martella, V., and Parashar, U.D. (2018). Viral gastroenteritis. *The Lancet* 392, 175–186. [https://doi.org/10.1016/S0140-6736\(18\)31128-0](https://doi.org/10.1016/S0140-6736(18)31128-0).
4. Corcoran, M.S., van Well, G.T.J., and van Loo, I.H.M. (2014). Diagnosis of viral gastroenteritis in children: interpretation of real-time PCR results and relation to clinical symptoms. *Eur J Clin Microbiol Infect Dis* 33, 1663–1673. <https://doi.org/10.1007/s10096-014-2135-6>.

### 13. INTERNATIONAL SYMBOLS

	Consult instructions for use or electronic instructions for use		Contains sufficient for <n> tests
	Catalogue number		<i>In vitro</i> diagnostic medical device
	Temperature limits		Use-by-date
	Manufacturer		Batch code
	Authorized representative in Switzerland		Do not use if package is damaged and consult the instructions for use
	Unique device identifier		Positive control
	Negative control		Control in Process
	Control Master mix		Keep away from sunlight
	Master mix		Cap strips
	Microplate		Qualitative detection by PCR
	Do not reuse		Importer

### 14. MANUFACTURER INFORMATION

**BIO SYNEX S.A.**  
 22 boulevard Sébastien Brant  
 67400 ILLKIRCH-GRAFFENSTADEN – France  
 Tel : +33 3 88 78 78 87

[www.biosynex.com](http://www.biosynex.com)

Contact France:

Tel.: +33 3 88 77 57 00

[service.clients@biosynex.com](mailto:service.clients@biosynex.com)

Contact other countries :

Tel.: +33 3 88 77 57 52

[sales@biosynex.com](mailto:sales@biosynex.com)

After-sales service:

Tel: +33 3 88 77 57 25

[Tech.support@biosynex.com](mailto:Tech.support@biosynex.com)

**BIO SYNEX SWISS S.A.**  
 Route de Rossemaison 100  
 2800 DELEMONT - Switzerland

**15. HISTORIC MODIFICATION**

<b>Record version</b>	<b>Amended paragraph</b>	<b>Change details</b>
V1	All	New product, creation of the IFU