**BIO SYNEX AMPLIQUICK LYSIS**

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## BIOSYNEX AMPLIQUICK LYSIS

NUCLEIC ACIDS RAPID RELEASE KIT FROM RESPIRATORY SAMPLES READY TO USE IN 96-WELL PRE-FILLED MICROPLATE.

*For professional in vitro diagnostic use only*

**REF** Ref: 3150059\_SEC01 / 3150059\_SEC02

### 1. INTENDED USE

BIOSYNEX AMPLIQUICK Lysis is a rapid lysis kit for the preparation of biological samples to access the nucleic acids. It comes in 96-well microplates pre-filled with ready-to-use lysis buffer. This kit is intended for use in *in vitro* molecular diagnostics by laboratory personnel.

This kit has been developed and validated for use only on nasopharyngeal swabs, nasopharyngeal/endonasal aspirates, or saliva samples collected with the BIOSYNEX AMPLIQUICK Sample Collection kit (ref 3150060 and related references). It is also validated for use on samples (urine, serum/plasma, oral, anal, vaginal or urethral swabs) treated with the BIOSYNEX AMPLIQUICK Sample Treatment kit. Samples prepared with the BIOSYNEX AMPLIQUICK Lysis kit can be used without nucleic acids extraction/purification steps with the following RT-PCR or PCR kits:

- BIOSYNEX AMPLIQUICK SARS-CoV-2 (ref 3150058 and related references) for the qualitative RT-PCR detection of SARS-CoV-2
- BIOSYNEX AMPLIQUICK Respiratory Triplex (ref 3150061\_SEC01 and related references) for the qualitative RT-PCR detection of SARS-CoV-2 (virus responsible for Covid-19), Influenza A and B, and Respiratory Syncytial Virus (RSV)
- BIOSYNEX AMPLIQUICK Bordetella (ref 3150068 and related references) for the qualitative PCR detection of *Bordetella pertussis*, *Bordetella parapertussis*, *Bordetella broncheseptica* or *Bordetella holmesii* in the context of whooping cough.
- BIOSYNEX AMPLIQUICK STI Panel I (ref 3150070 and related references) for the qualitative detection of *Chlamydia trachomatis*, *Neisseria gonorrhoeae* and *Mycoplasma genitalium* from a DNA extract obtained from urine samples, vaginal swabs urethral swabs, anal swabs or oral swabs.
- BIOSYNEX AMPLIQUICK Leptospira (ref 3150064 and related references) qualitative detection of *Leptospira interrogans* (the bacterium responsible for leptospirosis) from a DNA extract obtained from urine or blood samples (serum/plasma).

### 2. CLINICAL SUMMARY

COVID-19 is an acute respiratory illness caused by SARS-CoV-2, which belongs to the family of viruses responsible for Severe Acute Respiratory Syndrome (SARS). Approximately 2% of people with COVID-19 die.

Seasonal influenza A or B, caused by Influenza A or B viruses, is responsible for nearly half a million deaths worldwide each year.

Respiratory syncytial virus is a major cause of respiratory illness in young children in particular, causing upper and lower airway infections such as bronchiolitis.

The main four species of *Bordetella* responsible for respiratory diseases in humans also cause particularly severe respiratory symptoms in infants and children. The most common symptoms of these respiratory infections are the same: fever, fatigue, dry cough, runny or stuffy nose.

These diseases are spread from person to person by droplets from the nose or mouth when an infected person coughs or exhales. These droplets are deposited on objects and surfaces. Other people become infected when they come into contact with these objects or surfaces and then touch their eyes, nose or mouth, or when they breathe in the droplets from an infected person.

Leptospirosis is a bacterial disease that occurs worldwide. In humans, the disease is often benign with appropriate treatment, but it can lead to renal failure and even death in 5 to 20% of cases.

*Chlamydia trachomatis*, *Neisseria gonorrhoeae* and *Mycoplasma genitalium* are some of the most common sexually transmitted infections. Extra-genital infections and lesions are possible for all three bacteria.

For these infections, it is essential that the pathogen responsible for the symptoms be rapidly identified in order to take the appropriate therapeutic or even isolation measures necessary to control the spread and break the chain of infection. An accurate diagnosis is thus essential. The sensitivity of PCR tests allows for improved and faster diagnosis from different types of samples compared to traditional methods.

### 3. TEST PRINCIPLE

Tests based on the detection of the nucleic acids of specific pathogens allow the identification of infected persons. The BIOSYNEX AMPLIQUICK Lysis kit enables the release of RNA and DNA from biological samples in order to perform a molecular detection test to identify both asymptomatic and acute infections in patients.

The BIOSYNEX AMPLIQUICK Lysis kit is an *in vitro* assay based on the enzymatic digestion of protein elements in the samples collected with the BIOSYNEX AMPLIQUICK Sample Collection transport medium or treated with the BIOSYNEX AMPLIQUICK Sample Treatment kit. The Enzymatic digestion is performed by an endopeptidase. The proteolytic activity of the enzyme is performed during 10 minutes at 56°C in a thermocycler.

### 4. KIT CONTENT

Kit Reference	Material included	Material required but not included
3150059_SEC01 3150059_SEC02	<ul style="list-style-type: none"> <li>• 5 x 96-well divisible microplates prefilled with the lysis buffer</li> <li>• 5 bags containing cap strips</li> </ul>	<ul style="list-style-type: none"> <li>• Swabs &amp; Transport medium</li> <li>• RT-PCR or PCR Amplification Kit</li> <li>• Powder-free disposable gloves</li> <li>• Pipettes and filtered tips</li> <li>• Thermal Cycler (0.1 mL low-profile for ref 3150059_SEC01; 0.2mL high profile for ref 3150059_SEC02)</li> <li>• PCR microplates or tubes centrifuge</li> </ul>

### 5. PRECAUTIONS

- For *in vitro* diagnosis. For laboratory professional use only.
- Use the kit and all its components only before the expiry date.
- For optimal results, follow the procedure and storage conditions carefully.
- In the event of damage to the outer packaging only (no breakage or leakage), the kit remains usable. If the aluminum sealing foil is damaged (presenting some holes or unsticking parts), do not use the plate or strip concerned.
- Saliva samples for the detection of SARS-CoV-2 are validated in asymptomatic and symptomatic patients ≤6 days after appearance of the first symptoms only.
- Follow Good Laboratory Practices. Use disposable, powder-free laboratory gloves throughout the test procedure. Consider samples as potentially infectious and handle them with care, according to laboratory recommendations.
- Centrifuge the plates before opening and carefully remove the aluminum foil to avoid splashing of lysis buffer.
- The daily management of a large number of samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination, pre-handling PCR, post-PCR and DNA extraction should therefore be performed in separate rooms. The workflow in the laboratory should be unidirectional.
- Wear disposable gloves in each zone and change them before moving from one zone to another.

- Clean any splashes from the sample using an appropriate disinfectant.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste. Consider your local regulation on biological waste disposal.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Lysis, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. When in doubt, report.
- Safety data sheet available on request.

## 6. KIT STORAGE

- The kit is shipped at a temperature between 2 and 8°C. The components must be stored at this temperature at the time of arrival. Under these conditions the reagents are stable until the expiry date indicated.
- Remove from the refrigerator only the number of lysis buffer strips or plates required.
- Once the aluminum foil is removed, use the plate or strips immediately.

## 7. SAMPLE COLLECTION AND STORAGE

### When using the BIOSYNEX AMPLIQUICK Sample Collection kit

Use the BIOSYNEX AMPLIQUICK Sample Collection kit (ref 3150060 and related references) to collect nasopharyngeal or saliva samples with the provided sterile swabs and discharge into the transport medium tubes. Refer to § 7 and 8 of the instructions for use of BIOSYNEX AMPLIQUICK Sample Collection.

Perform nasopharyngeal/endonasal samples using a mucus aspirator in accordance with the manufacturer's recommendations. These aspirates must be used directly with the BIOSYNEX AMPLIQUICK lysis kit (ref. 3150059) without being unloaded into a transport medium.

Saliva samples can be collected in a clean and sterile pot directly, according to the requirements of the HAS (French National Authority for Health). In these conditions, saliva samples are stable 24h at 4°C. Samples discharged into the transport medium can be extracted immediately or within 4 hours if stored at room temperature, or within 24 hours if stored between 2 and 8°C, or in the long-term if frozen at -20°C or below.

**When using the BIOSYNEX AMPLIQUICK Sample Treatment kit** Use the BIOSYNEX AMPLIQUICK Sample Treatment kit (ref 3150063) to pre-treat urine, serum/plasma or swab samples (excluding respiratory infections). Please refer to §6 of the BIOSYNEX AMPLIQUICK Sample Treatment instructions.

The transport of clinical samples must comply with local regulations for the transport of infectious agents. If samples are transporting from the collection site to the processing site, centrifuge the sample collection tubes before opening them to ensure that there is no liquid in the caps.

## 8. PROTOCOL

1. Take a complete 96-well plate or the number of strips needed.
2. Centrifuge the plate or strips for 10 seconds to recover any droplets from the well edges or cap.
3. Carefully remove and discard the aluminum foil.

### From a sample discharged into the transport medium

- Vortex the sample tube.
- Pipet and deposit 50 µL of sample discharged in the transport medium per well and pipet up and down 3 times to homogenize the mix.
- Seal the wells with the provided cap strips or with a transparent PCR film.

### From pure saliva collected in a dry, clean and sterile sampling pot

- Pipet and deposit 50 µL of pure saliva per well and pipet up and down 3 times to homogenize the mix.
- Seal the wells with the provided cap strips or with a transparent PCR film.

#### From a sample pretreated with the BIOSYNEX AMPLIQUICK Sample kit Treatment

Once the pre-treatment steps have been carried out according to the kit protocol BIOSYNEX AMPLIQUICK Sample Treatment kit protocol:

- Vortex the sample tube.
- Take 50 µL of resuspended sample and place in a lysis well and pipette 3 times to homogenize the mixture.
- Close the wells with the caps provided or with transparent PCR film.

Briefly centrifuge the microplate or strips and place them in the thermal cycler using the following program:

Temperature	Duration
56°C	10 min
98°C	5 min
4°C	Until the microplate is removed from the system (if the device allows it)

Once the lysis step is completed, centrifuge the microplate or strips for 10 seconds to recover any condensation droplets. Gently remove the caps to avoid splashing.

If the use of the extracted nucleic acids is delayed, before adding it to the Master Mix, the microplate or strips can be stored between 2 and 8°C for 24h or frozen at -20°C or below for long term storage.

The treated samples can then be used with the BIOSYNEX AMPLIQUICK amplification kits indicated in paragraph 1 of this manual.

## 9. PERFORMANCE

This BIOSYNEX AMPLIQUICK Lysis kit was used to establish the performance of the BIOSYNEX AMPLIQUICK SARS-CoV-2 (ref 3150058 and related references), BIOSYNEX AMPLIQUICK Respiratory Triplex (ref 3150061\_SEC01 and related references) as well as the BIOSYNEX AMPLIQUICK Bordetella (ref 3150068 and related references).

## 10. INTERFERENCES

As for any protein, the enzyme is sensitive to chaotropic agents such as guanidium chloride, magnesium chloride or urea, among others. The use of buffers containing such solutes destabilizes the native state of the enzyme and may induce complete denaturation if present at high concentrations. In this case the enzymatic activity will be null and the release of RNA/DNA will not be possible.


## 11. LIMITS

The BIOSYNEX AMPLIQUICK Lysis kit has been developed and validated only for a use with samples collected with BIOSYNEX AMPLIQUICK Sample Collection kit or treated with BIOSYNEX AMPLIQUICK Sample Treatment. Saliva specimens can be added directly to the lysis buffer without use of the transport medium if sufficiently fluid to be pipetted (see §8). The quality and consistency of the samples play an important role in the efficiency of the enzymatic reaction. It is important, not to deposit mucus-rich or overly viscous samples into the lysis buffer. The volume and composition of the AMPLIQUICK Sample Collection transport medium or the resuspension buffer from the AMPLIQUICK Sample Treatment allow an optimal lysis reaction and subsequent amplification with one of the BIOSYNEX AMPLIQUICK amplification kits aforementioned; therefore, we do not recommend the use of other transport media or buffers with this kit.

## 12. BIBLIOGRAPHY

- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011;81:85-164.
- W. Guan et al. Clinical Characteristics of Coronavirus Disease 2019 in China NEngJMed (2020).
- Pittet LF, Emonet S, François P, Bonetti EJ, Schrenzel J, Hug M, Altwegg M, Siegrist CA, Posfay-Barbe KM. Diagnosis of whooping cough in Switzerland: Differentiating *Bordetella pertussis* from *Bordetella holmesii* by polymerase chain reaction. PLoS One. 2014; 9(2):e88936.
- Ñique AM, Coronado-Marquina F, Mendez Rico JA, García Mendoza MP, Rojas-Serrano N, Simas PVM, Cabezas Sanchez C, Drexler JF, Aboelhadid SM. A faster and less costly alternative for RNA extraction of SARS-CoV-2 using proteinase k treatment followed by thermal shock. PLoS One. 2021; 16(3)e0248885.
- Lever MA, Torti A, Eickenbusch P, Michaud AB, Santl-Temkiv T, Jørgensen BB. A modular method for the extraction of DNA and RNA, and the separation of DNA pools from diverse environmental sample types. Frontiers in Microbiology. 2015;6:476.

## 13. SYMBOLS

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Catalog number
	<i>In vitro</i> diagnostic medical device		Temperature limit		Do not reuse
	Manufacturer		Batch code		Use-by date
	Not for near patient testing		Do not use if package is damaged and consult instructions for use		Not for self-testing
	Importer		Unique Device Identifier		Authorized Representative in Switzerland
	Microplate		Bag of cap strips		

## 14. MANUFACTURER INFORMATION

**BIOSYNEX S.A.**  
 22 boulevard Sébastien Brant  
 67400 ILLKIRCH-GRAFFENSTADEN – France  
 Tel : +33 3 88 78 78 87  
[www.biosynex.com](http://www.biosynex.com)

Contact France:  
 Tel : +33 3 88 77 57 00  
[service.clients@biosynex.com](mailto:service.clients@biosynex.com)

Contact other countries :  
 Tel : +33 3 88 77 57 52  
[sales@biosynex.com](mailto:sales@biosynex.com)

After-sales service:  
 Tel : +33 3 88 77 57 25  
[Tech.support@biosynex.com](mailto:Tech.support@biosynex.com)

**BIOSYNEX SWISS S.A.**  
 Route de Rossemaison 100  
 2800 DELEMONT - Switzerland

Latest changes: Withdrawal of instructions for use from kit contents + modification §7

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