BIOSYNEX AMPLIQUICK® SAMPLE COLLECTION

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extraction of SARS-CoV-2 using proteinase k treatment followed by thermal shock. PLoS One. 2021; 16(3)e0248885.

- Lever MA, Torti A, Eickenbusch P, Michaud AB, Santl-Temkiv T, Jørgensen BB. A modular method for the extraction of DNA and RNA, and the separation of DNA pools from diverse environmental sample types. Frontiers in Microbiology. 2015;6:476.

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
BIOSYNEX AMPLIQUICK® SAMPLE COLLECTION

READY-TO-USE KIT FOR THE COLLECTION OF NASOPHARYNGEAL OR SALIVARY SAMPLES.

For professional *in vitro* diagnostic use only.**REF** Ref: 3150060, 3150060_ST**1. INTENDED PURPOSE**

BIOSYNEX AMPLIQUICK® Sample Collection is a non-automated nasopharyngeal and/or saliva sampling kit designed for the collection and conservation of specimens in the context of the of SARS-CoV-2, Influenza A & B, Respiratory Syncytial Virus (RSV) and Bordetella molecular diagnosis. It comes in tubes pre-filled with transport medium ready-to-use. This kit is intended for use in *in vitro* molecular diagnostics by laboratory personnel.

Samples collected with the kit can be used with BIOSYNEX AMPLIQUICK® Lysis kit (ref 3150059_SEC01 and related references) for nucleic acid release prior RT-PCR or PCR with one of the following kits:

- BIOSYNEX AMPLIQUICK® SARS-CoV-2 (ref 3150058 and related references) for the qualitative RT-PCR detection of SARS-CoV-2
- BIOSYNEX AMPLIQUICK® Respiratory Triplex (ref 3150061_SEC01 and related references) for the qualitative RT-PCR detection of SARS-CoV-2 (virus responsible for Covid-19), Influenza A and B, and Respiratory Syncytial Virus (RSV)
- BIOSYNEX AMPLIQUICK® Bordetella (ref 3150068 and related references) for the qualitative PCR detection of *Bordetella pertussis*, *Bordetella parapertussis*, *Bordetella broncheseptica* or *Bordetella holmesii* in the context of whooping cough.

2. INTRODUCTION

COVID-19 is an acute respiratory illness caused by SARS-CoV-2, which belongs to the family of viruses responsible for Severe Acute Respiratory Syndrome (SARS). Approximately 2% of people with COVID-19 die.

Seasonal influenza A or B, caused by Influenza A or B viruses, is responsible for nearly half a million deaths worldwide each year.

Respiratory syncytial virus is a major cause of respiratory illness in young children in particular, causing upper and lower airway infections such as bronchiolitis.

The main four species of Bordetella responsible for respiratory diseases like whooping cough in humans also cause particularly severe respiratory symptoms in infants and children.

The most common symptoms of these respiratory infections are the same: fever, fatigue, dry cough, runny or stuffy nose.

These diseases are spread from person to person by droplets from the nose or mouth when an infected person coughs or exhales. These droplets are deposited on objects and surfaces. Other people become infected when they come into contact with these objects or surfaces and then touch their eyes, nose or mouth, or when they breathe in the droplets from an infected person.


It is essential that the pathogen responsible for the symptoms be rapidly identified in order to take the appropriate therapeutic or even isolation measures necessary to control the spread and break the chain of infection.

3. PRINCIPLE

Tests based on the detection of the nucleic acids of specific pathogens such as SARS-CoV-2, Influenza A and B, RSV or Bordetella allow the identification of infected persons. The BIOSYNEX AMPLIQUICK® Sample Collection kit enables the collection and storage of respiratory and/or saliva specimens for a use with the BIOSYNEX AMPLIQUICK® Lysis prior to RT-PCR or PCR amplification.

The BIOSYNEX AMPLIQUICK® Sample Collection kit is composed of a buffer with no detectable DNAses, RNAses, proteases or phosphatases. It has been developed to ensure on one hand optimal results of the enzymatic digestion for RNA and DNA release when using the AMPLIQUICK® Lysis kit, and on the other hand solubilization, stabilization and protection of the nucleic acids from degradation.

4. KIT CONTENT

Kit Reference	Material included	Material required but not included
3150060 3150060_ST	<ul style="list-style-type: none"> •100 tubes of transport medium •100 swabs (ref. 2123-1029 - sterilised by irradiation - CE0197)  Citotest Labware Manufacturing Co., Ltd No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China  Wellkang Ltd Enterprise Hub, NW Business Complex, 1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland <ul style="list-style-type: none"> •1 Instructions for use 	<ul style="list-style-type: none"> •Powder-free disposable gloves •BIOSYNEX AMPLIQUICK® Lysis kit (ref 3150059_SEC01 and related references) •RT-PCR or PCR amplification kit

5. PRECAUTIONS

- For *in vitro* diagnosis. For laboratory professional use only.
- Use the kit and its components only before the expiry date. All components are for single use.
- For optimal results, follow the procedure and storage conditions carefully.
- In the event of damage to the outer packaging only (no breakage or leakage), the kit remains usable. If one of the tubes is damaged (cracked or broken) or badly capped, do not use it.
- Do not use the swab if there is a breakage of the pouch.
- Saliva samples for SARS-CoV-2 detection are validated in asymptomatic or symptomatic ≤6 days after onset of first symptoms only.
- Follow good laboratory practices. Wear a lab coat, disposable non-powdered laboratory gloves and eye protection throughout the nasopharyngeal or saliva sampling procedure.
- Consider samples as potentially infectious and handle them with care, according to laboratory recommendations.
- Wear disposable gloves in each zone and change them before moving from one zone to another.
- Clean any splashes from the sample using an appropriate disinfectant.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste. Consider your local regulation on biological waste disposal.
- If, in relation to the use of BIOSYNEX AMPLIQUICK® Sample Collection, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. When in doubt, report.
- Safety data sheet available on request.

6. CONSERVATION AND STABILITY

The kit is shipped at a temperature between 2°C and 30°C and the components must be stored at this temperature range at the time of arrival. Under these conditions the reagents are stable until the indicated expiry date.

7. PROTOCOL

Before opening the transport medium tube, centrifuge it to ensure that the entire volume is at the bottom, or alternatively, hold the tube vertically and tap it several times on a flat surface to bring any drops in the cap or on the walls to the bottom.

Nasopharyngeal sampling

1. Keep the patient's head tilted back.
2. Remove the swab from its packaging. Do not touch the head of the swab.
3. Carefully insert the swab horizontally into the nostril of the patient, reaching the surface of posterior nasopharynx.
4. Swab the posterior wall of the nasopharynx, turning the swab several times.
5. Remove the swab from the nasal cavity.
6. Then, immerse the swab in the transport medium contained in the tube.
7. Unload the sample by rotating the swab in the transport medium at least 15 times while pressing the head against the bottom and side of the tube.
8. Break the swab or remove it making sure to wring it out against the wall of the tube and dispose of it in a bin suitable for biological waste.
9. Close tightly the lid.

Nasopharyngeal/endonasal aspirates

Collect aspirates with a mucus aspirator according to the manufacturer's recommendations and discharge into a Sample Collection tube.

Saliva sampling using the material provided in this kit

1. Keep the patient's head tilted back.
2. Remove the swab from its packaging. Do not touch the head of the swab.
3. Carefully insert the swab on or under the patient's tongue.
4. Allow the swab to soak up the saliva for 1 minute.
5. Remove the swab from the oral cavity.
6. Then, immerse the swab in the transport medium contained in the tube.
7. Unload the sample by rotating the swab in the transport medium at least 15 times while pressing the head against the bottom and side of the tube.
8. Break the swab or remove it making sure to wring it out against the wall of the tube and dispose of it in a bin suitable for biological waste.
9. Close tightly the lid.

Saliva sampling in a dry, clean and sterile sampling pot

1. Ask the patient to salivate.
2. Remove the lid and collect the saliva in the sampling pot (no sputum).
3. Close tightly the lid.
4. At the laboratory, open the sampling pot under a microbiological safety cabinet.
5. If the saliva is fluid enough to be pipetted, pipet 50µL and add it directly into BIOSYNEX AMPLIQUICK® Lysis Buffer.
6. If the saliva is too thick to pipet, use the provided swab: remove the swab from its packaging. Do not touch the head of the swab.
7. Immerse the swab in the saliva and allow to soak for 1 minute.
8. Then, immerse the swab in the transport medium contained in the tube.
9. Unload the sample by rotating the swab in the transport medium at least 15 times while pressing the head against the bottom and side of the tube.
10. Break the swab or remove it making sure to wring it out against the wall of the tube and dispose of it in a bin suitable for biological waste.
11. Close tightly the lid.

N.B: The saliva sample should be taken more than 30 minutes after the last drink, food, cigarette / e-cigarette, tooth brushing or oral rinse, according to the recommendations of the HAS (French National Authority for Health).

8. SAMPLE STORAGE

Once collected, samples can be extracted immediately. Before treatment, they can also be stored up to 6 hours at room temperature (between 15 to 25°C), or 24h between 2°C and 8°C; alternatively, for long-term storage samples can be frozen at -20°C or below to ensure nucleic acid stability. Saliva samples

collected in a clean and sterile pot, according to the requirements of the HAS (French National Authority for Health), are stable 24h at 4°C.

The transport of clinical specimens must comply with local regulations for the transport of infectious agents.

The specimens discharged into the transport medium are intended to be used with the BIOSYNEX AMPLIQUICK® Lysis Kit (ref 3150059_SEC01 and related references).

9. PERFORMANCES

This BIOSYNEX AMPLIQUICK® Sample Collection kit was used to establish the performance of the BIOSYNEX AMPLIQUICK® SARS-CoV-2 (ref 3150058 and related references), BIOSYNEX AMPLIQUICK® Respiratory Triplex (ref 3150061_SEC01 and related references) as well as the BIOSYNEX AMPLIQUICK® Bordetella (ref 3150068 and related references) amplification kits on non-extracted samples treated with BIOSYNEX AMPLIQUICK® Lysis (ref 3150059_SEC01 and related references).

10. INTERFERENCES

There are no known interferences to report.

11. LIMITS

The BIOSYNEX AMPLIQUICK® Sample Collection kit has been developed and validated only for a use with BIOSYNEX AMPLIQUICK® Lysis kit (ref 3150059_SEC01 and related references). If following the diagnosis, other steps are required such as screening of variants or sequencing, we recommend adding 1mL of molecular grade water or saline to have the adequate volume.

The quality and consistency of the samples play an important role in the efficiency of the nucleic acid release. Therefore, we do not recommend leaving the swabs in the transport medium, if possible, as that increases the viscosity of the samples.

12. BIBLIOGRAPHY

- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011;81:85-164.
- W. Guan et al. Clinical Characteristics of Coronavirus Disease 2019 in China *N Engl J Med* (2020).
- Pittet LF, Emonet S, François P, Bonetti EJ, Schrenzel J, Hug M, Altwegg M, Siegrist CA, Posfay-Barbe KM. Diagnosis of whooping cough in Switzerland: Differentiating *Bordetella pertussis* from *Bordetella holmesii* by polymerase chain reaction. *PLoS One*. 2014; 9(2):e88936.
- Ñique AM, Coronado-Marquina F, Mendez Rico JA, García Mendoza MP, Rojas-Serrano N, Simas PVM, Cabezas Sanchez C, Drexler JF, Aboelhadid SM. A faster and less costly alternative for RNA extraction of SARS-CoV-2 using proteinase k treatment followed by thermal shock. *PLoS One*. 2021; 16(3):e0248885.
- Lever MA, Torti A, Eickenbusch P, Michaud AB, Santl-Temkiv T, Jørgensen BB. A modular method for the extraction of DNA and RNA, and the separation of DNA pools from diverse environmental sample types. *Frontiers in Microbiology*. 2015;6:476.

SYMBOLS

	Consult instructions for use or electronic instructions for use		Contains sufficient for <n> tests		Catalog number
	<i>In vitro</i> diagnostic medical device		Temperature limit		Do not reuse
	Manufacturer		Batch code		Use-by date
	Date of manufacture		Do not use if package is damaged and consult instructions for use		Not for self-testing
	Not for near patient testing		Unique Device Identifier		Authorized Representative in Switzerland
	Importer		Sterilized using irradiation		Authorized Representative in the European Community
	Caution		Sample collection tube		Swab
	Instructions for use				


MANUFACTURER INFORMATION

BIO SYNEX S.A.
 22 boulevard Sébastien Brant
 67400 ILLKIRCH-GRAFFENSTADEN – France
 Fax : +33 3 88 78 76 78
www.biosynex.com

Contacts for France:
 Tel.: +33 3 88 78 78 87
client.pro@biosynex.com

Contacts for other countries:
 Tel. : +33 3 88 77 57 52
export@biosynex.com

BIO SYNEX SWISS S.A.
 Route de Rossemaison 100
 2800 DELEMONT - Switzerland

Latest changes: Creation of multilingual instructions for use

BIO SYNEX AMPLIQUICK® SAMPLE COLLECTION

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- Pittet LF, Emonet S, François P, Bonetti EJ, Schrenzel J, Hug M, Altwegg M, Siegrist CA, Posfay-Barbe KM. Diagnosis of whooping cough in Switzerland: Differentiating *Bordetella pertussis* from *Bordetella holmesii* by polymerase chain reaction. *PLoS One*. 2014; 9(2):e88936.

- Nique AM, Coronado-Marquina F, Mendez Rico JA, García Mendoza MP, Rojas-Serrano N, Simas PVM, Cabezas Sanchez C, Drexler JF, Aboelhadid SM. A faster and less costly alternative for RNA extraction of SARS-CoV-2 using proteinase k treatment followed by thermal shock. 2021; 16(3)e0248885.

- Lever MA, Torti A, Eickenbusch P, Michaud AB, Santl-Temkiv T, Jørgensen BB. A modular method for the extraction of DNA and RNA, and the separation of DNA pools from diverse environmental sample types. *Frontiers in Microbiology*. 2015;6:476.

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BIOSYNEX S.A.
BIOSYNEX SWISS S.A.