## BIOSYNEX AMPLIQUICK *Leptospira*

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

- Urine sample with more than 2% blood may interfere with the test result. In case of blood in the urine sample, consider the result with caution. Serum and plasma samples should not contain blood.

### 13. BIBLIOGRAPHIE

- [www.pasteur.fr/fr/centre-medical/fiches-maladies/leptospirose](http://www.pasteur.fr/fr/centre-medical/fiches-maladies/leptospirose)
- Miotto BA, da Hora AS, Taniwaki SA, Brandão PE, Heinemann MB, Hagiwara MK. Development and validation of a modified TaqMan based real-time PCR assay targeting the lipL32 gene for detection of pathogenic *Leptospira* in canine urine samples. *Braz J Microbiol.* 49:3:584-590. July 2018.
- Woods K, Nic-Fhogartaigh C, Arnold C, Boutthasavong L, Phuklia W et al. A comparison of two molecular methods for diagnosing leptospirosis from three different sample types in patients presenting with fever in Laos. *Clin Microbiol Infect.* 24:9:1017.1-1017-7. Sept 2018.
- O Munch M, Chambers LC, Manhart LE, Domogala D, Lopez A, Fredricks DN, Srinivasan S. Optimizing bacterial DNA extraction in urine. *PLOS ONE.* 14:9. Sept 2019.
- Lucchesi PMA, Arroyo GH, Etcheverria AI, Parma AE, Seijo AC. Recommendations for the detection of *Leptospira* in urine by PCR. 37:2:131-134. May 2004.
- Stoddard RA, Gee JE, Wilkins PP, McCaustland K, Hoffmaster AR. Detection of pathogenic *Leptospira* spp. through TaqMan polymerase chain reaction targeting the LipL32 gene. *Diagn Microbiol Infect Dis.* 64:3:247-255. July 2009.
- Cyrille Goarant. Leptospirosis: risk factors and management challenges in developing countries. *Res Rep Trop Med.* 2016 Sep 28;7:49-62. doi: 10.2147/RRTM.S102543. PMID: 30050339; PMCID: PMC6028063.

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## BIOSYNEX AMPLIQUICK *Leptospira*

PCR FOR THE QUALITATIVE DETECTION OF *LEPTOSPIRA INTERROGANS* IN URINE, SERUM OR PLASMA

*In vitro* diagnostic medical device for professional use.

**REF** 3150064

### 1. INTENDED PURPOSE

BIOSYNEX AMPLIQUICK *Leptospira* is a molecular *in vitro* diagnostic test for the qualitative detection by qPCR of *Leptospira interrogans*. The test is an aid for the diagnosis of *Leptospirosis* performed using a DNA extract obtained from urine, serum or plasma samples. This test detects a molecular target of *Leptospira* located on the lipL32 gene encoding an outer membrane protein present only in pathogenic strains. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

### 2. CLINICAL SUMMARY

*Leptospirosis* is a bacterial disease that occurs worldwide. The disease is mainly caused by pathogenic bacteria of the species *Leptospira interrogans*. Its main reservoirs are rodents, especially rats, which excrete the bacteria in their urine. This bacterium is easily maintained in the environment (fresh water or muddy soils) which favors contamination. In humans, the disease is often benign with appropriate treatment.

Incubation lasts an average of 4 to 14 days. The disease can evolve into different clinical forms. In the moderate form, it resembles an influenza-like syndrome with fever, headache, joint and muscle pain. Without appropriate treatment, it can evolve into more serious forms with multivisceral damage that can lead to renal failure, neurological damage, haemorrhage and even death.

In France, the most exposed populations are agricultural workers, sewage workers, and more generally those who are regularly in contact with contaminated environments and surfaces such as water, soil or plants. There is also an incidence in people practicing freshwater leisure activities.

Clinical diagnosis can be confirmed by (real-time) PCR amplification of the bacterial genome. For serum or plasma samples, PCR is only of interest during the first week following the appearance of symptoms. After this period, serology should be performed to look for antibodies. For urine samples, PCR can be performed from the second week after the onset of fever.

Early diagnosis and antibiotic treatment reduce the risk of complications and shorten the course of the disease.

### 3. TEST PRINCIPLE

The BIOSYNEX AMPLIQUICK *Leptospira* kit is an *in vitro* diagnostic test based on real-time Polymerase Chain Reaction (PCR) technology by fluorescent probe hydrolysis. It allows the specific amplification and simultaneous detection of molecular targets carrying the sequences of interest: the lipL32 sequence of the bacterial DNA, an exogenous internal control CIEZ, as well as the endogenous internal control RNase P. This method uses probes labelled with the FAM fluorophore for the specific amplification of the bacterial gene, a probe labelled with the HEX fluorophore for the CIEZ exogenous control allowing to validate the amplification reaction, thus excluding false negative results, and a third probe labelled with the Cy5 fluorophore used for the amplification of the internal control RNase P region coding for highly conserved human RNA. Thus, this second procedural control allows to evaluate the quality of the biological sample collection, and to validate the DNA extraction and purification steps.

The increase of the fluorescence signal is only detected if the target sequence complementary to the amplified probe is present in the sample. The fluorescent signal is therefore directly proportional to the amplification of the target during the amplification phase. The C<sub>q</sub> (cycle quantification) value corresponds to the number of cycles at which the fluorescence starts to increase exponentially in contrast to the background. The kit can be used on three types of samples:

- serum
- plasma
- urine

As a first line of approach, we recommend that tests be performed using serum or plasma (Goarant, 2016). For urine samples, we recommend whenever possible to perform a DNA extraction/purification rather than using the Sample Treatment method as working from extracted DNA provides better sensitivity.

Prior to their use with the BIOSYNEX AMPLIQUICK Leptospira amplification kit (§9), samples must be first collected and treated according to the procedures detailed in §7 and §8 using BIOSYNEX AMPLIQUICK Sample Treatment kit (ref 3150063 and related references) or processed to obtain purified DNA extracts using a suitable kit.

## 4. KIT CONTENT

### Material included

- 1 Microtube containing 325µL of ready-to-use Master Mix (25 reactions)
- 1 Positive Control (CTRL +, red cap)
- 1 Negative Control (CTRL -, green cap)

### Material required but not included

- DNA extraction kit
- Powder-free disposable gloves
- Pipettes and filtered tips
- Real-time PCR Thermal Cycler
- PCR microplate or microtubes centrifuge
- PCR tubes/strips/plate

The PCR thermal cycler used for the test must have the following main characteristics:

- Open system
- Real-time quantitative PCR assays.
- Programmable Thermal Cycler block
- Excitation source : Leds, lamp or laser
- Filter sets (excitation/emission wavelengths) suitable for the detection of "reporter" fluorophores of the FAM, HEX and Cy5 probes.
- Connection with a computer using specific analysis software that allow the recovery of fluorescence data, absolute quantification assays, and the interpretation of results.

The kit has been developed and validated with the following real-time thermal cyclers:

- CFX96 Touch Real-Time PCR Detection System (Bio-Rad)
- CFX96 Opus Real-Time PCR Detection System (Bio-Rad)
- QuantGene 9600 (Bioer)
- LightCycler 480 (Roche)
- Amplix-DT lite 48 (DNA Technology).

If another thermal cycler is used, please perform your own validation of the BIOSYNEX AMPLIQUICK Leptospira kit using the provided controls and/or appropriate qualified samples before using the test.

## 5. PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of a large number of samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination. Pre-handling PCR, post-PCR and DNA extraction should therefore be performed in separate rooms. Wear disposable gloves in each zone and change them before moving from one zone to another.
- The test and cap are for single use only. Do not reuse. Do not open the PCR tubes at the end of the test.
- Do not use the test if aluminum foil is opened or damaged. Once the aluminum foil is removed, use the test immediately.
- Do not use the kit if it arrived unfrozen.
- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.

- Protect the kit from light.
- Centrifuge the tubes before opening, opening them one after the other and closing them well between each one to avoid any contamination.
- The CTRL+ contains significant amounts of DNA sequences. It can therefore potentially contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store this component outside the kit at the first opening of the kit.
- The negative and positive controls included in the kit mimic results obtained with negative or positive samples respectively. They must be used with each new run.
- To ensure that there is no contamination during handling, it is left to the user to include as part of their internal quality approach, at the frequency chosen, an experiment negative control well in which 8µL of molecular biology grade water are added to the Master mix.
- When using negative and positive controls with a series of patients, it is recommended to first deposit the negative control, then deposit the patient samples, and to finish with the deposit of the positive control.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste.
- The device contains material of bacterial or animal origin and may transmit infectious agents and should be handled with caution.
- Blood may interfere with test results. Sample should not contain blood.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Leptospira, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed

## 6. REAGENT STORAGE AND STABILITY

The kit is shipped frozen. Kit components must arrive frozen. Store the kit at a temperature of -20°C. Under these conditions, the reagents are stable until the expiry date indicated on the kit label. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light.

The kit can undergo 10 freeze/thaw cycles without the shelf-life and qualities of the product are altered.

## 7. TEST PRINCIPLE

BIOSYNEX AMPLIQUICK Leptospira allows the amplification of Leptospira DNA from serum, plasma, or urine samples. The test can be performed using purified DNA extract from the samples, or samples collected and treated with the BIOSYNEX AMPLIQUICK Sample Treatment kit (ref 3150063 and related references).

Urine samples:

Urine should be collected in a sterile, preservative-free dry collection container. Before analysis urine samples can be stored up to 8h at RT (room temperature, +21°C/+2°C) and 24h at +2-8°C. After the pre-treatment step, the pellet resuspended in resuspension buffer from the BIOSYNEX AMPLIQUICK Sample Treatment kit, can be stored up to 24h at RT or +2-8°C.

Serum/plasma samples:

The blood sample should be collected in an EDTA or heparinized tube.

Before analysis plasma and serum samples can be stored up to 8h at RT and 7 days at +2-8°C. After the centrifugation steps, plasma and serum pellets resuspended in the sample treatment buffer from the BIOSYNEX AMPLIQUICK Sample Treatment kit can be stored up to 8h at RT and 7 days at +2-8°C for plasma and 5 days at +2-8°C for serum.

Transportation of clinical specimens must comply with local regulations for the transport of infectious agents.

## 8. EXTRACTION OF NUCLEIC ACIDS

The nucleic acids extraction must be performed before the amplification protocol using a suitable extraction system for bacterial DNA. Please follow the manufacturer's instructions.

The BIOSYNEX AMPLIQUICK Leptospira kit has been validated with the following extraction kits:

- BIOSYNEX AMPLIQUICK Sample treatment (Biosynex - ref. 3150063 and related references)
- QIAamp DNA Mini kit (Qiagen – Ref 51306)
- NucleoMag Dx Pathogen (Macherey Nagel – Ref 744215)
- NucleoMag Pathogen (Macherey Nagel – Ref 744210)

- Nucleic Acid Extraction Reagents (DNA/RNA) kit (Singuway Biotech – Ref MTQM036).

If another DNA extraction method or kit is used, please perform your own method validation of the BIOSYNEX AMPLIQUICK Leptospira kit using qualified respiratory samples before using the test (do not use the controls provided).

It is not necessary to extract the positive and negative controls with the Nucleic Acid Extraction Kit.

## 9. AMPLIFICATION PROTOCOL

It is recommended that Master Mix deposits and sample additions be performed on a refrigerated rack or on ice.

1. Remove the Master mix tube from the freezer and let it thaw briefly at room temperature or on ice, then vortex.
2. Centrifuge the tube for 10 seconds to collect any droplets from the sides or cap.
3. Dispense 12  $\mu$ L of Master mix per PCR well.
4. Add 8  $\mu$ L of sample or control.
5. Close the wells with appropriate caps.
6. Centrifuge for 10 seconds.
7. Place them in the thermal cycler and run the following amplification program.

PCR Program :

Step	Repetitions	Temperature	Duration	Acquisition
Taq Activation	1x	95°C	3 min	-
Denaturation	11x	95°C	15 sec	-
Hybridation / Elongation		69.5°C	1 min	-
Denaturation	45x	95°C	15 sec	-
Hybridation / Elongation		59°C	1 min	Yes

Enter **20  $\mu$ L** of reaction volume into the Thermal Cycler program.

Please refer to the operating instructions of the Thermal Cycler used for programming information.

Setting the detection channels:

Target	Fluorochrome
LipL32 gene of Leptospira	FAM
Internal control CIEZ	HEX
Procedural Control gene RNase P	Cy5

## 10. DATA ANALYSIS AND INTERPRETATION OF RESULTS

### A. Trial Validation Criteria

#### Negative control :

The fluorescence emitted must be below the threshold. This is a non-specific amplification indicator. If the fluorescence is above the threshold, check for an atypical curve. In the case of an amplification curve, a contamination or a distribution error in the microtubes is to be considered. Only the internal and procedural controls signals should be amplified.

#### Positive control:

The value of the positive control should preferably be detected before 30 cycles ( $Cq \leq 30$ ). In the absence of amplification of the positive control, the existence of an amplification or fluorescence detection problem (defective thermal cycler or instrument non-adapted to the method) should be considered.

#### Internal amplification control:

The exogenous CIEZ amplification control ensures that the enzymes in the master mix are functional and validates the amplification reaction. The RNase P procedural control validates the quality of the biological sample collection and validates the DNA extraction and purification steps.

Indeed, an amplification curve of the exogenous internal control and the procedural control must be observed in the HEX and Cy5 channels.

Nevertheless, two situations of lack of amplification of the internal control can be observed:

- If the target gene is initially present in the sample with a high number of copies, the internal controls may not be amplified. This result is consistent and does not invalidate the test. It should be interpreted as a positive result despite the lack of signal from the internal controls. This phenomenon is the result of amplification competition between the internal controls and targets present at high copy numbers.
- If the target gene in the FAM channel is not amplified, as well as the internal controls in HEX and Cy5 channels, then no result can be rendered. This situation highlights the presence of PCR inhibitors. PCR should be repeated starting from the primary sample and preferably on DNA extract

#### B. Qualitative interpretation (positive or negative)

Signals above the threshold, **and visually consistent with a classical PCR amplification curve**, are considered positive results.

Some samples may show atypical curves that are not characteristic of amplification curves. In this case, the result should not be considered as interpretable and repeat the analysis of the sample with the controls.

Detection channels			Interpretation
FAM (gene lipL32)	HEX (internal ctrl CIEZ)	Cy5 (RNase P)	
-	+	+	Negative control
+	+	+	Positive control
+	+	+	Patient with <i>Leptospira</i> specific DNA
+	+	_*	
+	-	+	
+	_*	_*	
-	+	+	Patient with no <i>Leptospira</i> specific DNA
-	+	-	Invalid result: poor sample quality, inhibition of reaction, or problems during testing. Perform a new test or re-sampling.
-	-	-	

\*In case of positive signals in the channel (FAM) for pathogen target detection, the signal from the internal controls is not required to validate the result. A high pathogen load may result in a reduced or absent signal for the internal amplification control and/or the procedural control..

## 11. PERFORMANCES

### Analytical sensitivity

The detection limit (LoD) of the BIOSYNEX AMPLIQUICK *Leptospira* kit is defined as the concentration that can be detected at 95% at least on a specific *Leptospira interrogans* DNA sample.

DNA extracts from qualified samples were obtained after purification with the QIAamp DNA Mini kit (Qiagen).

The non-extracted samples (urine, serum and plasma) were prepared with the BIOSYNEX AMPLIQUICK Sample Treatment kit (ref. 3150063).

The LoD<sub>95</sub> has been statistically determined, based on results obtained by performing a dilution series of a reference sample with a known number of bacteria

Sample type	LoD <sub>95</sub> ( <i>Leptospira</i> /mL)
DNA extract (all sample types)	467.4
Urine (without DNA extraction)	1077.07
Serum (without DNA extraction)	561.26
Plasma (without DNA extraction)	450.094

### Analytical specificity

#### Specificity for the lipL32 gene

The design of oligonucleotides (primers and probes) was validated *in silico* by BLAST alignment. The comparison of the sequences obtained shows a specific detection of lipL32 gene from *Leptospira*.

No primers or probes detect bacterial DNA other than that of *Leptospira* targeted.

23 DNAs from pathogenic and saprophytic *Leptospira* strains expressing lipL32 including 3 in which the sequence is poorly conserved (-), were tested with the BIOSYNEX AMPLIQUICK *Leptospira* kit.

<i>L. mayottensis</i>	+	<i>L. gomenensis</i>	+
<i>L. alexanderi</i>	+	<i>L. putramalaysiae</i>	+
<i>L. licerasiae</i>	+	<i>L. tipperaryensis</i>	+
<i>L. fainei</i>	+	<i>L. borgpetersenii</i>	+
<i>L. kirschneri</i>	+	<i>L. interrogans</i>	+
<i>L. wolffii</i>	+	<i>L. noguchii</i>	+
<i>L. kmetyi</i>	+	<i>L. santarosai</i>	+
<i>L. alstonii</i>	+	<i>L. weilii</i>	+
<i>L. adleri</i>	+	<i>L. inadai</i> (-)	-
<i>L. barantonii</i>	+	<i>L. biflexa</i> (-)	-
<i>L. ellisii</i>	+	<i>L. wolbachii</i> (-)	-
<i>L. dzianensis</i>	+		

### Cross-reactivity

A panel of 79 DNA samples and 38 RNA samples from a biobank listed in the following tables were tested with the BIOSYNEX AMPLIQUICK Leptospira kit. For all these samples, no amplification of the Leptospira target (lipL32 gene) was observed.

DNA		
<i>Acanthamoeba castellanii</i>	<i>Enterococcus faecium</i> (VanA)	<i>Mycobacterium ulcerans</i>
<i>Adenovirus</i>	<i>Epstein-Barr Virus</i>	<i>Mycoplasma genitalium</i>
<i>Adenovirus 41</i>	<i>Escherichia coli</i> (EAEC)	<i>Mycoplasma hominis</i>
<i>Aspergillus fumigatus</i>	<i>Escherichia coli</i> (EIEC)	<i>Neisseria gonorrhoeae</i>
<i>Bacillus cereus</i>	<i>Escherichia coli</i> (ETEC)	<i>Neisseria meningitidis</i> Sg A
<i>Bartonella henselae</i>	<i>Escherichia coli</i> (VTEC)	<i>Neisseria meningitidis</i> Sg B
<i>Bartonella Quintana</i>	<i>Francisella tularensis</i>	<i>Neisseria meningitidis</i> Sg C
<i>Bk Virus</i>	<i>Gardnerella vaginalis</i>	<i>Papillomavirus type 16</i>
<i>Bordetella holmesii</i>	<i>Giardia intestinalis</i>	<i>Papillomavirus type 18</i>
<i>Bordetella parapertussis</i>	<i>Haemophilus ducreyi</i>	<i>Parvovirus B19</i> (Plasmid)
<i>Bordetella pertussis</i>	<i>Haemophilus influenzae</i>	<i>Rickettsia conorii</i>
<i>Borrelia afzelii</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>
<i>Borrelia garinii</i>	<i>Herpes simplex 1</i>	<i>Salmonella typhi</i>
<i>Borrelia burgdorferi</i>	<i>Herpes simplex 2</i>	<i>Staphylococcus aureus</i> (MecA-)
<i>Brucella abortus</i>	<i>Hhv-6</i>	<i>Staphylococcus aureus</i> (MecA+)
<i>Campylobacter jejuni</i>	<i>Hhv-8</i>	<i>Streptococcus agalactiae</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae</i> (NDM-1)	<i>Toxoplasma gondii</i>
<i>Candida auris</i>	<i>Legionella pneumophila</i>	<i>Treponema pallidum</i>
<i>Chlamydia trachomatis</i>	<i>Leishmania chagasi</i>	<i>Trichomonas vaginalis</i>
<i>Chlamydophila pneumoniae</i>	<i>Leishmania infantum</i>	<i>Trypanosoma cruzi</i>
<i>Chlamydophila psittaci</i>	<i>Listeria monocytogenes</i>	<i>Ureaplasma parvum</i>
<i>Clostridium difficile</i>	<i>Moraxella catarrhalis</i>	<i>Ureaplasma urealyticum</i>
<i>Coccidioides immitis</i>	<i>Mycobacterium avium</i>	<i>Varicella-Zoster Virus</i>
<i>Coxiella burnetii</i>	<i>Mycobacterium intracellulare</i>	<i>Vibrio cholerae</i>
<i>Cryptosporidium parvum</i>	<i>Mycobacterium kansasii</i>	<i>Yersinia enterocolitica</i>
<i>Cytomegalovirus</i>	<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pneumoniae</i>
<i>Enterococcus faecalis</i> (VanB)		

RNA	
<i>Coronavirus Oc43</i>	<i>Influenza A H5</i>
<i>Coronavirus</i>	<i>Influenza B</i>
<i>Coronavirus SARS (2003)</i>	<i>Measles</i>
<i>Coxsackie A6</i>	<i>MERS Coronavirus</i>
<i>Coxsackie B1</i>	<i>Mumps</i>
<i>Coxsackie B5</i>	<i>Norovirus</i>
<i>Dengue 1 Virus</i>	<i>Novel Influenza A H1n1</i>

<i>Dengue 2 Virus</i>	<i>Parainfluenza 1</i>
<i>Dengue 3 Virus</i>	<i>Parainfluenza 2</i>
<i>Dengue 4 Virus</i>	<i>Parainfluenza 3</i>
<i>Echovirus 5</i>	<i>Parainfluenza 4 A</i>
<i>Enterovirus 68</i>	<i>Respiratory Syncytial Virus (Subtype A)</i>
<i>Rhinovirus</i>	<i>Respiratory Syncytial Virus (Subtype B)</i>
<i>Rotavirus</i>	<i>West Nile Virus</i>
<i>Rubella</i>	<i>Yellow Fever Virus</i>
<i>Tick-Borne Encephalitis Virus</i>	<i>Zika Virus (Asian Lineage)</i>
<i>Human Parainfluenza 1</i>	<i>Zika Virus</i>
<i>Influenza A H1</i>	<i>Chikungunya Virus</i>
<i>Influenza A H3</i>	<i>SARS-CoV-2</i>

Various samples positive for pathogenic organisms' DNA were tested. No amplification signal was observed in the presence of the following micro-organisms:

- *Klebsiella pneumoniae*
- *Escherichia coli*
- *Pseudomonas aeruginosa*
- *Enterococcus faecalis*
- *Proteus mirabilis*
- *Citrobacter koseri*
- *Staphylococcus aureus*
- *Klebsiella oxytoca*
- *Staphylococcus saprophyticus*
- *Enterobacter aerogenes*
- *Candida Albicans*
- *Streptococcus B*

#### Interferences

The presence of PCR inhibitors in the sample may render the test result uninterpretable. The presence of various inhibitors in plasma, serum or urine samples was tested. Among the substances tested, only blood inhibited the amplification of target sequences with the BIOSYNEX AMPLIQUICK Leptospira kit.

Sample type	Substance tested	Test condition	Observed interference
Serum, Plasma	Blood	≥ 2%	Yes
	Heparin	N/A	None
	EDTA	N/A	
Urine	Albumin	5%	
		0.1%	
	Glucose	1%	
		0,1%	
	pH	pH4	
		pH9	
	Urea	5 mmol/L	
		1.4 mmol/L	
		0.4 mmol/L	
	Blood	2%	
5%			

#### Cut-off

No cut-off needs to be applied.

#### Clinical performances

The clinical performances were determined on 193 samples from 4 internal studies that were qualified positive or negative for *Leptospira* using CE marked PCR tests referenced by the French Health Ministry or bacterial culture paired with MALDI-TOF mass spectrometry.

The clinical performances were determined on 182 samples that were qualified positive or negative using CE marked PCR tests referenced by the French Health Ministry.

	Sample type	Leptospira qualification	
		Positive	Negative
Study 1	Urine	0	90
Study 2	Serum	0	30
	Plasma	0	30
Study 3	Urine	15	2*
	Serum	12	0
	Plasma	3	0
Study 4	Urine	6	0
	Serum	3	0
	Whole blood (only tested on DNA extract)	2	0

\*Among the 32 samples qualified as positive for *Leptospira* in Study 3, two were found negative using the BIOSYNEX AMPLIQUICK *Leptospira* kit and were confirmed negative using another CE marked kit on the market. Thus, those two samples were requalified as negative.

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK *Leptospira* kit on samples after DNA extraction/purification with a commercial kit:

		Qualification	
		Positive	Negative
BIOSYNEX AMPLIQUICK <i>Leptospira</i>	Positive	40	0
	Negative	0	150

**Sensitivity : 100%**      **Specificity : 100%**  
 (95%CI: 91.19% to 100.00%)      (95%CI: 97.60% to 100.00%)

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK *Leptospira* kit on samples treated with the BIOSYNEX AMPLIQUICK Sample Treatment kit:

		Qualification	
		Positive	Negative
BIOSYNEX AMPLIQUICK <i>Leptospira</i>	Positive	35	0
	Negative	4*	152

**Sensitivity : 89.74%**      **Specificity : 100%**  
 (95%CI: 75.78% to 97.13%)      (95%CI: 97.60% to 100.00%)

\*Among the 21 urine samples qualified as positive for *Leptospira*, 4 were found negative with the sample preparation method without DNA extraction. These samples had high Ct values associated and therefore low bacterial loads. LoD95 is much higher with urine sample than with other sample types.

### Precision

The precision data in the context of the BIOSYNEX AMPLIQUICK *Leptospira* kit was determined based on 5 conditions:

- Intra-assay variability (within a same experiment)
- Inter-assay variability (between different experiments)
- Inter-operator variability
- Inter-lot variability
- Inter-laboratory variability

The variability data are expressed in terms of mean value, standard deviation and coefficient of variation, based on the threshold cycle of quantification (Cq) values of the *Leptospira* DNA.

**Intra-assay variability data:**

Two sample dilutions were tested for each matrix: 1 high (+++) corresponding to 200 000 copies/μL and 1 low (+) corresponding to 5 copies/μL, as well as a negative sample (-). Each sample was tested 20 times

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**Samples pre-treated with the BIOSYNEX AMPLIQUICK Sample Treatment kit:**

GENE	Urine			Plasma			Serum		
	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %
Sample -	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Sample +++	27,6	0,4	1,5	27,5	0,3	1,0	27	0,3	1,3
Sample +	29,9	1,0	3,4	29,5	0,9	3,1	28,8	0,7	2,5

Internal control CIEZ	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %
Sample -	22,6	0,6	2,8	23,1	0,5	2,2	23,4	0,5	1,9
Sample +++	22,1	0,5	2,4	22,6	0,8	3,5	22,4	0,3	1,3
Sample +	21,6	0,5	2,4	21,8	0,1	0,6	21,6	0,2	1,0
Procedural control RNase P	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %	Average Cq value	Standard deviation	Coefficient of variation %
Sample -	19,9	0,1	0,7	17,2	0,2	0,9	18,0	0,2	1,0
Sample +++	18,3	0,1	0,7	18,0	0,1	0,8	15,3	0,1	0,7
Sample +	22,2	0,2	1,0	21,9	0,2	0,9	19,2	0,1	0,5

**Inter-assay, operators, and lots variability:**

Two sample dilutions were tested: 1 high (+++) corresponding to 200 000 copies/μL and 1 low (+) corresponding to 5 copies/μL. Negative sample give expected results.

**Inter-days and operators variability**

Each sample is tested on the same lot in duplicate twice each day by two technicians, during 5 days.

**Inter-lots variability**

Each sample is tested in triplicate on the first lot and in duplicate on the second and third lots.

**Inter-places variability**

Each sample is tested in duplicate in the first place and in duplicate in the second place.

GENE LipL32	Average Cq CTRL+	Standard deviation	Coefficient of variation %	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
Inter assay variability	24,3	1,0	4,3	11,4	0,2	2,1	29,5	0,7	2,4
Inter operator variability	24,3	0,0	0,0	11,2	0,2	2,1	29,2	0,5	1,7
Inter lot variability	24,2	0,7	2,9	10,9	0,5	4,4	29,5	1,0	3,3
Inter places variability	24,1	0,4	1,8	10,6	0,3	3,1	26,8	0,9	3,4
Internal control CIEZ	Average Cq CTRL+	Standard deviation	Coefficient of variation %	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
Inter assay variability	22,5	1,1	4,9	23,3	0,4	1,9	23,3	0,1	0,5
Inter operator variability	22,8	0,4	1,6	23,1	0,3	1,5	23,1	0,3	1,5
Inter lot variability	23,6	0,2	0,9	22,7	0,8	3,7	23,5	0,3	1,3
Inter places variability	25,9	0,1	0,5	25,9	0,5	2,1	25,7	0,7	2,7

Procedural control RNase P	Average Cq CTRL+	Standard deviation	Coefficient of variation %	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
Inter assay variability	21,4	0,6	2,9	22,2	0,3	1,5	25,0	0,3	1,2
Inter operator variability	21,8	0,4	2,0	21,9	0,5	2,2	25,0	0,1	0,2
Inter lot variability	22,7	0,2	0,9	21,8	0,3	1,4	24,0	0,4	1,9
Inter places variability	22,3	0,4	1,9	21,3	0,7	3,5	21,9	0,0	0,0

## 12. LIMITS

1. Failure to follow any instruction of the IFU may adversely affect the test performance and/or invalidate the test result.
2. As for any diagnostic test, results must be considered with other clinical information available to the physician. A negative result can never rule out the presence of *Leptospira* in the sample (e.g. genes may be present at a concentration below the minimum detection limit of the test, interference) and a positive result can never ensure the presence of *Leptospira* in the sample (e.g. contamination). A definitive diagnosis can only be made by the physician after evaluation of all clinical and laboratory data.
3. Urine sample with more than 2% blood may interfere with the test result. In case of blood in the urine sample, consider the result with caution. Serum and plasma samples should not contain blood.

## 13. BIBLIOGRAPHY

1. [www.pasteur.fr/fr/centre-medical/fiches-maladies/leptospirose](http://www.pasteur.fr/fr/centre-medical/fiches-maladies/leptospirose)
2. Miotto BA, da Hora AS, Taniwaki SA, Brandão PE, Heinemann MB, Hagiwara MK. Development and validation of a modified TaqMan based real-time PCR assay targeting the lipL32 gene for detection of pathogenic *Leptospira* in canine urine samples. *Braz J Microbiol.* 49:3:584-590. July 2018.
3. Woods K, Nic-Fhogartaigh C, Arnold C, Boutthasavong L, Phuklia W et al. A comparison of two molecular methods for diagnosing leptospirosis from three different sample types in patients presenting with fever in Laos. *Clin Microbiol Infect.* 24:9:1017.1-1017-7. Sept 2018.
4. O Munch M, Chambers LC, Manhart LE, Domogala D, Lopez A, Fredricks DN, Srinivasan S. Optimizing bacterial DNA extraction in urine. *PLOS ONE.* 14:9. Sept 2019.
5. Lucchesi PMA, Arroyo GH, Etcheverria AI, Parma AE, Seijo AC. Recommendations for the detection of *Leptospira* in urine by PCR. *37:2:131-134.* May 2004.
6. Stoddard RA, Gee JE, Wilkins PP, McCaustland K, Hoffmaster AR. Detection of pathogenic *Leptospira* spp. through TaqMan polymerase chain reaction targeting the LipL32 gene. *Diagn Microbiol Infect Dis.* 64:3:247-255. July 2009.
7. Cyrille Goarant. Leptospirosis: risk factors and management challenges in developing countries. *Res Rep Trop Med.* 2016 Sep 28;7:49-62. doi: 10.2147/RRTM.S102543. PMID: 30050339; PMCID: PMC6028063..

## 14. INTERNATIONAL SYMBOLS

	Consult instructions for use or electronic instructions for use		Contains sufficient for <n> tests		Catalog number
	In vitro diagnostic medical device		Temperature limit		Batch code
	Qualitative detection by PCR		Keep away from sunlight		Use-by date
	Negative control		Positive control		Master MIX
	Do not use if package is damaged and consult instructions for use		Manufacturer		Authorized Representative in Switzerland
	Importer				

## 15. MANUFACTURER INFORMATION

**BIO SYNEX S.A.**  
 22 boulevard Sébastien Brant  
 67400 ILLKIRCH-  
 GRAFFENSTADEN – France  
 Standard:  
 Tel: +33 3 88 78 78 87  
[www.biosynex.com](http://www.biosynex.com)  
 Contacts France:  
 Phone : +33 3 88 77 57 00  
[service.clients@biosynex.com](mailto:service.clients@biosynex.com)  
 Contacts other countries:  
 Phone: +33 3 88 77 57 52  
[sales@biosynex.com](mailto:sales@biosynex.com)  
 Customer service  
 Tel. : +33 3 88 77 57 25  
[tech.support@biosynex.com](mailto:tech.support@biosynex.com)

**BIO SYNEX SWISS S.A.**  
 Route de Rossemaison 100  
 2800 DELEMONT - Switzerland

**Latest changes:** IVDR (EU) 2017/746 compliance