BIOSYNEX

BIOSYNEX AMPLIFLASH MALARIA

REF 3150053L

FR (FRANÇAIS / FRENCH)	1
EN (ENGLISH)	12
ES (SPANISH)	22
IT (ITALIAN)	33

CONTENTS

MANUFACTURER INFORMATION

BIOSYNEX S.A.
22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN – France
Standard : +33 3 88 78 78 87
www.biosynex.com

Contact France : +33 3 88 77 57 00
service.clients@biosynex.com
Contact other countries : +33 3 88 77 57 52
sales@biosynex.com
After-sales service: +33 3 88 77 57 25
tech.support@biosynex.com

CH REP

BIOSYNEX SWISS S.A.
Route de Rossemaison 100
2800 DELEMONT - Switzerland

NON-STANDARD SYMBOLS

	FR	EN	ES	IT
	Détection qualitative par PCR	Qualitative detection by PCR	Detección cualitativa mediante PCR	Rilevazione qualitativa mediante PCR
	Tampon d'activation	Activation buffer	Tampón de activación	Tampone di attivazione
	Master mix	Master mix	Master mix	Master mix
	Barrettes de bouchons	Cap strips	Tapones	Tappi

CONTENTS

BIOSYNEX AMPLIFLASH MALARIA

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13. LITERATURE

1 : ECDC Surveillance Report: Malaria, Annual Epidemiology Report for 2021. <https://www.ecdc.europa.eu/en/publications-data/malaria-annual-epidemiological-report-2021>

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





BIOSYNEX AMPLIFLASH MALARIA

REAL-TIME ISOTHERMAL PCR FOR THE QUALITATIVE DETECTION OF PLASMODIUM SPP AND SPECIFIC IDENTIFICATION OF PLASMODIUM FALCIPARUM FROM EXTRACTED WHOLE BLOOD SAMPLE.

In vitro diagnostic medical device for professional use.

REF 3150053L

1. INTENDED PURPOSE

BIOSYNEX AMPLIFLASH Malaria is a molecular *in vitro* diagnostic test for the qualitative detection by real-time isothermal PCR of 5 *Plasmodium* species - *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale* and *Plasmodium knowlesi* (PAN) and the specific identification of *Plasmodium falciparum* (FALCI). The test is an aid for the diagnosis of Malaria performed using a DNA extract obtained from whole blood sample. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

2. CLINICAL SUMMARY

Malaria is a life-threatening disease mostly spreads to people through the bites of some infected female *Anopheles* mosquitoes. Blood transfusion and contaminated needles may also transmit malaria. It is mostly endemic in tropical countries but can be imported through travel and migrations. Malaria can be prevented by avoiding mosquito bites and with adapted medical treatment.

The most common early symptoms of malaria are fever, headache and chills. Symptoms usually start within 10–15 days of getting bitten by an infected mosquito.

Symptoms may be mild for some people, especially for those who have had a malaria infection before. Because some malaria symptoms are not specific, getting tested early is crucial.

Some types of malaria can cause severe illness and death. Infants, children under 5 years, pregnant women, travelers and people with AIDS or immunosuppressed in general are at higher risk. Severe symptoms include: extreme tiredness and fatigue, impaired consciousness, multiple convulsions, difficult breathing, dark or bloody urine, jaundice (yellowing of the eyes and skin), abnormal bleeding, death in extreme cases.

People with severe symptoms should get emergency care right away. Getting treatment early for mild malaria can stop the infection from becoming severe.

Malaria infection during pregnancy can also cause premature delivery or delivery of a baby with low birth weight.

The ECDC Annual Epidemiology Report for 2021¹ gives European Epidemiology data in the European Union.

Of 4 482 confirmed cases for which the pathogen was reported, 3 771 (84.1%) had *P. falciparum*, 96 (2.1%) had *P. vivax*, 181 (4.0%) had *P. ovale*, 137 (3.1%) had *P. malariae*, one had *P. knowlesi*, and 296 (6.6%) cases had infections with unspecified *Plasmodium* species. The case fatality was 1.1% among all 1 944 malaria cases with a known outcome. Among 4 257 cases with known importation status, 99.7% were travel related.

Nearly all malaria cases reported by EU/EEA countries for 2021 were imported. The countries reporting the highest numbers of cases have historical, economic, linguistic, and cultural links with endemic areas, particularly in Africa and the Americas. Most of the imported malaria cases in France and Spain are linked to travel routes from West Africa. Seasonality and age distribution of cases in Europe most likely reflect travel patterns to malaria-endemic countries. The literature suggests that a substantial proportion of imported malaria cases in the EU/EEA occur among recent immigrants from malaria-endemic countries and more settled migrants who have travelled to visit friends and relatives in malaria-endemic countries.

In Europe, malaria chemoprophylaxis is only recommended for travelers to malaria-endemic countries, which are classified into three risk groups (very limited risk of malaria transmission, risk of non-*P. falciparum* malaria and risk of *P. falciparum* malaria) to determine the most effective drug regimen (see WHO's requirements and recommendations¹ for international travelers, including a list of affected countries, as of 2022).

3. TEST PRINCIPLE

BIOSYNEX AMPLIFLASH Malaria kit contains reagents allowing the real-time isothermal LAMP (Loop-Mediated Amplification) amplification to determine positive or negative parasitemia for *Plasmodium spp* (PAN) or *Plasmodium falciparum* specifically (FALCI). The targeted DNA sequence of the PAN version is a conserved sequence of *cox* gene from *Plasmodium spp* genomes. The targeted DNA sequence of the FALCI version is a specific sequence of *cox* gene from *Plasmodium falciparum* genome. The kit contains a well with an endogenous Internal Control (RNAse P encoding gene) to ensure the correct running of the extraction and amplification and assess the quality of sample.

LAMP is a method of specific isothermal amplification of a nucleic sequence with incorporation of fluorescent nucleic acid stain during amplification. It amplifies a target sequence with three to four sets of primers at a constant temperature of +65°C.

BIOSYNEX AMPLIFLASH Malaria allows the test of one sample at a time.

An analysis with the BIOSYNEX AMPLIFLASH Malaria consists of carrying out the following steps:

- Quick DNA extraction from blood with BIOSYNEX Blood Lysis kit;
- Isothermal amplification (LAMP) with the reaction buffer and the primers distributed in the microtubes (assay tube), as shown on below picture.


Figure 1: 1= PAN; 2= Falci; 3= Internal Control; 4= Empty (red bead)

The test allows the amplification of 2 targeted sequences and one internal control:

- PAN: target conserved in *Plasmodium spp* – position 1 on the strip of microtubes (see Figure 1);
- FALCI: target specific of *Plasmodium falciparum* – position 2 on the strip of microtubes (see Figure 1);
- IC: human RNase P encoded gene used as internal control for the lysis and amplification processes – position 3 on the strip of microtubes (see Figure 1).

The tongue (the black strip in the figure) and the well containing the red bead act as polarizers, ensuring the correct orientation of the strip. For correct orientation, the tongue must be positioned to the left and the well containing the red bead to the right when the samples are deposited and when the strip is inserted into the reader.

The kit contains an activation buffer (ratio 10µL sample for 30µL activation buffer) containing all the salts allowing the activation of the reaction at +65°C.

These amplifications are performed in parallel, and the amplification result is obtained after a period of between 3 and 25 minutes depending on the parasitemia of the sample being analyzed.

The isothermal amplification test is based on the detection, by the BIOSYNEX AMPLIFLASH Reader, of fluorescence generated during amplification of the target DNA. The BIOSYNEX AMPLIFLASH Reader device indicates in real time for each tube whether the amplification is positive (detection of the DNA target) or negative (no detection of the DNA target).

4. REQUIRED MATERIAL

Material provided:

- 2 plastic bags each containing 12 strips of 4 wells (1 well containing 1 red bead and 3 containing the master mixes)
- 1 plastic bag containing 24 strips of 4 caps
- 24 activation buffer tubes
- 1 tube of positive control (CONTROL+, red cap)
- 1 tube of negative control (CONTROL-, green cap)

Note: 5 control reactions are possible with the tubes provided.

Material required but not provided:

- Sample pre-treatment kit BIOSYNEX Blood Lysis (Ref. 3150077)
- Powder-free single-use gloves
- Micropipette & filter-tips
- Bench top vortex
- BIOSYNEX AMPLIFLASH Reader (Ref. 5150001)

5. PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of many samples and the high sensitivity of the LAMP technique can, in the absence of precaution, generate false positive results by contamination. If possible, pre-handling LAMP, post-LAMP and DNA extraction should therefore be performed in separate rooms. Wear disposable gloves in each zone and change them before moving from one zone to another.
- The test and caps are for single use only. Do not reuse. Do not open the tubes at the end of the test.
- Do not use the test if the tongue is opened or damaged.
- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.
- Protect the kit from light.

- Centrifuge the tubes before opening, opening them one after the other by closing them well between each one to avoid any contamination.
- The CTRL+ contains significant amounts of DNA sequences. It can therefore potentially contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store this component outside the kit at the first opening of the kit.
- Negative and positive controls are included in the kit. The negative test serves as a sign of possible contamination. The positive control is used to verify the validity of the reagents and the reader. Each laboratory shall define its own quality control system.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste. Consider your local regulation on biological waste disposal.
- The device contains material of bacterial or animal origin and should be handled with caution.
- If, in relation to the use of BIOSYNEX AMPLIFLASH Malaria, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed.

6. RECOMMENDATIONS FOR USE

In view of the significant risk of contamination of the equipment and ambient air associated with isothermal amplification, it is important to carefully observe the recommendations and precautions indicated in the instructions for use.

- Always make sure that the caps of the reaction tubes are tightly closed before starting a heating protocol in the incubator used.
- At the end of the program, when the strips are removed from the reader, place your fingers over the caps to prevent them from opening during removal. Immediately take off gloves used to remove the used strips from the reader by wrapping the microtubes in them before disposing of it in the appropriate waste bin.
- Clean the equipment and the workstation daily with the appropriate detergents and clean monthly with isopropanol according to the protocol provided by Biosynex (refer to the user manual for the BIOSYNEX AMPLIFLASH Reader or the instrument used).

7. KIT STORAGE

The kit is shipped frozen. Kit components must arrive frozen. Upon receipt, keep the kit at a temperature below or equal to -20°C. Under these conditions, the reagents are stable until the expiration date. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light. Do not thaw the kit components except when needed and only before usage. Control tubes can be stored at +2/8°C until the expiration date.

Note: the positive and negative control tubes can undergo 5 freeze – thaw cycles.

8. SAMPLE COLLECTION AND STORAGE

DNA of whole blood samples collected in EDTA or heparin collection tubes extracted with BIOSYNEX Blood Lysis kit should be used with the BIOSYNEX AMPLIFLASH Malaria kit.

Whole blood samples can either be extracted with BIOSYNEX Blood Lysis immediately or within 7 hours if stored at room temperature (+21°C - +30°C) or stored at +4°C for 5 days.

Repeated freeze/thaw cycles of the sample should be avoided (maximum of 1 cycle is allowed).

Samples treated with the BIOSYNEX Blood Lysis kit (step 1 of the test procedure described in section 9) can either be tested immediately or within 7 days if stored between +4°C and +30°C or stored between -20°C and -80°C for 3 months.

Samples mixed with the activation buffer (step 2 of the test procedure described in section 9) are stable for 2 hours at room temperature (+21°C - +30°C). After this time, restart the test procedure.

Transportation of clinical specimens must comply with local regulations for the transport of infectious agents.

9. TEST PROCEDURE

For a test or control, take a tube of activation buffer, a test strip and caps. Use of the instrument to perform a test is described in appendix in paragraph 15 of these instructions for use.

1. Sample lysis and DNA release


For rapid lysis of samples with the BIOSYNEX Blood Lysis kit, please refer to the package insert for this kit, reference 3150077. Positive and negative controls are ready to use and do not need any previous extraction step.

2. Preparation of the reaction mixture

- a. Open the tube of activation buffer.
- b. Using a micropipette fitted with a filter tip, place 10µL of lysate obtained in the previous step in the 30µL of activation buffer. Homogenize by pipetting up and down or by vortexing the tube.

3. Addition of the reaction mixture to the strip

- a. Take the 4 microtubes of the strip without separating them and let them thaw.
- b. Orient the wells of the strip according to the position of the polarizers (tongue to the left and well containing the red bead to the right); well #1 is the proximal well.
- c. Remove the aluminum foil and discard it.

Microtubes from the left to the right: 1= PAN; 2= FALCI; 3= Internal control; 4= red bead (Do not add anything to the well)

- d. Pipet 10µL of reaction mixture with a micropipette and filtered tip and dispense the 10µL into each microtube of the strip.
 - e. Homogenize by pipetting.
 - f. Close the wells with the provided caps.
 - g. Hold the strip vertically and tap it on a flat surface to remove any air bubbles.
 - h. Without delay, place the strip in the AMPLIFLASH Reader in the correct direction.
- The well containing the red bead allows to maintain the correct orientation of the strip.

Note 1: The controls provided allow 5 reactions to be performed.

10. DATA ANALYSIS AND INTERPRETATION OF RESULTS

A. Study validation criteria

A positive control and a negative control are supplied in the kit. Each laboratory must define the frequency of the quality controls it applies. The results expected with each control are defined in the following table:

Control	Target		
	PAN	FALCI	RNaseP
Positive control	+	+	+
Negative control	-	-	-

B. Qualitative interpretation (positive or negative)

The interpretation of blood samples is given in the following table:

Target			Interpretation
PAN	FALCI	RNaseP	
+	+	+	Patient with <i>Plasmodium falciparum</i> specific DNA and may also have co-infection status
+	+	-	
+	-	+	Patient with <i>Plasmodium spp</i> DNA other than <i>P. falciparum</i> DNA
+	-	-	
-	-	+	Patient without specific DNA of any of the 5 detectable <i>Plasmodium spp</i> species
-	-	-	Invalid result: poor sample quality, inhibition of PCR reaction, or problems during run or test process - perform a new test or re-sampling

+ Positive results - Negative results

Note: When using the AMPLIFLASH Reader, wells in positions 1, 2 and 3 can appear orange in the Gene Series Ampliflash software or with the "!" pictogram on the reader. In this case, the test is invalid. Repeat the test with a new sample. Positions 5 to 8 of the reader do not accommodate tubes, so an invalid result corresponds to normal functioning.

11. PERFORMANCES

Analytical sensitivity

The analytical sensitivity of the kit's isothermal amplification reaction was determined from quantified DNA from the 1st WHO International Standard for *Plasmodium falciparum* and quantified genomic DNA of *Plasmodium falciparum* form ATCC (American Type Culture Collection). The limit of detection (LoD) of the BIOSYNEX AMPLIFLASH Malaria kit is defined as the concentration, in copy number/µL, that can be detected at 95% at least, on a specific *Plasmodium spp* DNA sample.

Each dilution was tested 24 times. The statistically determined results are presented in the following table:

Target	LoD ₉₅	
	<i>P. falciparum</i> WHO Standard	<i>P. falciparum</i> ATCC gDNA
FALCI	3.09 IU/reaction (0.127 IU/µL)	2.75 copies/reaction (1.1 copies/µL)
PAN	3.19 IU/reaction (0.123 IU/µL)	2.3 copies/reaction (0.92 copy/µL)

Analytical specificity

The design of oligonucleotides was validated *in silico* by BLAST alignment. The comparison of the sequences obtained shows a specific detection of the BIOSYNEX AMPLIFLASH Malaria on its *Plasmodium spp* target.

No primers detect DNA other than the targeted *Plasmodium* spp.

Cross-reactivity

A panel of 99 DNA samples and 42 RNA samples from a biobank listed in the following tables were tested with the BIOSYNEX AMPLIFLASH Malaria kit. For all these samples, no amplification of the targets (PAN or FALCI) was observed.

DNA	
<i>Acanthamoeba castellanii</i>	<i>Helicobacter pylori</i>
Adenovirus 41	Herpès simplex 1
Adenovirus	Herpès simplex 2
<i>Aspergillus fumigatus</i>	HHV-6
<i>Bacillus cereus</i>	HHV-8
<i>Bartonella henselae</i>	<i>Klebsiella pneumoniae (ndm-1)</i>
<i>Bartonella quintana</i>	<i>Legionella pneumophila</i>
Virus BK	<i>Leishmania chagasi</i>
<i>Bordetella holmesii</i>	<i>Leishmania infantum</i>
<i>Borrelia afzelii</i>	<i>Listeria monocytogenes</i>
<i>Borrelia burgdorferi</i>	<i>Moraxella catarrhalis</i>
<i>Borrelia garinii</i>	<i>Mycobacterium avium</i>
<i>Brucella abortus</i>	<i>Mycobacterium intracellulare</i>
<i>Campylobacter jejuni</i>	<i>Mycobacterium kansasii</i>
<i>Candida albicans</i>	<i>Mycobacterium tuberculosis</i>
<i>Candida auris</i>	<i>Mycobacterium ulcerans</i>
<i>Chlamydia trachomatis</i>	<i>Mycoplasma genitalium</i>
<i>Chlamydophila pneumoniae</i>	<i>Mycoplasma hominis</i>
<i>Chlamydophila psittaci</i>	<i>Neisseria gonorrhoeae</i>
<i>Clostridium difficile</i>	<i>Neisseria meningitidis sg a</i>
<i>Coccidioides immitis</i>	<i>Neisseria meningitidis sg b</i>
<i>Coxiella burnetii</i>	Papillomavirus type 16
<i>Cryptosporidium parvum</i>	Papillomavirus type 18
Cytomegalovirus	Parvovirus b19 (plasmide)
<i>Enterococcus faecalis (vanb)</i>	<i>Vibrio cholerae</i>
<i>Enterococcus faecium (vana)</i>	<i>Yersinia enterocolitica</i>
Epstein-barr virus	<i>Aeromonas caviae</i>
<i>Escherichia coli (eaec)</i>	<i>Aeromonas dhakensis</i>
<i>Escherichia coli (eiec)</i>	<i>Vibrio vulnificus</i>
<i>Escherichia coli (epec)</i>	<i>Aeromonas veronii</i>
<i>Escherichia coli (etec)</i>	<i>Aeromonas hydrophila</i>
<i>Escherichia coli (vtcc)</i>	<i>Shigella sonnei</i>
<i>Rickettsia conorii</i>	<i>Shigella boydii</i>
<i>Salmonella enteritidis</i>	<i>Shigella flexneri</i>
<i>Salmonella typhi</i>	<i>Shigella dysenteriae</i>
<i>Shigella flexneri</i>	<i>Campylobacter upsaliensis</i>
<i>Staphylococcus aureus (mecA+)</i>	<i>Plesiomonas shigelloides</i>
<i>Staphylococcus aureus (mecA-)</i>	<i>Campylobacter coli</i>
<i>Streptococcus pneumoniae</i>	<i>Escherichia coli o111/nm</i>
<i>Toxoplasma gondii</i>	<i>Campylobacter fetus</i>
<i>Treponema pallidum</i>	<i>Yersinia pseudotuberculosis</i>
<i>Trichomonas vaginalis</i>	Amplirun adenovirus 41
<i>Trypanosoma cruzi</i>	Adenovirus type 40 (species f) (strain: dugan)
<i>Ureaplasma parvum</i>	Adenovirus type 04 (species e)
<i>Ureaplasma urealyticum</i>	Adenovirus type 20 (species d)
<i>Varicella-zoster virus</i>	Adenovirus type 11 (species b)
<i>Francisella tularensis</i>	Adenovirus type 02 (species c)
<i>Gardnerella vaginalis</i>	Adenovirus type 31 (species a)
<i>Haemophilus ducreyi</i>	-
<i>Haemophilus influenzae</i>	-

RNA		
Coronavirus oc43	Rotavirus	Parainfluenza 2
Coronavirus 229-e	Rubella	Parainfluenza 3
Coronavirus sars (2003)	Sars-cov-2	Parainfluenza 4 a
Coronavirus nl63	Tick-borne encephalitis virus	Respiratory syncytial virus (subtype a)
Coxsackie a6	Human parainfluenza 1	Respiratory syncytial virus (subtype b)
Coxsackie b1	Influenza a h1	West Nile virus
Coxsackie b5	Influenza a h3	Yellow fever virus
Dengue 1 virus	Influenza a h5	Zika virus (asian lineage)
Dengue 2 virus	Influenza b	Zika virus
Dengue 3 virus	Measles	Chikungunya virus
Dengue 4 virus	Mers coronavirus	Norovirus g1
Echovirus 5	Mumps	Norovirus g2
Enterovirus 68	Novel influenza a h1n1	Sapovirus
Rhinovirus	Parainfluenza 1	Astrovirus

• **Interference study**

The presence of PCR inhibitors in the sample may induce a positive or negative interference on test results. Whole blood samples, who are usually collected in classical EDTA tubes (EDTA blood samples used to develop and validate the BIOSYNEX Blood Lysis = prototype), can, in some cases, be collected in heparin tubes and heparin could be a source of inhibition. Whole blood with abnormal haematological and biochemical parameters could also be a source of inhibition.

The interference tests of Heparin on the BIOSYNEX AMPLIFLASH Malaria kit were performed using negative blood samples spiked or not with positive blood sample of *P. falciparum* from the WHO STD. No positive or negative interference was demonstrated. Whole blood interference tests with abnormal haematological or biochemical parameters on the BIOSYNEX AMPLIFLASH Malaria kit were performed using negative blood samples spiked or not with quantified *P. falciparum* gDNA from the ATCC. No positive or negative interference was demonstrated for the following haematological and biochemical parameters:

- Bilirubin ≤ 158.6 mg/L
- Triglycerides ≤ 8.29 g/L
- Leukocytes ≤ 13,880 cells/m³
- Erythrocytes ≤ 6,300,000 cells/m³

To this date, our clinical studies on samples with blood have not shown any interference on the results.

• **Precision**

Repeatability and Reproducibility

Repeatability data for the BIOSYNEX AMPLIFLASH Malaria kit was determined on replicate tests using positive and negative blood samples, each run 30 times. As no positive blood sample of *P. knowlesi* is available, results for this species were obtained with dilution of purified genomic DNA from BEI Resources. A 100% agreement was shown on each hypothesis.

Results (min:sec)					
Sample		PAN		FALCI	
		Mean	SD	Mean	SD
Negative sample	-	N/A	-	N/A	-
<i>Plasmodium falciparum</i>	+++	03:18	00:10	04:08	00:17
	++	04:49	00:13	06:17	01:18
	+	06:11	00:27	07:33	00:36
<i>Plasmodium malariae</i>	++	07:31	00:20	N/A	-
	+	10:19	00:53	N/A	-
<i>Plasmodium ovale</i>	+++	04:50	00:09	N/A	-
	++	05:24	00:30	N/A	-
	+	07:59	00:37	N/A	-
<i>Plasmodium vivax</i>	++	06:01	00:12	N/A	-
	+	07:56	00:24	N/A	-

<i>Plasmodium knowlesi</i>	+++	05:32	00:17	N/A	-
	++	06:47	00:40	N/A	-
	+	06:53	00:55	N/A	-

Inter-day, inter-batch, inter-laboratory, and inter-operator reproducibility data for the BIOSYNEX AMPLIFLASH Malaria kit were determined on replicate tests from positive and negative blood or genomic DNA samples. A 100% agreement was shown on each hypothesis.

Inter-days, operators, lots and places variability:

Positive and Negative Controls were tested. Two dilutions of gDNA from the ATCC were tested: 1 medium (++) at 100 copies/μL and 1 low (+) at 2 copies/μL) and two dilutions of the WHO Standard : 1 high (+++) at 25000 IU/reaction and 1 low (+) at 10 IU/reaction). Negative sample was also tested.

Inter-days variability

Each sample is tested on the same lot in duplicate each day by the same technician, during 5 days.

Inter-operators variability

Each sample is tested on the same lot in duplicate by two technicians.

Inter-lots variability

Each sample is tested in duplicate on three lots.

Inter-places variability

Each sample is tested in duplicate in the first place and in duplicate in the second place.

PAN	CTRL+		CTRL-		ATCC gDNA ++		ATCC gDNA +		WHO Standard +++		WHO Standard +		Negative Sample	
	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD
Inter days variability	03:59	00:20	N/A	-	03:56	00:04	05:32	00:20	04:00	00:32	06:46	02:15	N/A	-
Inter operators variability	03:42	00:10	N/A	-	04:02	00:10	07:25	01:24	04:05	00:42	06:22	00:17	N/A	-
Inter lots variability	04:15	00:18	N/A	-	04:11	00:24	06:38	01:22	03:21	00:12	06:01	00:02	N/A	-
Inter places variability	03:55	00:07	N/A	-	04:17	00:24	05:12	00:10	03:35	00:07	05:52	00:24	N/A	-
FALCI	CTRL+		CTRL-		ATCC gDNA ++		ATCC gDNA +		WHO Standard +++		WHO Standard +		Negative Sample	
	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD
Inter days variability	04:35	00:12	N/A	-	04:48	00:05	06:47	00:49	04:34	00:13	08:16	01:17	N/A	-
Inter operators variability	04:30	00:21	N/A	-	04:50	00:00	07:30	01:17	04:42	00:17	07:27	00:17	N/A	-
Inter lots variability	04:50	00:32	N/A	-	04:51	00:58	06:56	00:59	03:46	00:40	07:08	00:59	N/A	-
Inter places variability	04:47	00:17	N/A	-	05:05	00:21	07:15	01:03	04:27	00:31	06:57	00:17	N/A	-
RNaseP - Internal Control	CTRL+		CTRL-		ATCC gDNA ++		ATCC gDNA +		WHO Standard +++		WHO Standard +		Negative Sample	
	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD	Mean	SD
Inter days variability	05:35	00:31	N/A	-	N/A	-	N/A	-	07:48	01:12	07:47	01:48	07:04	01:32
Inter operators variability	05:37	00:10	N/A	-	N/A	-	N/A	-	06:20	00:07	06:17	00:03	06:15	00:07
Inter lots variability	05:36	00:22	N/A	-	N/A	-	N/A	-	06:48	00:45	07:13	00:51	06:43	00:23
Inter places variability	06:05	00:00	N/A	-	N/A	-	N/A	-	07:22	00:45	06:45	00:07	06:32	00:17

• **Clinical performances**

The clinical performance of the BIOSYNEX AMPLIFLASH Malaria kit was evaluated on clinical samples from various studies:

- Study on negative blood bank samples (Study 1)
- Study on positive and negative samples from a hospital center (Study 2)
- Study on positive and negative samples from an African hospital (Study 3)

Total samples from all studies for samples treated with BIOSYNEX Blood Lysis kit:

	Number of Samples	Malaria Qualification	
		Negative Samples	Positive Samples (PAN well)
Study 1	30	30	0
Study 2	47	10	37
Study 3	56	20	36
Total	133	60	73

Samples with distinction between *P. falciparum* and other *Plasmodium* species:

	Number of Samples	Malaria Qualification			
		Negative Samples	<i>P. falciparum</i>	<i>P. malariae</i>	<i>P. ovale</i>
Study 1	30	30	0	0	0
Study 2	47	10	33	2	2
Study 3	55	20	35*	0	0
Total	132	60	68	2	2

*One sample was qualified positive for *Plasmodium spp* and negative for *Plasmodium falciparum*, therefore not considered in this table (unknown qualification against malaria species).

Then, data from all studies are merged.

Results obtained with the 133 blood samples from all the studies are presented below.

The contingency table below shows the overall performance on whole blood samples.

BIOSYNEX AMPLIFLASH Malaria		Sample Qualification	
		Positive	Negative
	Positive	73	0
	Negative	0	60

Sensitivity : 100% (95%CI: 95.07-100%)

Specificity : 100% (95%CI: 94.04-100%)

Precision : 100% (95%CI: 97.26-100%)

(CI: Confidence Interval)

12. LIMITATIONS

1. Failure to follow any instructions of the IFU may adversely affect the test performance and/or invalidate the test result.
2. The possibility to obtain false positive or false negative results cannot be ruled out (e.g. in case of accidental contamination, in case of poor quality of the sample or in case the test procedure is not correctly followed). As with all diagnostic tests, results must be considered with other clinical information available to the physician.

13. LITERATURE

1-ECDC Surveillance Report: Malaria, Annual Epidemiology Report for 2021. <https://www.ecdc.europa.eu/en/publications-data/malaria-annual-epidemiological-report-2021>.

14. HISTORY

Record version	Amended paragraph	Change details
V1	All	Creation of the document.


15. ANNEX §9 TEST PROTOCOL

Alternative procedure without the use of the computer.


1. Unplug the USB connection between the AMPLIFLASH Reader instrument and the computer and start the instrument using the ON/OFF button on the left side. The screen opens to the list of protocols stored on the instrument.


2. Select the "AMPLIFLASH Malaria" protocol on the touch screen. The "New protocol" window opens.



3. Manually enter the patient ID on the "Name" line and press "OK".

4. The list of targets detected by the test is displayed. Select "Run" at the top right of the screen.

5. A list of instructions appears on the screen corresponding to "Open the cover", "Insert the strip into the machine", "Close the cover", and "Click OK". Insert the test strip in the correct orientation, close the cover and select "OK".

6. The results can be observed in real time in the "Amp." and "Results" tabs, allowing you to see in real time the appearance of the amplification curves and the detection or not of the targets in the wells as well as the detection time of the positive signal.
7. At the end of the run, the results can be retrieved in .csv format using a USB stick inserted into the port on the left side of the instrument and selecting "Export" on the bottom right of the screen. It is also possible to retrieve old tests, including those carried out via the control software by going to the "Data" tab at the bottom of the home screen. All the tests done on the instrument are recorded there and can also be retrieved on a USB key according to the same procedure. Note that exporting test results from the screen, one by one, enables them to be retrieved in spreadsheet format, whereas exporting several results together from the internal memory will be in "text" format.

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