BIOSYNEX AMPLIQUICK Bordetella

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BIOSYNEX AMPLIQUICK Bordetella

PCR FOR THE QUALITATIVE DETECTION AND DIFFERENTIATION OF BORDETELLA PERTUSSIS, BORDETELLA PARAPERTUSSIS, BORDETELLA BRONCHISEPTICA AND BORDETELLA HOLMESII IN NASOPHARYNGEAL SWABS OR NASAL ASPIRATES.

In vitro diagnostic medical device for professional use.

REF 3150068_SEC01 / 3150068_SEC02

1 | INTENDED PURPOSE

BIOSYNEX AMPLIQUICK Bordetella is a molecular *in vitro* diagnostic test for the qualitative detection by qPCR and differentiation of 4 bacterial species of the genus *Bordetella* encountered in human medicine: *B. pertussis*, *B. parapertussis*, *B. bronchiseptica* and *B. holmesii*. The test is an aid for the diagnosis of whooping cough performed using a DNA extract obtained from nasopharyngeal swabs or nasal aspirates. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

2 | CLINICAL SUMMARY

The genus *Bordetella* is composed of 9 different species. The main species responsible for respiratory diseases in humans are:

- *B. pertussis* which is the agent responsible for whooping cough and whose reservoir is exclusively human,
- *B. parapertussis*, a species close to *B. pertussis* but which does not secrete *pertussis* toxin, and which may be responsible for a whooping cough syndrome that is usually less severe,
- *B. bronchiseptica*, affecting a large number of mammals and responsible in humans for respiratory infections mainly in the immunocompromised,
- and *B. holmesii* which can also be responsible for respiratory symptoms similar to whooping cough and even septicemia.

In unvaccinated patients or those whose whooping cough vaccination status is unknown, presenting with a persistent cough, *Bordetella* testing should be performed as soon as possible to initiate antibiotic treatment and isolate the patient in order to break the chain of infection.


Diagnosis of whooping cough by PCR should be performed within the first three weeks after the onset of symptoms.

3 | TEST PRINCIPLE

The BIOSYNEX AMPLIQUICK Bordetella kit is an *in vitro* diagnostic test based on real-time Polymerase Chain Reaction (PCR) technology by fluorescent probe hydrolysis. It is composed of 96-wells microplates pre-filled with two ready-to-use Master mixes containing dNTPs, MgCl₂, primers and fluorescent probes, Taq polymerase enzyme, and reaction buffer.

The assay consists of a multiplex PCR step carried out in a real-time thermocycler, allowing the specific amplification and simultaneous detection of molecular targets carrying the sequences of interest: the TOX operon (sequence IS481), the pertussis toxin promoter region (ptxA-Pr), the insertion sequence h-IS1001 specific to *B. holmesii*, and the flagellin gene specific to *B. parapertussis* and which allows to differentiate it from *B. bronchiseptica* in two steps (FLA1 and FLA2 sequences), as well as the internal procedural control – (human RNase P gene) which allows to evaluate the quality of the biological sample collection, and to validate the DNA extraction and purification steps.

For this purpose, two Master mixes containing probes labelled with FAM, HEX and Cy5 fluorophores are used.

In the blue colored Master mix contained in the wells marked "1" on the pre-filled plate, the fluorophore FAM is used for the sequence-specific amplification of IS481, HEX for the sequence ptxA-Pr and Cy5 for the RNase P gene. In the clear Master mix contained in the wells marked "2" on the plate, the fluorophore FAM is used for the sequence-specific amplification of IS1001, HEX for the sequence FLA1 and Cy5 for the sequence FLA2. Determination of the Bordetella species depends on the combination of detection of these genes as described in §10.

The increase in fluorescence signal is only detected if the target sequence complementary to the amplified probe is present in the sample. The fluorescent signal is therefore directly proportional to the amplification of the target during the amplification phase. The Cq (cycle quantification) value corresponds to the number of cycles at which the fluorescence starts to increase exponentially in contrast to the background. The kit can be used on two types of samples:

- nasopharyngeal swabs
- nasopharyngeal/endonasal aspirates

Prior to their use with the BIOSYNEX AMPLIQUICK Bordetella amplification kit (§9), samples must be first collected and treated according to the procedures detailed in §7 and §8 using BIOSYNEX AMPLIQUICK Sample Collection kit (ref 3150060 and related references) and BIOSYNEX AMPLIQUICK Lysis kit (ref 3150059 and related references), or processed to obtain purified DNA extracts using a suitable kit.

4 | KIT CONTENT

Material included

- 2 divisible 96-wells microplates prefilled with the Master mixes
- 1 tube of positive control (CONTROL +, red cap)
- 1 tube of negative control (CONTROL -, green cap)
- 2 bags of optical cap strips, for microplate sealing

Material required but not provided

- DNA extraction kit
- Powder-free disposable gloves
- Pipettes and filtered tips
- Real-time PCR Thermal Cycler
- PCR microplate or microtubes centrifuge

The PCR thermal cycler used for the test must have the following main characteristics:

- Open system
- Real-time quantitative PCR assays.
- Programmable Thermal Cycler block

Reference	Thermal cycling block
3150068_SEC01	0.1mL low profile
3150068_SEC02	0.2mL high profile

- Excitation source: Leds, lamp or laser
- Filter sets (excitation/emission wavelengths) suitable for the detection of "reporter" fluorophores of the FAM, HEX and Cy5 probes.
- Connection with a computer using specific analysis software that allow the recovery of fluorescence data, absolute quantification assays, and the interpretation of results.

The kit has been developed and validated with the following real-time PCR thermal cyclers:

Reference	Thermal cycler
3150068_SEC01	CFX96™ Touch Real-Time PCR Detection System (Bio-Rad)
	CFX96™ Opus Real-Time PCR Detection System (Bio-Rad)
	QuantGene 9600 (Bioer)
	QuantStudio 5 System (Applied Biosystems)
	LightCycler 480 (Roche)
	Amplix-DT lite 48 (DNA Technology)

3150068_SEC02	CFX96™ Touch Real-Time PCR Detection System (Bio-Rad) CFX96™ Opus Real-Time PCR Detection System (Bio-Rad) QuantGene 9600 (Bioer) QuantStudio 5 System (Applied Biosystems) Amplix-DT lite 48 (DNA Technology)
---------------	--

If another thermal cycler is used, please perform your own validation of the BIOSYNEX AMPLIQUICK Bordetella kit using the provided controls and/or qualified respiratory samples before using the test.

5 | PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of a large number of samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination. Pre-handling PCR, post-PCR and DNA extraction should therefore be performed in separate rooms. Wear disposable gloves in each zone, and change them before moving from one zone to another.
- The test and cap are for single use only. Do not reuse. Do not open the PCR tubes at the end of the test.
- Do not use the test if aluminum foil is opened or damaged. Once the aluminum foil is removed, use the test immediately.
- Do not use the kit if it arrived unfrozen.
- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.
- Protect the kit from light.
- Centrifuge the tubes before opening, opening them one after the other by closing them well between each one to avoid any contamination.
- The CTRL+ contains significant amounts of DNA sequences. It can therefore potentially contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store this component outside the kit at the first opening of the kit.
- The negative and positive controls included in the kit mimic results obtained with negative or positive samples respectively. They must be used with each new run.
- To ensure that there is no contamination during handling, it is left to the user to include as part of their internal quality approach, at the frequency chosen, an experiment negative control well in which 8µL of molecular biology grade water are added to the Master mix.
- When using negative and positive controls with a series of patients, it is recommended to first deposit the negative control, then deposit the patient samples, and to finish with the deposit of the positive control.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste.
- The device contains material of bacterial or animal origin and may transmit infectious agents and should be handled with extreme caution.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Bordetella, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed.

6 | REAGENT STORAGE AND STABILITY

The kit is shipped frozen. Kit components must arrive frozen. Store the kit at a temperature of -20°C. Under these conditions, the reagents are stable until the expiry date indicated on the kit label. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light.

The positive and negative controls may undergo up to 15 freeze/thaw cycles.

Remove and thaw only the required number of well strips. As the strips are ready to use, there is no reason to subject them to repeated freeze-thaw cycles.

7 | SAMPLE COLLECTION AND STORAGE

BIOSYNEX AMPLIQUICK Bordetella allows the amplification of DNA of all 4 Bordetella species. The test can be performed using purified DNA extract from nasopharyngeal samples (swabs) or nasopharyngeal/endonasal

aspirates; or the same samples collected in the BIOSYNEX AMPLIQUICK Sample Collection (ref 3150060 and related references) and treated with the BIOSYNEX AMPLIQUICK Lysis kit (ref 3150059 and related references).

- Collect nasopharyngeal specimens using sterile swabs discharged into sterile tubes containing transport medium, or those available in the BIOSYNEX AMPLIQUICK Sample Collection kit (ref 3150060 and related references).
- Collect nasopharyngeal/endonasal aspirates with a mucus aspirator according to the manufacturer's recommendations and discharge into sterile tubes containing transport medium, or those available in the BIOSYNEX AMPLIQUICK Sample Collection kit (ref 3150060 and related references). When using the BIOSYNEX AMPLIQUICK Lysis kit (ref. 3150059), aspirated samples must be used directly without being unloaded into a transport environment.

Samples discharged into the BIOSYNEX AMPLIQUICK Sample Collection transport medium can either be extracted immediately or within 4 hours if stored at room temperature or stored between 2°C and 8°C for 24 hours.

Transportation of clinical specimens must comply with local regulations for the transport of infectious agents.

8 | EXTRACTION OF NUCLEIC ACIDS

The nucleic acids extraction must be performed before the amplification protocol using a suitable extraction system for bacterial DNA. Please follow the manufacturer's instructions.

The BIOSYNEX AMPLIQUICK Bordetella kit has been validated with the following extraction kits:

- BIOSYNEX AMPLIQUICK Lysis (ref. 3150059 and related references). **BIOSYNEX AMPLIQUICK Lysis can only be used with nasopharyngeal sample collected in BIOSYNEX AMPLIQUICK Sample Collection. Aspirated samples can be used without unloading into a transport environment.**
- QIAamp DNA Mini kit (Qiagen – Ref 51306)
- NucleoMag Dx Pathogen (Macherey Nagel – 744215.4)
- Amplix Bacterial DNA Extraction kit on Amplix Automated System (DNA Technology - Ref OP05006)

If another DNA extraction method or kit is used, please perform your own method validation of the BIOSYNEX AMPLIQUICK Bordetella kit using qualified respiratory samples before using the test (do not use the controls provided).

It is not necessary to extract the positive and negative controls with the Nucleic Acid Extraction Kit. If use of purified DNA is delayed, before addition to the Master Mix, store at 4°C or on ice on the day of testing.

9 | AMPLIFICATION PROTOCOL

For the plate design, please refer to §3.

It is recommended that the Master mix strips be placed on a refrigerated rack or on ice while samples are being added.

1. Take one breakable plate in strips of 8 wells, and take the number of strips required. If less than a whole strip of 8 wells is needed, the strip can be cut horizontally using a pair of scissors.
2. Centrifuge the strips for a few seconds to collect any droplets on the well edges or on the seal.
3. Carefully remove and discard the aluminum foil. Once the aluminum foil is removed, use the strip immediately.
4. Add 8 µL of sample or control according to the plate design.
5. Close the wells with the transparent caps provided.
6. Centrifuge the strips for a few seconds.
7. Place them in the thermal cycler and run the following amplification program:

PCR Program :

Step	Repetitions	Temperature	Duration	Acquisition
Taq Activation	1x	95°C	3 min	-
Denaturation	50x	95°C	10 sec	-
Hybridation / Elongation		58°C	45 sec	yes

Enter **20 µL** of reaction volume into the Thermal Cycler program.
Please refer to the operating instructions of the Thermal Cycler used for programming information.

Detection channels setting:

Blue Master mix	
Target	Fluorochrome
IS481	FAM
ptxA-Pr	HEX
RNAse P (procedural control)	Cy5
Clear Master mix	
Target	Fluorochrome
IS1001	FAM
FLA1	HEX
FLA2	Cy5

10 | DATA ANALYSIS AND INTERPRETATION OF RESULTS

A. Trial Validation Criteria

Negative control:

The fluorescence emitted must be below the threshold with the exception of the Cy5 channel of the blue Master mix. This is an indicator of non-specific amplification. If the fluorescence is above the threshold, check for an atypical curve. In the case of an amplification curve, a contamination or a distribution error in the microtubes is to be considered. Only the internal control RNAse P should be amplified.

Positive control:

The value of the positive control should preferably be detected before 30 cycles ($C_q \leq 30$). In the absence of amplification of the positive control, the existence of an amplification or fluorescence detection problem (defective thermal cycler or instrument non-adapted to the method) should be considered.

Internal amplification control:

The exogenous internal control RNAse P ensures that the enzymes in the Master mix are functional and validates the nucleic acid collection and extraction steps. Indeed, the amplification curve of the RNAse P procedural control should be observed in the Cy5 channel of the wells containing the blue colored Master mix.

Nevertheless, two situations of lack of amplification of the internal control can be observed:

- If the target genes are initially present in the sample with a high number of copies, the internal control may not be amplified. This result is consistent and does not invalidate the test. It should be interpreted as a positive result despite the lack of signal from the internal control. This phenomenon is the result of amplification competition between the Internal Control and targets present at high copy numbers.
- If the target genes in the FAM and HEX channels are not amplified, as well as the internal control in the Cy5 channel in the blue Master mix, then no result can be rendered. This situation highlights the presence of PCR inhibitors. PCR should be repeated starting from the primary sample and preferably on DNA extract.

B. Qualitative interpretation (positive or negative)

Signals above the threshold, **and visually consistent with a classical PCR amplification curve**, are considered positive results.

Some samples may show atypical curves that are not characteristic of amplification curves. In this case, the result should not be considered as interpretable and the analysis of the sample repeated with the controls.

Detection channels						Interpretation
Blue Master mix 1			Clear Master mix 2			
FAM (IS481)	HEX (ptxA-Pr)	Cy5 (RNaseP)	FAM (IS1001)	HEX (FLA1)	Cy5 (FLA2)	
-	-	+	-	-	-	Negative control
+	+	+	+	+	+	Positive control
+	+	+	-	-	-	Patient with <i>B. pertussis</i> specific DNA
+	+	-	-	-	-	
+	-	+	-	-	-	Patient with <i>B. holmesii</i> specific DNA
+	-	-	-	-	-	
+ If Cq >30	-	+/-	-	-	-	Patient with <i>Bordetella spp</i> DNA
-	-	+	+	+	+	Patient with <i>B. parapertussis</i> specific DNA
-	-	-	+	+	+	
-	-	+/-	+ If Cq >30	-	-	Patient with <i>B. parapertussis</i> or <i>B. bronchiseptica</i> DNA
+/-*	-	+	+/-*	-	+	Patient with <i>B. bronchiseptica</i> specific DNA
+/-*	-	-	+/-*	-	+	
-	-	-	-	-	-	Invalid result, retest

*The result is strain dependent. Indeed, depending on the strain, some species may or may not carry the gene. Thus, the *B. bronchiseptica* strain FR3523 is positive for IS1001 while the reference strain RB50 is negative for this gene.

The Master mix 1 (blue) allows the detection of *B. pertussis* or *B. holmesii* and Master mix 2 (colorless) the detection of *B. parapertussis*. The combination of the expression of the different genes detected in the two wells allows the identification of the other specie *B. bronchiseptica*.

11 | PERFORMANCES

• Analytical sensitivity

The detection limit (LoD) of the BIOSYNEX AMPLIQUICK Bordetella kit is defined as the concentration that can be detected at 95% at least on a specific *Bordetella* DNA sample.

Detection limits on extracted DNA samples per specie:

The LoD₉₅ on extracted DNA is expressed in number of copies of DNA/μL. Quantified DNA extracts were purchased from Vircell (Amplirun®). LoDs were assessed for *Bordetella holmesii* on the master mix 1 as it is positive for the gene IS481, *Bordetella pertussis* on the master mix 1 as it is positive for both targets IS481 and ptxA-Pr and for *Bordetella parapertussis* on the master mix 2 as it is positive for IS1001 and both FLA genes. Quantified DNA of *B. bronchiseptica* is not available on the market, therefore the LoD₉₅ could not be assessed for this specie.

The LoD₉₅ of the BIOSYNEX AMPLIQUICK Bordetella on extracted DNA have been statistically determined for each strain and gene sequence. Results are given in the table below.

Bordetella strain	LoD ₉₅ (copies/μL)
<i>Bordetella pertussis</i>	IS481: 0.043 ptxA-Pr: 4.017
<i>Bordetella parapertussis</i>	IS1001: 0.265 FLA 1: 18.418 FLA 2: 16.033
<i>Bordetella holmesii</i>	IS481: 0.253

Detection limits for samples collected with AMPLIQUICK Sample Collection transport medium and treated with the BIOSYNEX AMPLIQUICK Lysis kit:

The LoDs were determined by performing a series of dilutions of a reference sample with a known concentration CFU/mL in qualified *Bordetella* negative nasopharyngeal swabs discharged into the BIOSYNEX AMPLIQUICK Sample Collection transport medium or negative nasal aspirate samples and treated with the BIOSYNEX AMPLIQUICK Lysis kit.

The LoD₉₅ of the BIOSYNEX AMPLIQUICK *Bordetella* on DNA treated with BIOSYNEX AMPLIQUICK Lysis has been statistically determined for each strain and gene sequence. Results are given in the tables below.

Bordetella strain	LoD ₉₅ for nasopharyngeal swabs				
	IS481	ptxA-Pr	IS1001	FLA Publi	FLA 4
<i>Bordetella pertussis</i>	22.837 CFU/mL	451.291 CFU/mL	-	-	-
<i>Bordetella parapertussis</i>	-	-	38.264 CFU/mL	528.032 CFU/mL	429.806 CFU/mL
<i>Bordetella holmesii</i>	113.354 CFU/mL	-	-	-	-
<i>Bordetella bronchiseptica</i>	-	-	-	-	1791.123 CFU/mL

Bordetella strain	LoD ₉₅ for nasal aspirates				
	IS481	ptxA-Pr	IS1001	FLA Publi	FLA 4
<i>Bordetella pertussis</i>	36.679 CFU/mL	2758.238 CFU/mL	-	-	-
<i>Bordetella parapertussis</i>	-	-	52.0 CFU/mL	1535.5 CFU/mL	445.297 CFU/mL
<i>Bordetella holmesii</i>	87.674 CFU/mL	-	-	-	-
<i>Bordetella bronchiseptica</i>	-	-	-	-	1101.8 CFU/mL

- Analytical specificity**

The design of oligonucleotides (primers and probes) was validated *in silico* by BLAST alignment. The comparison of the sequences obtained shows a specific detection of the BIOSYNEX AMPLIQUICK *Bordetella* targets.

No primers or probes detect bacterial DNA other than that of the 4 species of *Bordetella* targeted.

18 DNA extracts from qualified and referenced strains were tested.

11 DNA extracts from the French National Reference Center positive for whooping cough and other *Bordetellosis*, 3 control DNA sequences from a biobank (Amplirun®, Vircell) and 4 cultures of standard bacterial strain (Zeptomatrix®).

		Control DNA	IS481	ptxA-Pr	IS1001	FLA1	FLA2	
Biobank		<i>Bordetella pertussis</i> DNA	+	+	-	-	-	
		<i>Bordetella parapertussis</i> DNA	-	-	+	+	+	
		<i>Bordetella holmesii</i> DNA	+	-	-	-	-	
		Species	Strain	IS481	ptxA-Pr	IS1001	FLA1	FLA2
Standard culture		<i>Bordetella pertussis</i>	A639	+	+	-	-	-
		<i>Bordetella parapertussis</i>	A747	-	-	+	+	+
		<i>Bordetella holmesii</i>	F061	+	-	-	-	-
		<i>Bordetella bronchiseptica</i>	Z341	-	-	-	+	-
NRC	<i>Bordetella bronchiseptica</i>	RB50	+	-	-	-	+	
		LAR	+	-	-	+	+	
		FR3523	-	-	+	-	+	
	<i>Bordetella Petrii</i>	KMBW	-	-	-	-	-	
		FR5141	-	-	-	-	-	

<i>Bordetella avium</i>	CIP103348T	-	-	-	-	-
	FR6062	-	-	-	-	-
<i>Bordetella hinzii</i>	CIP104527T	-	-	-	-	-
	FR5948	-	-	-	-	-
<i>Bordetella trematum</i>	CIP105351T	-	-	-	-	-
	FR5853	-	-	-	-	-

Cross-reactivity

A panel of 76 DNA samples and 38 RNA samples from a biobank listed in the following tables were tested with the BIOSYNEX AMPLIQUICK Bordetella kit. For all these samples, no amplification of the targets of interest was observed.

DNA		
<i>Acanthamoeba castellanii</i>	<i>Escherichia coli</i> (ETEC)	<i>Neisseria meningitidis</i> Sg A
<i>Adenovirus</i>	<i>Escherichia coli</i> (VTEC)	<i>Neisseria meningitidis</i> Sg B
<i>Adenovirus 41</i>	<i>Francisella tularensis</i>	<i>Neisseria meningitidis</i> Sg C
<i>Aspergillus fumigatus</i>	<i>Gardnerella vaginalis</i>	<i>Papillomavirus</i> type 16
<i>Bacillus cereus</i>	<i>Giardia intestinalis</i>	<i>Papillomavirus</i> type 18
<i>Bartonella henselae</i>	<i>Haemophilus ducreyi</i>	<i>Parvovirus B19</i> (Plasmid)
<i>Bartonella Quintana</i>	<i>Haemophilus influenzae</i>	<i>Rickettsia conorii</i>
<i>Bk Virus</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>
<i>Borrelia afzelii</i>	<i>Herpes simplex 1</i>	<i>Salmonella typhi</i>
<i>Borrelia burgdorferi</i>	<i>Herpes simplex 2</i>	<i>Staphylococcus aureus</i> (MecA-)
<i>Brucella abortus</i>	Hhv-6	<i>Staphylococcus aureus</i> (MecA+)
<i>Campylobacter jejuni</i>	Hhv-8	<i>Streptococcus agalactiae</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae</i> (NDM-1)	<i>Toxoplasma gondii</i>
<i>Candida auris</i>	<i>Legionella pneumophila</i>	<i>Treponema pallidum</i>
<i>Chlamydia trachomatis</i>	<i>Leishmania chagasi</i>	<i>Trichomonas vaginalis</i>
<i>Chlamydophila pneumoniae</i>	<i>Leishmania infantum</i>	<i>Trypanosoma cruzi</i>
<i>Chlamydophila psittaci</i>	<i>Listeria monocytogenes</i>	<i>Ureaplasma parvum</i>
<i>Clostridium difficile</i>	<i>Moraxella catarrhalis</i>	<i>Ureaplasma urealyticum</i>
<i>Coccidioides immitis</i>	<i>Mycobacterium avium</i>	<i>Varicella-Zoster Virus</i>
<i>Coxiella burnetii</i>	<i>Mycobacterium intracellulare</i>	<i>Vibrio cholerae</i>
<i>Cryptosporidium parvum</i>	<i>Mycobacterium kansasii</i>	<i>Yersinia enterocolitica</i>
<i>Enterococcus faecalis</i> (VanB)	<i>Mycobacterium tuberculosis</i>	<i>Borrelia garinii</i>
<i>Enterococcus faecium</i> (VanA)	<i>Mycobacterium ulcerans</i>	<i>Cytomegalovirus</i>
<i>Epstein-Barr Virus</i>	<i>Mycoplasma genitalium</i>	<i>Streptococcus pneumoniae</i>
<i>Escherichia coli</i> (EAEC)	<i>Mycoplasma hominis</i>	
<i>Escherichia coli</i> (EIEC)	<i>Neisseria gonorrhoeae</i>	

RNA		
<i>Coronavirus</i> Oc43	<i>Human Parainfluenza 1</i>	<i>Parainfluenza 4 A</i>
<i>Chikungunya Virus</i>	<i>Influenza A H1</i>	<i>Respiratory Syncytial Virus</i> (Subtype A)
<i>Coronavirus</i>	<i>Influenza A H3</i>	<i>Respiratory Syncytial Virus</i> (Subtype B)
<i>Coronavirus SARS (2003)</i>	<i>Influenza A H5</i>	<i>Rhinovirus</i>
<i>Coxsackie A6</i>	<i>Influenza B</i>	<i>Rotavirus</i>
<i>Coxsackie B1</i>	<i>Measles</i>	<i>Rubella</i>
<i>Coxsackie B5</i>	<i>MERS Coronavirus</i>	<i>SARS-CoV-2</i>
<i>Dengue 1 Virus</i>	<i>Mumps</i>	<i>Tick-Borne Encephalitis Virus</i>
<i>Dengue 2 Virus</i>	<i>Norovirus</i>	<i>West Nile Virus</i>
<i>Dengue 3 Virus</i>	<i>Novel Influenza A H1n1</i>	<i>Yellow Fever Virus</i>
<i>Dengue 4 Virus</i>	<i>Parainfluenza 1</i>	<i>Zika Virus</i>
<i>Echovirus 5</i>	<i>Parainfluenza 2</i>	<i>Zika Virus</i> (Asian Lineage)
<i>Enterovirus 68</i>	<i>Parainfluenza 3</i>	

- Interference study**

The presence of PCR inhibitors in the sample may induce a positive or negative interference on test results. The presence of various inhibitors in samples was tested.

Nasal sprays

No positive or negative interference was found by loading up to 10% in the BIOSYNEX AMPLIQUICK Sample Collection transport medium from a matrix of nasopharyngeal secretions treated with the BIOSYNEX AMPLIQUICK Lysis kit, with the following substances:

Product	Active Substance
Nasacort	Triamcinolone acetonide
Aturgyl	Oxymetazoline
Béconase	Beclometasone dipropionate
Budésonide	Budesonide
Mométasone	Mometasone
Avamys	Fluticasone furoate
Fixorinox	Fluticasone propionate

Blood

The presence of blood in nasopharyngeal swabs can induce a positive or negative interference on the PCR test result. Our study shows that blood concentrations up to 1% do not interfere with the PCR reaction when the BIOSYNEX AMPLIQUICK Bordetella kit is used with the BIOSYNEX AMPLIQUICK Sample Collection and BIOSYNEX AMPLIQUICK Lysis kits. For blood concentrations from 2% and above, there is a decrease in sensitivity. Therefore, blood concentration in the sample should not exceed 1%.

• Clinical performances

The clinical performances were determined on 146 samples from 4 different laboratories that were qualified positive or negative for *B. pertussis* or *B. parapertussis* using CE marked PCR tests referenced by the French Health Ministry.

These tests allow the identification of these two strains through the amplification of the IS481 and IS1001 genes, which are also present in the *B. holmesii* and *B. bronchiseptica* strains.

The samples were qualified with the BIOSYNEX AMPLIQUICK Bordetella kit on purified DNA extracts or after treatment with the BIOSYNEX AMPLIQUICK Lysis kit.

	Number of samples	<i>B. pertussis</i>		<i>B. parapertussis</i>		Co-infections	Negative for Bordetella
		+	-	+	-		
Laboratory 1	63	58	1	8	51	8*	5
Laboratory 2	12	12	0	2	10	2	0
Laboratory 3	31	0	0	0	0	0	31
Laboratory 4	40	0	0	0	0	0	40

*Among the 8 samples qualified as positive for *B. pertussis* and *B. parapertussis* from laboratory 1, no co-infection was confirmed using the BIOSYNEX AMPLIQUICK Bordetella kit as well as another CE marked kit on the market.

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Bordetella kit on samples after DNA extraction/purification with a commercial kit:

BIOSYNEX AMPLIQUICK Bordetella		Reference PCR	
		Positive	Negative
	Positive	69	0
	Negative	1	76

Sensitivity: 98.57% **Specificity: 100%**
(95%CI: 92.30% to 99.96%) (95%CI: 95.26% to 100%)

The contingency table below shows the performance of the BIOSYNEX AMPLIQUICK Bordetella on samples treated with BIOSYNEX AMPLIQUICK Lysis kit:

		Reference PCR	
		Positive	Negative
BIOSYNEX AMPLIQUICK Bordetella	Positive	66	0
	Negative	4	76

Sensitivity : 94.29% (95%CI: 86.01% to 98.42%) **Specificity : 100%** (95%CI: 95.26% to 100%)

For this study, we did not have access to clinical samples from patients positive for *B. bronchiseptica* and *B. holmesii* (rare samples). The performances of the kit for these two species were validated using strains provided by the French National Reference Center for whooping cough and other Bordetellosis.

• **Precision**

The precision data in the context of the BIOSYNEX AMPLIQUICK Bordetella kit was determined based on 4 conditions:

- Intra-assay variability (intra-day, within a same experiment)
- Inter-days variability
- Inter-operators variability
- Inter-lots variability

The variability data are expressed in terms of mean value, standard deviation and coefficient of variation, based on the threshold cycle of quantification (Cq) values of the Bordetella DNA.

Intra-assay variability data:

Two sample dilutions were tested for each Bordetella specie detected by the kit: 1 high (+++) corresponding to 1000 copies/μL and 1 low (+) corresponding to 10 copies/μL. Two different low dilutions were used for *B. pertussis* and *B. parapertussis*, as their two positive targets have different sensitivities and the IS481 and IS1001 genes are multicopy and the ptxA-Pr and FLA genes are monocopy.

Each sample was tested 30 times.

Extracted DNA

	Sample dilution	Target	Mix 1			Sample dilution	Target	Mix 2					
			Mean Cq	SD	CV %			Mean Cq	SD	CV %			
<i>B pertussis</i>	+++	IS481	13.3	0.1	0.9	<i>B parapertussis</i>	+++	IS1001	15.8	0.3	1.7		
	+		(ptxA-Pr)	27.4	1.1		4.0		+	(FLA1&2)	29.1	0.1	0.5
	-		(IS481)	34.7	0.3		1.0		+	(IS1001)	35.1	0.7	1.9
	+++	ptxA-Pr	19.6	0.2	0.8	+++	FLA1	21.5	0.2	0.8			
	+		(ptxA-Pr)	34.3	0.7	2.2		+	(FLA1&2)	35.2	0.4	1.2	
<i>B holmesii</i>	+++	IS481	21.2	0.1	0.7	<i>B bronchiseptica</i>	+++	FLA 2	21.8	0.1	0.6		
	+		36.2	0.5	1.3		+		36.6	0.6	1.7		

Pre-treatment with BIOSYNEX AMPLIQUICK Lysis kit

	Sample dilution	Target	Mix 1			Sample dilution	Target	Mix 2					
			Mean Cq	SD	CV %			Mean Cq	SD	CV %			
<i>B pertussis</i>	+++	IS481	14.7	0.2	1.0	<i>B parapertussis</i>	+++	IS1001	13.7	0.2	1.2		
	+		(ptxA-Pr)	28.8	1.0		3.4		+	(FLA1&2)	31.3	0.2	0.8
	+		(IS481)	34.4	1.2		3.4		+	(IS1001)	34.5	0.8	2.4
	+++	ptxA-Pr	22.4	0.1	0.6	+++	FLA1	19.2	0.2	1.1			

	+		35.5	1.1	3.0		+		37.1	0.6	1.7
	(ptxA-Pr)						(FLA1&2)				
							+++	FLA2	18.9	0.2	1.1
							(FLA1&2)		34.9	0.5	1.4
<i>B. holmesii</i>	+++	IS481	21.5	0.0	0.2	<i>B. bronchiseptica</i>	+++	FLA 2	18.6	0.2	0.8
	+		35.4	1.4	4.0		+		34.8	0.3	0.8

Inter-assays variability data:

Two sample dilutions were tested for each Bordetella specie detected by the kit: 1 high (+++) corresponding to 1000 copies/μL and 1 low (+) corresponding to 10 copies/μL.

Each sample is tested in duplicate twice each day, during 5 days.

Negative sample give negative results.

Gene	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
IS481	27.7	0.1	0.4	33.9	0.2	0.6
ptxA-Pr	27.3	0.1	0.3	33.4	0.2	0.5
IS1001	27.6	0.1	0.4	33.6	0.3	0.8
FLA1	29.5	0.2	0.6	35.7	0.1	0.4
FLA2	28.6	0.2	0.8	34.7	0.3	0.7

Inter-operators variability data:

Two sample dilutions were tested for each Bordetella specie detected by the kit: 1 high (+++) corresponding to 1000 copies/μL and 1 low (+) corresponding to 10 copies/μL.

twice each day by two different operators, during 5 days.

Negative sample give negative results.

Gene	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
IS481	27.6	0.1	0.4	33.9	0.0	0.1
ptxA-Pr	27.2	0.1	0.4	33.4	0.0	0.1
IS1001	27.5	0.1	0.2	33.6	0.0	0.1
FLA1	29.4	0.1	0.4	35.6	0.0	0.1
FLA2	28.6	0.0	0.0	34.8	0.0	0.1

Inter-lots variability data:

Two sample dilutions were tested for each Bordetella specie detected by the kit: 1 high (+++) corresponding to 1000 copies/μL and 1 low (+) corresponding to 10 copies/μL.

Each sample was tested in triplicate or duplicate on two different lots.

Negative sample give negative results.

Gene	+++			+		
	Average Cq value sample +++	Standard deviation	Coefficient of variation %	Average Cq value sample +	Standard deviation	Coefficient of variation %
IS481	26.9	0.2	0.6	33.4	0.4	1.1
ptxA-Pr	27.4	0.1	0.5	33.7	0.7	2.0
IS1001	26.9	0.0	0.0	33.0	0.3	0.9
FLA1	29.1	0.1	0.5	35.1	0.0	0.0
FLA2	28.5	0.0	0.1	34.6	0.1	0.3




12 I LIMITS

1. Failure to follow any instruction of the IFU may adversely affect the test performance and/or invalidate the test result.
2. As for any diagnostic test, results must be considered with other clinical information available to the physician. A negative result can never rule out the presence of Bordetella in the sample (e.g. genes may be present at a concentration below the minimum detection limit of the test, interference) and a positive result can never ensure the presence of Bordetella in the sample (e.g. contamination). A definitive diagnosis can only be made by the physician after evaluation of all clinical and laboratory data.
3. Samples with more than 1% blood may interfere with the test result. In case of blood in the sample, consider the result with caution.

13 I BIBLIOGRAPHY

- Pittet LF, Emonet S, Schrenzel J, Siegrist CA, Posfay-Barbe KM. Bordetella holmesii : an under-recognized Bordetella species. The Lancet Infectious Diseases. 2014; 14(6):510-9.
- Pittet LF, Emonet S, François P, Bonetti EJ, Schrenzel J, Hug M, Altwegg M, Siegrist CA, Posfay-Barbe KM. Diagnosis of whooping cough in Switzerland: Differentiating Bordetella pertussis from Bordetella holmesii by polymerase chain reaction. PLoS One. 2014; 9(2): e88936
- Sanz JC, Abad R, Sanz C, Miguel A. Diagnóstico diferencial de Bordetella bronchiseptica por RT-PCR en un niño con tos paroxística sin antecedentes patológicos previos. Enferm Infecc Microbiol Clin. 2019. <https://doi.org/10.1016/j.eimc.2018.09.016>

14 I INTERNATIONAL SYMBOLS

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Catalog number
	<i>In vitro</i> diagnostic medical device		Temperature limit		Do not reuse
	Manufacturer		Batch code		Use-by date
	Keep away from sunlight		Master Mix		Microplate
	Negative control		Positive control		Unique Device Identifier
	Bag of cap strips		Do not use if package is damaged and consult instructions for use		Authorized Representative in Switzerland
	Importer		Qualitative detection by PCR		

15 I MANUFACTURER INFORMATION

BIO SYNEX S.A.
 22 boulevard Sébastien Brant
 67400 ILLKIRCH-
 GRAFFENSTADEN – France
 Standard:
 Phone: +33 3 88 78 78 87
www.biosynex.com
 Contacts France:
 Phone: +33 3 88 77 57 00
service.clients@biosynex.com
 Contacts other countries:
 Phone: +33 3 88 77 57 52
sales@biosynex.com
 Customer service:
 Phone : +33 3 88 77 57 25
tech.support@biosynex.com

BIO SYNEX S.A.
 Route de Rossemaison 100
 2800 DELEMONT - Switzerland

Latest changes :
 Update of the table of symbols and
 manufacturer's contact details.
 §7: correction on suction sampling

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Et panel med 76 DNA-prøver og 38 RNA-prøver fra en biobank som er oppført i følgende tabeller, ble testet med BIOSYNEXAMPLIQUICK © Bordetella-settet. For alle disse prøvene ble det ikke observert noen forsterkning av målene av interesse.
BIO SYNEX SWISS S.A.