BIOSYNEX AMPLIQUICK Fecal Bacteriology

REF 3150046_SEC01 / 3150046_SEC02 / 3150046_BLK96

FR (FRANÇAIS / FRENCH)	1
EN (ENGLISH)	26

CONTENTS

BIOSYNEX AMPLIQUICK Fecal Bacteriology

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12. BIBLIOGRAPHIE

1. Graves, N.S. (2013). Acute Gastroenteritis. Primary Care: Clinics in Office Practice 40, 727–741. 10.1016/j.pop.2013.05.006.
2. Bányai, K., Estes, M.K., Martella, V., and Parashar, U.D. (2018). Viral gastroenteritis. The Lancet 392, 175–186. 10.1016/S0140-6736(18)31128-0.
3. Centers for Disease Control and Prevention (CDC). "Outbreak of Salmonella Infections Linked to Backyard Poultry." 2020.
4. European Food Safety Authority (EFSA). "The European Union summary report on trends and sources of zoonoses, zoonotic agents and food-borne outbreaks in 2017." EFSA Journal, 2018; 16(12):5500.
5. Brenner FW, Villar RG, Angulo FJ, Tauxe R, Swaminathan B. 2000. Salmonella Nomenclature. J Clin Microbiol 38: https://doi.org/10.1128/jcm.38.7.2465-2467.2000.

6. European Centre for Disease Prevention and Control. Healthcare-associated infections acquired in intensive care units. In: ECDC. Annual Epidemiological Report for 2020. Stockholm: ECDC; 2024.
7. Fukushima H, Shimizu S, Inatsu Y. Yersinia enterocolitica and Yersinia pseudotuberculosis Detection in Foods. J Pathog. 2011;2011:735308. doi: 10.4061/2011/735308. Epub 2011 Oct 5. PMID: 22567341; PMCID: PMC3335482.
8. Kaakoush NO, Castaño-Rodríguez N, Mitchell HM, Man SM. Global Epidemiology of Campylobacter Infection. Clin Microbiol Rev. 2015 Jul;28(3):687-720. doi: 10.1128/CMR.00006-15. PMID: 26062576; PMCID: PMC4462680.
9. World Health Organization (WHO). "Campylobacter." Fact sheet, 2020.
10. Janda JM, Abbott SL. The genus Aeromonas: taxonomy, pathogenicity, and infection. Clin Microbiol Rev. 2010 Jan;23(1):35-73. doi: 10.1128/CMR.00039-09. PMID: 20065325; PMCID: PMC2806660.
11. Yuwono C, Wehrhahn MC, Liu F, Zhang L. 2023. Enteric Aeromonas Infection: a Common Enteric Bacterial Infection with a Novel Infection Pattern Detected in an Australian Population with Gastroenteritis. Microbiol Spectr 11:e00286-23. <https://doi.org/10.1128/spectrum.00286-23>.
12. Guh AY, Kuty PK. Ann Intern Med. 2018 October 02; 169(7)
13. Spigaglia P *et al*. Clostridium difficile causing pediatric infections: New findings from a hospital-based study in Italy. Anaerobe. 2017 Dec;48:262-268. doi: 10.1016/j.anaerobe.2017.10.008.
14. European Centre for Disease Prevention and Control. Antimicrobial resistance surveillance in Europe 2011. Annual Report of the European Antimicrobial Resistance Surveillance Network (EARS-Net). Stockholm: ECDC; 2012.
15. Janda JM, Abbott SL, McIver CJ. Plesiomonas shigelloides Revisited. Clin Microbiol Rev. 2016 Apr;29(2):349-74. doi: 10.1128/CMR.00103-15. PMID: 26960939; PMCID: PMC4786884.
16. World Health Organization (WHO). "Cholera." Fact sheet, 2020.
17. Baker-Austin, C., Stockley, L., Rangdale, R. and Martinez-Urtaza, J. (2010), Environmental occurrence and clinical impact of Vibrio vulnificus and Vibrio parahaemolyticus: a European perspective. Environmental Microbiology Reports, 2: 7-18. <https://doi.org/10.1111/j.1758-2229.2009.00096.x>
18. Pakbin B, Brück WM, Rossen JWA. Virulence Factors of Enteric Pathogenic Escherichia coli: A Review. Int J Mol Sci. 2021 Sep 14;22(18):9922. doi: 10.3390/ijms22189922. PMID: 34576083; PMCID: PMC8468683.
19. Shigella species and characteristics: "Shigellose : symptômes, traitement, prévention." Institut Pasteur.
20. Schmidt, H., Scheef, J., Huppertz, H.I., Frosch, M., and Karch, H. (1999). Escherichia coli O157:H7 and O157:H Σ Strains That Do Not Produce Shiga Toxin: Phenotypic and Genetic Characterization of Isolates Associated with Diarrhea and Hemolytic-Uremic Syndrome. J. CLIN. MICROBIOL. 37.
21. Stephan, R., Zhang, W., Bielaszewska, M., Mellmann, A., and Karch, H. (2009). Phenotypic and Genotypic Traits of Shiga Toxin–Negative *E. coli* O157:H7/H Σ - Bovine and Porcine Strains. Foodborne Pathogens and Disease 6, 235–243. <https://doi.org/10.1089/fpd.2008.0205>.

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BIOSYNEX AMPLIQUICK Fecal Bacteriology

PCR FOR THE QUALITATIVE DETECTION OF GASTROINTESTINAL BACTERIA IN STOOL SAMPLES.

In vitro diagnostic medical device for professional use.

REF 3150046_SEC01 / 3150046_SEC02 / 3150046_BLK96

1. INTENDED PURPOSE

BIOSYNEX AMPLIQUICK Fecal Bacteriology is a molecular *in vitro* diagnostic test for the qualitative detection of 15 pathogenous agents commonly associated with gastroenteritis : *Campylobacter spp.*, *Yersinia enterocolitica* and *Yersinia pseudotuberculosis*, *Salmonella spp.*, *Clostridioides difficile*, *Clostridioides difficile* hypervirulent strains, *Aeromonas spp.*, Cholera toxin, *Vibrio spp.*, *Escherichia coli* O157 strain, Enteroinvasive *Escherichia coli* (EIEC)/*Shigella spp.*, Enteropathogenic *Escherichia coli* (EPEC), Enterotoxigenic *Escherichia coli* (ETEC), Enteroaggregative *Escherichia coli* (EAEC), Enterohemorrhagic *Escherichia coli* (EHEC) and *Plesiomonas shigelloides*.

The test is an aid in the diagnosis of bacterial gastroenteritis performed using a purified DNA extract obtained from stool samples collected with BIOSYNEX AMPLIQUICK Fecal Pretreatment or with FecalSwab™ (Copan). For a syndromic diagnosis, the PCR can be performed simultaneously with the BIOSYNEX AMPLIQUICK Fecal Virology test (ref. 3150047 and related references) using the same samples and the same amplification program. This non-automated test is intended for *in vitro* diagnostic use in laboratory by professionals only.

2. CLINICAL SUMMARY

Gastroenteritis is the inflammation of the gastrointestinal tract that causes diarrhea, abdominal pain, nausea, fever and sometimes vomiting. This inflammation affects nearly 1.7 billion children every year and is a leading cause of death in children under five years old mainly due to severe dehydration (World Health Organization). It can be provoked by many different infectious agents such as parasites (10-15%), bacteria (15-20%) and viruses (50-70%) that are spread by the fecal-oral route through the contamination of food and water.¹

The disease symptoms usually occur 6 hours to 5 days after the infection and can last from 2-7 days. Diarrhea is the most common clinical symptom of gastrointestinal infections. Other symptoms may appear along with diarrhea such as abdominal pain, nausea, vomiting, hemorrhagic diarrhea (in the case of EHEC infection for instance) and fever. These symptoms can be severe, and in the lack of accurate and prompt diagnosis, it can lead to serious consequences, including death. The diagnosis of such disease should be carried out on patients suffering from acute diarrhea and/or any of the symptoms listed above, to allow the administration of specific treatment in patients with severe gastroenteritis especially for those with bloody stool, prolonged illness or with chronic conditions.²

The most common pathogenic bacteria that trigger gastrointestinal infections include *Salmonella spp.*, *Yersinia spp.*, *Campylobacter spp.*, *Aeromonas spp.*, *Clostridioides difficile*, *Plesiomonas shigelloides*, *Vibrio spp.*, *Escherichia coli*, and *Shigella spp.*

- *Salmonella* remains a major cause of acute gastroenteritis, with notable outbreaks in Europe (2017) and the US (2020).^{3,4} The *Salmonella* genus, part of the Enterobacteriaceae family, comprises two species: *Salmonella bongori* and *Salmonella enterica*, both of which are Gram-negative rods. *Salmonella enterica* is the most common species and is divided into six subspecies, comprising more than 2,500 different serotypes, the best known of which are *Salmonella Typhimurium* and *Salmonella Enteritidis*.⁵ *Salmonella bongori* is a distinct species of *Salmonella*, less common than *Salmonella enterica*. Unlike *Salmonella enterica*, *Salmonella bongori* lacks certain islands of pathogenicity (SPI2), which makes it less common than *Salmonella enterica*. *Salmonella* is mainly transmitted through the consumption of contaminated

food, such as meat, eggs, raw milk and certain dairy products. Direct contact with infected animals, particularly reptiles, can also be a source of contamination.

- The *Yersinia* genus, also within the Enterobacteriaceae family, comprises 11 Gram-negative bacillus species including three human pathogens *Y. enterocolitica*, *Y. pseudotuberculosis*, and *Y. pestis*. *Y. enterocolitica* and *Y. pseudotuberculosis* are associated with gastrointestinal infections, particularly in winter. *Yersinia* is a leading cause of gastroenteritis, with the European Center for Disease Prevention and Control (ECDC) reporting 7,048 cases of yersiniosis in 2019.⁶ *Y. enterocolitica* is more common in cooler climates and can grow at refrigeration temperatures, making it a concern for food safety.⁷ *Y. pseudotuberculosis* causes a disease like tuberculosis in animals, hence its name; *Y. pseudotuberculosis* infections are less common than those caused by *Y. enterocolitica*. Both species share several virulence factors, such as the ability to invade intestinal cells and evade the host immune response. They are both zoonotic pathogens, meaning they can be transmitted from animals to humans.
- *Campylobacter* infections are among the most widespread infectious diseases, increasing in both developed and developing countries over the past decade.⁸ According to the WHO, *Campylobacter* is a global cause of diarrheal diseases.⁹ The genus includes several Gram-negative species that are spiral, curved, or rod-shaped. *C. jejuni* and *C. coli* are major causes of human gastroenteritis, while species like *C. upsaliensis* and *C. fetus* are emerging pathogens. The main reservoirs are animals, particularly birds, pigs, and domestic pets.
- *Aeromonas* species are emerging pathogens causing various human diseases, including septicemia, wound infections, and gastroenteritis.^{10,11} They are commonly found in aquatic environments. The incidence of *Aeromonas* in gastroenteritis ranges from 2% to 10%, depending on the country. The genus includes 14 Gram-negative rod-shaped species, with *A. caviae*, *A. veronii*, *A. dhakensis*, and *A. hydrophila* being most associated with human infections. *A. caviae* is known for its ability to form biofilms, which can enhance its survival in various environments and contribute to its pathogenicity. *A. veronii* produces various virulence factors and it is a significant cause of gastroenteritis, especially in children. *A. dhakensis* can cause severe infections in humans and it is considered more virulent than *A. hydrophila*, with which it is often confused. Molecular tools are essential for accurate detection.¹¹ *A. hydrophila* is widely distributed in freshwater environments and is known for its ability to survive in various conditions, including low temperatures; moreover, it resists to chlorination and can persist in contaminated water.
- *Clostridioides difficile* is a Gram-positive, spore-forming bacterium that can cause severe gastrointestinal infections, particularly after antibiotic use.¹² Its virulence is primarily mediated by two enterotoxins, TcdA and TcdB, which destroy the actin cytoskeleton of intestinal cells. TcdB is considered more potent and widespread than TcdA, and it is primarily responsible for the cytotoxic effect seen in *C. difficile* infections. Moreover, it is possible to distinguish a hypervirulent strain which secretes more TcdA and TcdB thanks to a deletion in the regulatory sequence of the operon. In addition of those two toxins, some strains produce a third toxin, the binary toxin *C. difficile* transferase (CDT), which can also contribute to the pathogen virulence and disease. Interestingly, some *C. difficile* strains, such as RT033 and RT442, do not possess the classical TcdA/TcdB toxins, produce the binary toxin, provoke gastroenteritis and are mostly missed during diagnosis.¹³ *C. difficile* caused over 453,000 gastroenteritis cases in the US in 2011 and is listed in the highest priority category of antimicrobial resistance threats by the ECDC.¹⁴ The infections are challenging to treat for several reasons: highly resistant and long-lived spores; resistance to antibiotics; proliferation after disruption of the normal gut microbiota; biofilm formation and recurrent infections. When antibiotic treatment fails, Fecal Microbiota Transplantation is an alternative approach to get rid of recurrent infections. Probiotics and immunotherapy are also being explored as adjunct therapies to enhance treatment efficacy and prevent recurrence.
- *Plesiomonas shigelloides* is a Gram-negative, rod-shaped bacterium that belongs to the Enterobacteriaceae family. It is globally distributed but predominantly found in tropical areas, where it primarily inhabits freshwater, especially during warmer months.¹⁵ It is often associated

with aquatic animals like fish, shellfish and amphibians but also some terrestrial animals. *P. shigelloides* can cause gastroenteritis, often leading to diarrhea, especially after consuming contaminated water or seafood.¹⁵

- The *Vibrio* genus includes Gram-negative, rod-shaped bacteria that naturally inhabit freshwater and marine environments. The most common pathogenic species are *V. cholerae*, *V. vulnificus*, and *V. parahaemolyticus*. *V. cholerae* causes severe diarrheal illness, with an estimated 3-5 million infections worldwide annually. The species has over 200 serogroups, with O1 being the most predominant.¹⁶ *V. parahaemolyticus* is ubiquitous in temperate and tropical coastal areas, causing around 30,000 infections annually. *V. vulnificus* is an opportunistic pathogen affecting patients with chronic diseases.¹⁷
- *Escherichia coli* is a Gram-negative, rod-shaped bacterium belonging to the Enterobacteriaceae family. From the first hours after birth, it colonizes the gastrointestinal tract of infants, becoming part of a normal and healthy microbiota through a commensal relationship. However, *E. coli* can cause disease under certain conditions, such as in immunocompromised hosts or when the gastrointestinal barrier is breached. Pathogenic *E. coli* is responsible for approximately 2 million deaths annually, particularly in children under five years old. Each *E. coli* pathotype has its characteristic pathogenic mechanism and a specific profile of virulence factors.¹⁸ We can distinguish four main virulence classes: colonization, fitness, toxins and effectors, each of which consists of specific virulence factor with a defined function and activity. Seven pathotypes of *E. coli* have been identified, including the five most associated with gastroenteritis: enteroaggregative (EAEC), enteroinvasive (EIEC), shiga-toxin secreting *E. coli* (STEC, including enterohemorrhagic *E. coli* (EHEC), with the O157 serotype being the most common in recent European outbreaks, enteropathogenic (EPEC) and enterotoxigenic (ETEC). These pathotypes can differ in colonization and/or cellular invasion strategies; moreover, they are distinguished by the presence or absence of various virulence genes, such as the stx1/stx2 toxins in EHEC.¹⁸ These virulence genes can be encoded by plasmids or in chromosomal regions. It is crucial to distinguish the pathotype of *E. coli* in the context of gastroenteritis, as different pathotypes can cause a variety of symptoms and require specific therapeutic approaches. Some pathotypes of *E. coli*, such as enterohaemorrhagic *E. coli* (EHEC) or enterotoxigenic *E. coli* (ETEC), can cause symptoms ranging from watery diarrhea to bloody diarrhea and severe abdominal cramps. For example, EHEC can lead to hemorrhagic colitis and serious extraintestinal complications such as hemolytic uremic syndrome. In addition, knowing the pathotype helps to identify potential sources of contamination and to implement appropriate preventive measures. For example, EHEC is often associated with the consumption of undercooked meat and unpasteurized dairy products. Treatments may vary as well depending on the pathotype: for example, the use of antibiotics may be contraindicated for certain pathotypes such as EHEC, as this can worsen symptoms. Finally, it is also of epidemiological interest, as the surveillance and the control of epidemics require precise identification of the pathotype to understand modes of transmission and implement effective public health interventions.
- *Shigella* is a leading cause of mortality due to diarrheal diseases, accounting for 212,000 deaths in 2016. The *Shigella* genus includes Gram-negative, rod-shaped bacteria that are genetically closely related to *E. coli*. There are four *Shigella* species: *S. dysenteriae*, *S. flexneri*, *S. boydii*, and *S. sonnei*.¹⁹ *Shigella* invades the epithelial cells of the intestinal mucosa, leading to severe inflammation and tissue destruction. It also produces several toxins contributing to its pathogenicity. *Shigella* and EIEC are so closely related genetically that they are often considered to be part of the same pathovar within the genus *Escherichia*; they also share the same virulence plasmid. *Shigella*, particularly *S. flexneri*, tends to cause more severe disease with a pronounced inflammatory response compared to EIEC. In France and industrialized countries, the most common serogroup is *S. sonnei*, but it is also emerging worldwide, including in developing countries. In the past, *Shigella dysenteriae* serotype 1 (or *Shiga bacillus*) caused deadly epidemics. Since 2011, no more strains of this serotype have been reported.

3. TEST PRINCIPLE

The assay consists of a multiplex PCR allowing the specific amplification and the simultaneous detection of bacterial genes and/or toxins specific to *Campylobacter*, *Salmonella*, *Vibrio* and *Aeromonas* species, *Yersinia enterocolitica*, *Yersinia pseudotuberculosis*, cholera toxin genes, *Plesiomonas shigelloides*, 5 different pathotypes of *E. coli* EPEC, EAEC, ETEC, EHEC, EIEC and *Shigella spp.*, as well as an indication of the presence of the *E. coli* O157 strain.

The assay is composed of 5 different master mixes. The targets and the detection channels for each master mix are detailed in the table below.

Master Mixes	Channel	Bacterial pathogens	Target genes/sequences
A	FAM	<i>Campylobacter spp.</i>	16S (<i>C. jejuni/C. coli</i>) lpxA (<i>C. upsaliensis</i>) gyrB (<i>C. fetus</i>)
	HEX	<i>Yersinia enterocolitica</i> <i>Yersinia pseudotuberculosis</i>	YST ail
	TxRed	<i>Salmonella spp.</i>	invA
	Cy5	Internal amplification control	CIEZ
B	FAM	<i>Clostridioides difficile</i> Hypervirulent	tcdC ^{A117}
	HEX	<i>Clostridioides difficile</i> ToxB/Binary toxin	tcdB, CDT
	TxRed	<i>Aeromonas spp.</i> (<i>A. hydrophila</i> , <i>A. veronii</i> , <i>A. dhakensis</i> and <i>A. caviae</i>)	gyrB
	Cy5	Internal amplification control	CIEZ sequence
C	FAM	Cholera toxins	ctxAB, tcpA
	HEX	<i>Vibrio spp.</i> (<i>V. cholerae</i> , <i>V. vulnificus</i> , <i>V. Parahaemolyticus</i> , <i>V. alginolyticus</i> and <i>V. fluvialis</i>)	ATPsynthase
	TxRed	<i>E. coli</i> O157 strain	z3276
	Cy5	Internal amplification control	CIEZ sequence
D	FAM	<i>Shigella spp./E. coli</i> EIEC	ipaH
	HEX	<i>E. coli</i> EPEC	eaeA, bfp
	TxRed	<i>E. coli</i> ETEC	lt1/sta
	Cy5	Internal amplification control	CIEZ sequence
E	FAM	<i>Plesiomonas shigelloides</i>	HugA
	HEX	<i>E. coli</i> EAEC	aggR, aap
	TxRed	<i>E. coli</i> EHEC	stx1, stx2
	Cy5	Control in process	Rp49

An exogenous internal control CIEZ, as well as an exogenous control in process Rp49 are included: The Internal Control CIEZ is detected in the Cy5 channel of Mix A, B, C and D while the exogenous Control in Process Rp49 is detected in the Cy5 channel of Mix E.


The test is based on real-time Polymerase Chain Reaction (PCR) technology by fluorescent probe hydrolysis. This method uses probes labeled with FAM, HEX, Texas Red and Cy5 for the amplification of bacterial targets and Cy5-labeled probes for the detection of the sequences used as controls the CIEZ and the Rp49. It should be noted that the detection of the CIEZ internal control validates the PCR amplification step and rules out a possible false negative result due to PCR inhibition. Rp49 is a known DNA sequence from *Drosophila* that is absent in the human and bacterial genomes and is therefore added to stool samples as an extraction control, also known as a control in process. The

supplementation of the Rp49 exogenous sequence in the fecal samples and its detection allow to control the quality of the biological sample and to validate the genomic material purification step.

The increase in the fluorescence signal is only detected if the target sequence is complementary to the amplified probe present in the sample. The fluorescent signal is therefore directly proportional to the amplification of the target during the amplification phase. The Cq value (quantification cycle) corresponds to the number of cycles during which the fluorescence increases exponentially above the background noise.

BIO SYNEX AMPLIQUICK Fecal Bacteriology is available in two formats :

- In bulk tubes (ref. 3150046_BLK96): each master mix is provided in a microtube for the user to dispense in PCR plates.
- In prefilled microplates (ref. 3150046_SEC01 and 3150046_SEC02): each master mix is distributed in one well of each strip according to order shown below.

Regardless the packaging, the Master mixes are ready-to-use and contains dNTPs, MgCl₂, primers and fluorescent probes, Taq polymerase enzyme, and reaction buffer.

BIO SYNEX AMPLIQUICK Fecal Bacteriology must be used with stool sample collected and treated with BIO SYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or collected with FecalSwab™ (Copan). Regardless the collection system used, the nucleic acids must be extracted and purified before performing the amplification step.

4. REQUIRED MATERIAL

3150046_SEC01 or 3150046_SEC02	3150046_BLK96
Kit composition	
<ul style="list-style-type: none"> - 5 dividable 96-wells microplates prefilled with the five master mixes (Mmix) - 1 tube of positive control (CONTROL +, red cap) - 1 tube of negative control (CONTROL -, green cap) - 1 tube of control in process (CONTROL IP, yellow cap) - 1 tube of master mix control (CONTROL Mmix, blue cap) - 5 bags of optical cap strips 	<ul style="list-style-type: none"> - 1 tube of Master Mix A (Mmix A, white cap marked B1) - 1 tube of Master Mix B (Mmix B, white cap marked B2) - 1 tube of Master Mix C (Mmix C, white cap marked B3) - 1 tube of Master Mix D (Mmix D, white cap marked B4) - 1 tube of Master Mix E (Mmix E, white cap marked B5) - 1 tube of positive control (CONTROL +, red cap) - 1 tube of negative control (CONTROL -, green cap) - 2 tubes of control in process (CONTROL IP, yellow cap)

	- 1 tube of master mix control (CONTROL Mmix, blue cap)
Material required but not provided	
- BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or FecalSwab™ transport medium (Copan, ref. 470CE) - DNA extraction kit - Powder-free disposable gloves - Pipettes and filtered tips - PCR microplate or microtubes centrifuge - qPCR Thermal Cycler	- BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or FecalSwab™ transport medium (Copan, ref. 470CE) - DNA extraction kit - Powder-free disposable gloves - Pipettes and filtered tips - PCR microplate or microtubes with optical caps - PCR microplate or microtubes centrifuge - qPCR Thermal Cycler

The PCR thermal cycler used for the test must have the following main characteristics:

- Open system
- Real-time quantitative PCR assays
- Programmable Thermal Cycler block

Reference	Thermal cycling block
3150046_SEC01	0.1mL low profile
3150046_SEC02	0.2mL high profile

- Excitation source: Leds, lamp or laser
- Filter sets (excitation/emission wavelengths) suitable for the detection of "reporter" fluorophores of the FAM, HEX, TxRed and Cy5 probes.
- Connection with a computer using specific analysis software that allow the recovery of fluorescence data, absolute quantification assays, and the interpretation of results.

The kit has been developed and validated with the following real-time PCR thermal cycler:

- CFX96 Touch™ Real-Time PCR Detection System (BioRad)
- CFX96 OPUS™ Real-Time PCR Detection System (BioRad)
- QuantGene 9600 (Bioer)
- QuantStudio™ 5 System (Applied Biosystems™)

If another thermal cycler is used, please perform your own validation of the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit before using the test. We recommend using qualified samples and the negative and positive controls provided in the kit.

5. PRECAUTIONS

- Carefully follow these instructions for use. Failure to follow any instruction of this IFU may adversely affect the test performance and have harmful consequences.
- Follow Good Laboratory Practice and wear powder-free laboratory disposable gloves throughout the test procedure.
- The daily management of many samples and the high sensitivity of the PCR technique can, in the absence of precaution, generate false positive results by contamination. Pre-handling PCR, post-PCR and DNA extraction should therefore be performed in separate rooms. Wear disposable gloves in each room and change them before moving between different rooms.
- The provided strips of caps are for single use only. Do not reuse them. Do not open the PCR tubes or microplates at the end of the test.
- For the prefilled kit formats (ref. 3150046_SEC01 and 3150046_SEC02): Do not use the test if aluminum foil is opened or damaged. Once the aluminum foil is removed, use the test immediately.
- Do not use the kit if it is received defrosted.

- Do not use the kit in case of breakage or leakage. In the event of damage to the packaging only (no breakage or leakage), the kit remains usable.
- Protect the kit from light.
- Centrifuge each tube before opening and open them one after the other, closing them tightly between each one to avoid contamination.
- The positive control (CONTROL +) contains significant amounts of DNA sequences. It can therefore potentially contaminate the other components of the kit if good molecular biology practices are not followed. To limit this risk of contamination, it is recommended to store the positive control separately from the other components of the kit, preferably outside the kit on the first opening.
- The negative and positive controls included in the kit mimic results obtained with negative or positive samples respectively. They must be used with each run.
- To ensure that there is no contamination during handling, it is left to the user to include as part of their internal quality approach, at the frequency chosen, an additional negative control well in which 5µL of molecular biology grade water are added to the Master mix.
- When using negative and positive controls with a series of patients, it is recommended to first deposit the negative control, then deposit the patient samples, and to finish with the deposition of the positive control.
- Dispose of soiled parts or empty kit components in a trash bin suitable for biological waste.
- The device contains material of bacterial or animal origin and should be handled with caution.
- If, in relation to the use of BIOSYNEX AMPLIQUICK Fecal Bacteriology, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country. If in doubt, report it.
- Safety data sheet available on request. Summary of safety and performance will be available online on Eudamed.

6. KIT STORAGE

The kit is shipped frozen. Kit components must arrive frozen. Store the kit at a temperature below or equal to -20°C. Under these conditions, the reagents are stable until the expiry date indicated on the kit label. Do not use the kit and any of its components after the expiry date. Protect the kit and the reagents from direct light.

The master mixes in bulk format (3150046_BLK96) and the Control Master mix in the kit can undergo 16 freeze/thaw cycles without the shelf-life and qualities of the product being altered. The controls (including the Control In Process) of the kit can undergo 30 freeze/thaw cycles without the shelf-life and qualities of the product being altered.

NB: Only thaw the number of well strips or plates needed for the prefilled kit format. As the plates and strips are ready to use, there is no need for repeated thaw/freeze cycles.

7. SPECIMENS/CONTROLS COLLECTION AND STORAGE

The stool sample should be collected with BIOSYNEX AMPLIQUICK Fecal Pretreatment (ref. 3150065) or with FecalSwab™ (Copan) from feces collected in a preservative-free container.

For sample collection and storage, please refer to the instructions for use of BIOSYNEX AMPLIQUICK Fecal Pretreatment or FecalSwab™ (Copan).

Repeated freeze/thaw cycles of samples should be avoided.

The transport of clinical samples must comply with local regulations for the transport of infectious agents.

8. TEST PROCEDURE

1. Sample preparation:

Add 5µL of the control in process (CONTROL IP, yellow cap) in the sample preparation (i.e. FecalSwab™ or BIOSYNEX AMPLIQUICK Fecal Pretreatment tube with the stool sample).

2. Extraction of nucleic acids:

The nucleic acids extraction must be performed before the amplification protocol using a suitable extraction system for bacterial DNA. Please follow the manufacturer's instructions for use.

The positive and negative controls provided in the kit should not undergo extraction and should be used as provided.

The BIOSYNEX AMPLIQUICK Fecal Bacteriology kit has been validated with the following extraction kits:

- Macherey-Nagel NucleoMag® Dx Pathogen (ref. 744215.4) with the MagnetaPure 32 Automated System (Dutscher)
- QIAGEN QIAampFAST DNA Stool Mini Kit (ref. 51604).


If another DNA extraction method or kit is used, please perform your own validation method of the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit using positive and negative qualified samples supplemented with the BIOSYNEX AMPLIQUICK Fecal Bacteriology control in process.

3. Amplification protocol:

It is recommended that the Master mixes be placed on a refrigerated rack or on ice while samples are being added.

For microplates prefilled with ready-to-use master mixes (REF 3150046_SEC01 & 3150046_SEC02):


1. Take from the freezer the number of strips or plates needed.
2. Centrifuge for 10 seconds to remove any droplets from the walls of the wells or the aluminum sealing foil.
3. Add 15µL of Control master mix to wells G and H of a single strip. If you are using multiple strips at the same time, it is not necessary to deposit control master mix into each strip.
4. Add 5 µL of negative control to well G containing the control mix added in the previous step.
5. Add 5µL of extraction eluate from the same sample into the wells containing the master mixes (wells from line A to line E as shown on the picture below).
6. Finally, add 5µL of positive control in well H containing the Master mix control added in step 4
7. Close the plate or strips using the provided caps and centrifuge once more for 10 seconds to eliminate potential bubbles at the surface.

With tubes of bulk ready-to-use master mixes (REF 3150046_BLK96):

1. Remove the Master mixes tubes from the freezer and allow them to thaw on ice or refrigerated rack.
2. Centrifuge the tubes for 10 seconds to remove any droplets from the sides or caps.
3. Dispense 15 μ L of each Master mix into PCR tubes or plate. Each Master Mix must be deposited in a separate well. So, for each patient, 5 wells enable a diagnosis over the entire panel.
4. Dispense 15 μ L of control master mix in two empty wells for the controls.
5. Add 5 μ L of negative control to one of the tube/well containing the control Master mix added in the previous step.
6. Add 5 μ L of extraction eluate from the same sample to the various dedicated reaction wells (Master mixes A to E).
7. Finally, add 5 μ L of positive control in the last tube/well containing the Master mix control added in step 4.
8. Close the wells with optical caps or film.
9. Centrifuge for 10 seconds.
10. Return the master mixes tubes to the freezer.

N.B. To limit the loss of reagent, it is important to pipette the Master Mix without plunging the tip completely into the volume of Master Mix. Alternatively, you could use the “low-binding” or “low retention” tips to limit reagent loss.

Place the plate, strips or tubes in the thermal cycler and run the following amplification program:

Step	Repetitions	Temperature	Duration	Acquisition
Reverse transcription	1x	45°C	15 min	-
Taq Activation	1x	95°C	2 min	-
Denaturation	50x	95°C	5 sec	-
Hybridation / Elongation		60°C	20 sec	yes

If you use :

- CFX96 Touch™ Real-Time PCR Detection System (BioRad)
- CFX96 OPUS™ Real-Time PCR Detection System (BioRad)

To analyze the amplification data, please select **Apply fluorescent drift correction**.

Enter 20 μ L as the reaction volume into the Thermal Cycler program.

Please refer to the operating instructions of the Thermal Cycler used for programming information.

If a syndromic diagnosis is performed using the BIOSYNEX AMPLIQUICK Fecal Virology kit in simultaneously, for the detailed BIOSYNEX AMPLIQUICK Fecal Virology kit test procedure, please refer to the kit instructions for use.

Detection channels settings :

Master mixes	Channel	Bacterial target/genes
A	FAM	<i>Campylobacter spp.</i>
	HEX	<i>Yersinia enterocolitica/Yersinia pseudotuberculosis</i>
	TxRed	<i>Salmonella spp.</i>
	Cy5	Internal amplification control
B	FAM	<i>Clostridioides difficile</i> Hypervirulent
	HEX	<i>Clostridioides difficile</i> ToxB/Binary toxin
	TxRed	<i>Aeromonas spp.</i> (<i>A. hydrophila</i> , <i>A. veronii</i> , <i>A. dhakensis</i> and <i>A. caviae</i>)
	Cy5	Internal amplification control
C	FAM	Cholera toxins
	HEX	<i>Vibrio spp.</i> (<i>V. cholerae</i> , <i>V. vulnificus</i> , <i>V. parahaemolyticus</i>)
	TxRed	<i>E. coli</i> O157 strain
	Cy5	Internal amplification control
D	FAM	<i>Shigella spp./</i> EIEC
	HEX	EPEC
	TxRed	EPEC
	Cy5	Internal amplification control
E	FAM	<i>Plesiomonas shigelloides</i>
	HEX	EAEC
	TxRed	EHEC
	Cy5	Control in process

9. INTERPRETATION OF RESULTS

The threshold of the real-time PCR reaction corresponds to the level of signal that reflects a statistically significant increase in fluorescence emission compared to the baseline signal calculated; the value is set to distinguish a significant amplification signal from background noise. We recommend automatic setting of the threshold by the real-time PCR instrument software, rather than manual setting. If further settings should be done for the test, please contact the thermal cycler manufacturer's technical support team.

Samples:

Signals above the threshold, **and visually consistent with a classical PCR amplification curve**, are considered positive results. Please take into consideration the amplification curves of the positive control from the kit as reference.

Some samples may show atypical curves that are not characteristic of amplification curves. In this case, the result should not be considered as interpretable, and the analysis of the sample should be repeated in a new PCR run that includes positive and negative controls. Cq cut-off values should also be considered when interpreting the results.

Result interpretations and Cq cut-off values are given per target/channel in the table below.

Samples can be positive for one or more pathogens detected in the same well or in different wells.

Quality control:

Internal control

The internal control CIEZ ensures that the enzymes in the Master mix are functional. Indeed, an amplification curve needs to be observed in Cy5 channel of Master mixes A, B, C and D. The internal control value must preferentially be detected before 35 cycles ($Cq \leq 35$). Nevertheless, 2 two situations of lack of amplification of the internal control can be observed:

- If the target genes are initially present in the sample with a high number of copies, the internal control may not be amplified. This result is consistent and does not invalidate the test. It should be interpreted as a positive result despite the lack of signal from the internal control. This phenomenon is the result of amplification competition between the Internal Control and targets that are present at high copy numbers.
- If the target genes in the FAM, HEX and Texas Red channels are not amplified, as well as the internal control in the Cy5 channel, then no result can be rendered. This situation highlights the presence of PCR inhibitors or a technical problem occurring during test realization. The PCR should be repeated:
 - after diluting the obtained DNA extract at 1/10
 - or after realizing a new extraction / purification of the nucleic acids.

Control in process (CONTROL IP, yellow cap)

The exogenous Control in Process Rp49 should be detected in the Cy5 channel of Mix E. This procedural control validates the quality of the biological sample collection and validates the DNA extraction and purification steps.

In the absence of an amplification signal for the procedural control, the test is considered invalid when no target has been detected. In this case, the nucleic acids extraction or the sample treatment must be repeated.

Negative control (CONTROL -, green cap)

The fluorescence emitted must be below the threshold except for the Cy5 channel. This is an indicator of non-specific amplification. If the fluorescence is above the threshold, check for an atypical curve. In the case of an amplification curve, a contamination or a distribution error in the microtubes should be considered.

Positive control (CONTROL +, red cap)

The value of the positive control should preferably be detected before 30 cycles ($Cq \leq 30$) in the four channels (FAM, HEX, TxRed and Cy5). In the absence of amplification of the positive control, the existence of an amplification problem can be suggested. Some samples show atypical curves that are not characteristic of sigmoid amplification curves. In this case, the result should not be considered interpretable, and the sample should be analyzed again.

Master mix A				Sample interpretation
<i>Campylobacter</i> spp (16S/pxA/gyrB)	<i>Yersinia enterocolitica</i> (yst) / <i>Yersinia pseudotuberculosis</i> (ail)	<i>Salmonella</i> spp (invA)	Internal control (CIEZ)	
FAM	HEX	TxRed	Cy5	
+	-	-	+/- $Cq \leq 35$	Positive for <i>Campylobacter</i> species (<i>C. jejuni</i> , <i>C. coli</i> , <i>C. upsaliensis</i> , <i>C. fetus</i>)
-	+	-	+/- $Cq \leq 35$	Positive for <i>Yersinia enterocolitica</i> / <i>pseudotuberculosis</i>
-	-	+	+/- $Cq \leq 35$	Positive for <i>Salmonella</i> species
-	-	-	+ $Cq \leq 35$	Negative for the tested pathogens

Master mix B				Sample interpretation
<i>C. difficile</i> Hypervirulent (<i>TcdC^{Δ117}</i>)	<i>C. difficile</i> (<i>TcdB/CDT</i>)	<i>Aeromonas spp</i> (<i>gyrB</i>)	Internal control (CIEZ)	
FAM	HEX	TxRed	Cy5	
+	+	-	+/- Cq ≤ 35	Positive for <i>Clostridioides difficile</i> Hypervirulent strain
-	+	-	+/- Cq ≤ 35	Positive for <i>Clostridioides difficile</i> Toxin B or Binary Toxin
-	-	+ Cq<41	+/- Cq ≤ 35	Positive for <i>Aeromonas</i> species (<i>A. hydrophila</i> , <i>A. caviae</i> , <i>A. veronii</i> , <i>A. dhakensis</i>)
+/-	-	-	+ Cq ≤ 35	Negative for the tested pathogens

Master mix C				Sample Interpretation
Cholera toxins (<i>ctxAB/TcpA</i>)	<i>Vibrio spp</i> (<i>ATPsynthase</i>)	<i>E. coli</i> O157 strain (<i>z3276</i>)	Internal Control (CIEZ)	
FAM	HEX	TxRed	Cy5	
+	+	-	+/- Cq ≤ 35	Positive for <i>Vibrio cholerae</i> with cholera toxin genes
-	+ Cq<40	-	+/- Cq ≤ 35	Positive for <i>Vibrio</i> species (<i>V. cholerae</i> , <i>V. vulnificus</i> , <i>V. parahaemolyticus</i> , <i>V. fluvialis</i> , <i>V. alginolyticus</i>)
-	-	+	+/- Cq ≤ 35	Positive for <i>E. coli</i> O157 strain
+/-	-	-	+ Cq ≤ 35	Negative for the tested pathogens

Master mix D				Sample Interpretation
<i>E. coli</i> EIEC/ <i>Shigella spp</i> (<i>ipaH</i>)	<i>E. coli</i> EPEC (<i>bfp/eaеA</i>)	<i>E. coli</i> ETEC (<i>It1/sta</i>)	Internal Control (CIEZ)	
FAM	HEX	TxRed	Cy5	
+ Cq<40	-	-	+/- Cq ≤ 35	Positive for <i>E. coli</i> EIEC or <i>Shigella</i> species (<i>S. dysenteriae</i> , <i>S. boydii</i> , <i>S. flexneri</i> , <i>S. sonnei</i>)
-	+	-	+/- Cq ≤ 35	Positive for <i>E. coli</i> EPEC
-	-	+	+/- Cq ≤ 35	Positive for <i>E. coli</i> ETEC
-	-	-	+ Cq ≤ 35	Negative for the tested pathogens

Master mix E				Sample interpretation
<i>Plesiomonas shigelloides</i> (<i>HugA</i>)	<i>E. coli</i> EAEC (<i>aggR/aap</i>)	<i>E. coli</i> EHEC (<i>stx1/stx2</i>)	Control in process (Rp49)	
FAM	HEX	TxRed	Cy5	
+	-	-	+/- Cq ≤ 35	Positive for <i>Plesiomonas shigelloides</i>
-	+	-	+/- Cq ≤ 35	Positive for <i>E. coli</i> EAEC
-	-	+	+/- Cq ≤ 35	Positive for <i>E. coli</i> EHEC
-	-	-	+ Cq ≤ 35	Negative for the tested pathogens

As the E coli pathotypes are detected in different wells, result interpretation in wells C, D and E is required. The table below is an aid for a prompt interpretation.

Master mix C			Master mix D			Master mix E			Sample interpretation
Cholera toxins (ctxAB/TcpA)	<i>Vibrio spp</i> (ATPsynthase)	<i>E. coli</i> O157 strain (z3276)	<i>E. coli</i> EIEC / <i>Shigella spp</i> (ipaH)	<i>E. coli</i> EPEC (bfp/eaecA)	<i>E. coli</i> ETEC (lt1/sta)	<i>Plesiomonas shigelloides</i> (HugA)	<i>E. coli</i> EAEC (aggR/aap)	<i>E. coli</i> EHEC (stx1/stx2)	
FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	
-	-	-	+	-	-	-	-	-	EIEC and/or <i>Shigella spp</i> (<i>S. boydii</i> , <i>S. sonnei</i> , <i>S. dysenteriae</i> , <i>S. flexneri</i>)
-	-	-	+	-	-	-	-	+	<i>Shigella dysenteriae</i> / <i>Shigella flexneri</i> or co-infection EIEC and EHEC or <i>Shigella spp</i> and EHEC
-	-	-	-	-	-	-	+	-	Typical or atypical EAEC
-	-	-	-	-	+	-	-	-	ETEC
-	-	-	-	+	-	-	-	-	EPEC
-	-	-	-	-	-	-	-	+	EHEC
-	-	+	-	-	-	-	-	+	EHEC O157 strain
-	-	-	-	+	-	-	-	+	STEC (EHEC) or co-infection EPEC and EHEC
-	-	+	-	+	-	-	-	+	STEC (EHEC O157) or co-infection EPEC and EHEC O157 strain
-	-	+	-	-	-	-	-	-	Negative*

*The BIOSYNEX AMPLIQUICK Fecal Bacteriology kit can detect the O157 serotype of *E. coli* EHEC that is responsible of major outbreaks. The virulence of the EHEC pathotype is mediated by the secretion of shigatoxins (stx1/stx2). Several findings have documented the presence of non-pathogenic *E. coli* O157 strains in feces. These non-pathogenic O157 strains do not carry stx1/stx2 genes and are not associated with gastroenteritis.^{20,21} Therefore, a positive detection of *E. coli* O157 should be interpreted as negative results in the case of negative detection of EHEC (stx1/stx2). The detection results of O157 characterization should be interpreted when only the EHEC is found positive.

10. PERFORMANCES

Analytical sensitivity

Target sequence detection limits:

The detection limit of the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit is defined as the concentration, in number of copies/μL, that can be detected at 95% in a DNA sample specific to each target. It has been determined by performing serial dilution of reference samples with a known number of copies.

The LoD₉₅ of the BIOSYNEX AMPLIQUICK Fecal Bacteriology on DNA samples was statistically determined for each pathogen and targeted gene. Results are given in the following table:

Target	Tested strains	LoD ₉₅ (DNA copy/μL CFU/μL)	LoD ₉₅ (DNA copy/rxn CFU/rxn)
<i>Campylobacter jejuni</i>	Vircell MBC088-R	0.579	2.895
<i>Campylobacter coli</i>	ATCC RM 2228	0.245	1.225
<i>Campylobacter upsaliensis</i>	ATCC RM 3195	1.264	6.32
<i>Campylobacter fetus</i>	ATCC NCTC 10842	1.978	9.89
<i>Yersinia enterocolitica</i>	Vircell MBC027	2.221	11.105
<i>Yersinia pseudotuberculosis</i>	ATCC NCTC 2476	1.167	5.835
<i>Salmonella typhimurium</i>	Vircell MBC044	3.894	19.47
<i>Salmonella enteritidis</i>	Vircell MBC003	3.076	15.38
<i>Clostridium difficile</i> Strain 4118 (Hypervirulent detection)	ATCC 4118	1.092	5.46
<i>Clostridium difficile</i> Strain 4118 (Toxin detection)	ATCC 4118	0.432	2.16
<i>Clostridium difficile</i>	Vircell MBC043	3.607	18.035
<i>Aeromonas veronii</i>	DSMZ 7386	2.562	12.81
<i>Aeromonas dhakensis</i>	DSMZ 18362	4.372	21.86

<i>Aeromonas hydrophilia</i>	DSMZ 30187	2.523	12.615
<i>Aeromonas caviae</i>	DSMZ 7323	1.122	5.61
<i>Vibrio cholerae</i> (Cholera toxin detection)	Vircell MBC118	3.416	17.08
<i>Vibrio cholerae</i> (<i>Vibrio spp</i> detection)	Vircell MBC118	4.345	21.725
<i>Vibrio parahaemolyticus</i>	ATCC EB101	2.981	14.905
<i>Vibrio vulnificus</i>	DSMZ 10143	3.328	16.64
<i>E. coli</i> VTEC (O157 serotype detection)	Vircell MBC022	3.570	17.85
<i>E. coli</i> EIEC	Vircell MBC122	2.461	12.305
<i>Shigella flexneri</i>	DSMZ 4782	0.811	4.055
<i>E. coli</i> EPEC	Vircell MBC123	0.722	3.61
<i>E. coli</i> ETEC	Vircell MBC124	5.428	27.14
<i>Plesiomonas shigelloides</i>	ATCC GNI14	9.221	46.105
<i>E. coli</i> EAEC	Vircell MBC121	2.225	11.125
<i>E. coli</i> VTEC (EHEC <i>stx1/stx2</i> detection Mix E/TxRed)	Vircell MBC022	1.820	9.10

Analytical specificity

The design of oligonucleotides was validated *in silico* by BLAST alignment. The comparison of the primer and probe sequences shows that the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit is specific to the detection of the specified targets.

Cross reactivity

A panel of 42 RNA and 72 DNA samples from a biobank, listed in the following tables, were tested using the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit. For all these samples, no amplification of targets in the various master mixes was observed.

RNA		
Astrovirus	Human Parainfluenza 1	Parainfluenza 3
Chikungunya virus	Influenza A H1	Parainfluenza 4 A
Coronavirus	Influenza A H3	Respiratory syncytial virus (subtype A)
Coronavirus Oc43	Influenza A H5	Respiratory syncytial virus (subtype B)
Coronavirus SARS	Influenza B	Rhinovirus
Coxsackie A6	Measles	Rotavirus
Coxsackie B1	MERS Coronavirus	Rubella
Coxsackie B5	Mumps	Sapovirus
Dengue 1	Norovirus	SARS-CoV-2
Dengue 2	Norovirus G1	Tick-Borne Encephalitis virus
Dengue 3	Norovirus G2	West Nile virus
Dengue 4	Novel Influenza A H1N1	Yellow Fever virus
Echovirus 5	Parainfluenza 1	Zika virus
Enterovirus 68	Parainfluenza 2	Zika virus (asian lineage)

DNA		
<i>Acanthamoeba castellanii</i>	<i>Coccidioides immitis</i>	<i>Mycobacterium kansasii</i>
Adenovirus	<i>Coxiella burnetii</i>	<i>Mycobacterium tuberculosis</i>
Adenovirus type 02	<i>Cryptosporidium parvum</i>	<i>Mycobacterium ulcerans</i>
Adenovirus type 04	Cytomegalovirus	<i>Mycoplasma genitalium</i>
Adenovirus type 11	<i>Enterococcus faecalis</i> (vanB)	<i>Mycoplasma hominis</i>
Adenovirus type 20	<i>Enterococcus faecium</i> (vanA)	<i>Neisseria gonorrhoeae</i>
Adenovirus type 31	Epstein-Barr virus	<i>Neisseria meningitidis</i> sg A
Adenovirus type 40	<i>Francisella tularensis</i>	Papillomavirus type 16
Adenovirus 41	<i>Gardnerella vaginalis</i>	Papillomavirus type 18
<i>Aspergillus fumigatus</i>	<i>Giardia intestinalis</i>	Parvovirus b19 (plasmid)
<i>Bacillus cereus</i>	<i>Haemophilus ducreyi</i>	<i>Plasmodium falciparum</i>
<i>Bartonella henselae</i>	<i>Haemophilus influenzae</i>	<i>Rickettsia conorii</i>
<i>Bartonella quintana</i>	<i>Helicobacter pylori</i>	<i>Staphylococcus aureus</i> (mecA-)
Bk virus	Herpes simplex 1	<i>Staphylococcus aureus</i> (mecA+)
<i>Bordetella holmesii</i>	Herpes simplex 2	<i>Streptococcus pneumoniae</i>
<i>Borrelia afzelii</i>	HHV-6	<i>Toxoplasma gondii</i>

<i>Borrelia burgdorferi</i>	HHV-8	<i>Treponema pallidum</i>
<i>Borrelia garinii</i>	<i>Klebsiella pneumoniae</i> (NDM-1)	<i>Trichomonas vaginalis</i>
<i>Brucella abortus</i>	<i>Legionella pneumophila</i>	<i>Trypanosoma cruzi</i>
<i>Candida albicans</i>	<i>Leishmania chagasi</i>	<i>Ureaplasma parvum</i>
<i>Candida auris</i>	<i>Leishmania infantum</i>	<i>Ureaplasma urealyticum</i>
<i>Chlamydia trachomatis</i>	<i>Listeria monocytogenes</i>	Varicella-zoster virus
<i>Chlamydophila pneumoniae</i>	<i>Moraxella catarrhalis</i>	-
<i>Chlamydophila psittaci</i>	<i>Mycobacterium avium</i>	-
<i>Clostridioides difficile</i> ribotype 010	<i>Mycobacterium intracellulare</i>	-

Analytical reactivity

The analytical reactivity of the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit was confirmed by testing a panel of 51 DNA of different bacterial strains targeted by this kit. The results are given in Annex 1.

Chemical interference

Various substances that may be present in patients' stool samples may cause positive or negative interference on PCR results. Interference tests on the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit were performed by adding the potential interferent to fecal swab positive samples and tested on the BIOSYNEX AMPLIQUICK Fecal Bacteriology:

Tested substances	Tested concentration	Interference observed
Betamethasone	2.5% v/v	None
Hemorrhoid cream (Titanoréine crème)	0.25% v/v	
Mucin	0.8% w/v	
Imodium (Lopéramide-Lyoc)	5% w/v	
Ampicillin sodium salt	152 µmol/L	
Fatty acid supplement	4.8%	
Pepto-Bismol	5%	
Starch	1.5% w/v	
Cellulose	4% w/v	
Pectin	0.01% w/v	
Lubricant	5% v/v	
Ethanol	0.2% v/v	
Urine	20% v/v	
Blood	40% v/v	

Precision data

The precision data for the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit were determined based on 5 conditions:

- Intra-assay variation (within one test run)
- Inter-day variation (between different test runs)
- Inter-operator variation (between different operators)
- Inter-laboratory variation (between different locations)
- Inter-batch variation (between different batches)

The precision data are expressed in terms of mean Cq value, standard deviation (SD) and coefficient of variation (CV), based on the threshold quantification cycle (Cq) values for the DNA of the various targets. Intra-assay variability data for the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit was determined on 30 replicate tests using negative and positive bacterial samples for each bacterial target. The SD and CV determinations showed that the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit is repeatable with CV<5%. Inter-days, inter-operators, inter-laboratories and inter-batches variabilities were determined on duplicate tests from extracted bacterial DNA. The BIOSYNEX AMPLIQUICK Fecal Bacteriology results are reproducible with CV<5%.

The variability per target is given in the tables below:

<i>C. jejuni</i> variability	Master mix A / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	24.71	0.09	0.37%	36.47	0.98	2.70%
Inter-days	19.77	0.50	2.51%	36.10	0.85	2.35
Inter-operators	19.43	0.57	2.93%	35.53	0.55	1.54
Inter-laboratories	19.85	0.77	3.87%	36.13	1.43	3.95%
Inter-batches	19.93	0.06	0.3%	35.31	0.18	0.5%

<i>C. fetus</i> variability	Master mix A / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	24.31	0.08	0.33%	33.85	0.22	0.64%
Inter-days	20.28	0.35	1.72%	36.53	0.84	2.30%
Inter-operators	20.42	0.34	1.64%	36.98	0.17	0.45%
Inter-laboratories	20.46	0.31	1.54%	35.46	0.20	0.57%
Inter-batches	20.34	0.35	1.72%	36.03	0.63	1.76%

<i>C. upsaliensis</i> variability	Master mix A / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	24.17	0.11	0.46%	34.55	0.31	0.89%
Inter-days	19.75	0.28	1.41%	35.73	0.52	1.44%
Inter-operators	19.75	0.18	0.91%	36.34	0.35	0.95%
Inter-laboratories	19.33	0.12	0.60%	35.20	1.32	3.75%
Inter-batches	20.59	0.80	3.88%	35.70	0.35	0.98%

<i>Y. enterocolitica</i> variability	Master mix A / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	26.38	0.06	0.21%	35.78	0.64	1.79%
Inter-days	20.43	0.41	2.03%	35.58	0.33	0.94%
Inter-operators	20.28	0.25	1.22%	35.75	0.31	0.86%
Inter-laboratories	20.85	0.33	1.59%	35.58	0.28	0.79%
Inter-batches	20.24	0.52	2.59%	35.03	0.26	0.75%

<i>Y. pseudotuberculosis</i> variability	Master mix A / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	23.57	0.06	0.26%	34.51	0.90	2.61%
Inter-days	18.82	0.41	2.19%	35.21	0.20	0.56%
Inter-operators	18.82	0.11	0.60%	35.09	0.10	0.27%
Inter-laboratories	19.17	0.51	2.64%	35.54	0.01	0.02%
Inter-batches	18.45	0.23	1.23%	35.13	0.00	0.01%

<i>S. enteritidis</i> variability	Master mix A / Texas Red					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.31	0.08	0.30%	35.98	0.47	1.31%
Inter-days	19.91	0.72	3.60%	35.94	0.89	2.47%
Inter-operators	19.54	0.59	3.02%	36.55	0.28	0.77%
Inter-laboratories	20.30	0.16	0.78%	36.14	0.11	0.31%
Inter-batches	19.95	0.6	3.03%	35.28	0.65	1.85%

<i>C. difficile</i> hypervirulent variability	Master mix B / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	25.24	0.06	0.22%	34.17	0.59	1.74%
Inter-days	19.43	0.54	2.76%	35.37	0.28	0.79%
Inter-operators	19.10	0.56	2.93%	35.47	0.20	0.56%
Inter-laboratories	19.62	0.82	4.16%	34.77	1.24	3.56%
Inter-batches	19.57	0.12	0.60%	35.01	0.14	0.39%

<i>C. difficile</i> toxins variability	Master mix B / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	25.47	0.07	0.29%	34.89	0.50	1.44%
Inter-days	19.79	0.40	2.04%	35.51	0.45	1.28%
Inter-operators	19.64	0.29	1.46%	35.28	0.25	0.71%
Inter-laboratories	20.22	0.23	1.14%	36.55	0.74	2.02%
Inter-batches	19.84	0.24	1.23%	34.84	0.35	0.99%

<i>A. veronii</i> variability	Master mix B / Texas Red					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	25.78	0.06	0.22%	34.86	0.37	1.07%
Inter-days	20.29	0.57	2.81%	36.72	0.97	2.63%
Inter-operators	19.96	0.55	2.76%	36.81	1.40	3.79%
Inter-laboratories	20.81	0.05	0.25%	36.79	0.27	0.74%
Inter-batches	20.44	0.64	3.11%	35.79	0.20	0.55%

Cholera toxin variability	Master mix C / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.24	0.09	0.34%	36.34	0.88	2.41%
Inter-days	20.15	0.54	2.70%	36.54	0.74	2.03%
Inter-operators	19.78	0.65	3.27%	36.20	0.71	1.96%
Inter-laboratories	20.29	0.77	3.78%	36.10	0.72	2.01%
Inter-batches	20.18	0.05	0.25%	36.20	0.41	1.14%

<i>Vibrio spp</i> variability	Master mix C / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.21	0.08	0.30%	35.58	0.85	2.38%
Inter-days	20.36	0.38	1.84%	35.62	0.64	1.80%
Inter-operators	20.22	0.25	1.24%	36.39	1.01	2.77%
Inter-laboratories	20.76	0.25	1.21%	36.25	1.52	4.20%
Inter-batches	20.34	0.28	1.39%	35.36	0.42	1.19%

EHEC O157 variability	Master mix C / Texas Red					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.00	0.07	0.23%	35.23	0.42	1.19%
Inter-days	21.40	0.63	2.94%	37.62	0.47	1.26%
Inter-operators	21.10	0.68	3.22%	37.82	0.41	1.09%
Inter-laboratories	21.89	0.07	0.31%	37.24	0.46	1.22%
Inter-batches	21.45	0.41	1.93%	37.00	0.13	0.36%

EIEC/ <i>Shigella spp</i> variability	Master mix D / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	26.00	0.45	1.73%	37.31	1.23	3.29%
Inter-days	20.27	0.35	1.74%	36.63	0.74	2.03%
Inter-operators	20.44	0.48	2.37%	36.04	0.68	1.88%
Inter-laboratories	20.03	0.63	3.16%	36.29	1.62	4.46%
Inter-batches	20.00	0.03	0.16%	36.01	0.30	0.84%

EPEC variability	Master mix D / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	26.45	0.07	0.26%	37.88	0.74	1.94%
Inter-days	20.16	0.35	1.72%	35.45	0.78	2.21%
Inter-operators	19.99	0.29	1.47%	35.11	0.68	1.94%
Inter-laboratories	20.51	0.22	1.05%	35.89	0.72	2.00%
Inter-batches	20.16	0.27	1.35%	34.94	0.46	1.31%

ETEC variability	Master mix D / Texas Red					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.62	0.18	0.66%	35.86	0.86	2.41%
Inter-days	19.60	0.42	2.15%	35.53	0.48	1.34%
Inter-operators	19.19	0.09	0.48%	36.12	0.92	2.55%
Inter-laboratories	19.31	0.16	0.81%	36.25	0.07	0.19%
Inter-batches	18.90	0.30	1.57%	35.48	0.01	0.02%

<i>P. shigelloides</i> variability	Master mix E / FAM					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	26.55	0.06	0.23%	32.24	0.12	0.38%
Inter-days	19.33	0.47	2.43%	35.79	0.78	2.17%
Inter-operators	19.53	0.57	2.93%	36.03	0.92	2.54%
Inter-laboratories	19.05	0.62	3.23%	35.84	0.92	2.57%
Inter-batches	18.79	0.29	1.52%	35.14	0.44	1.26%

EAEC variability	Master mix E / HEX					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	27.20	0.07	0.27%	35.61	0.52	1.47%
Inter-days	20.13	0.34	1.67%	35.15	0.79	2.25%
Inter-operators	19.99	0.34	1.68%	35.87	0.08	0.22%
Inter-laboratories	20.51	0.19	0.95%	35.72	0.09	0.25%
Inter-batches	19.98	0.35	1.77%	34.83	0.41	1.18%

EHEC O157 variability	Master mix E / Texas Red					
	"High" Sample Concentration			"Low" Sample Concentration		
	Mean Cq	SD	CV	Mean Cq	SD	CV
Intra-assay	28.05	0.06	0.23%	36.39	0.46	1.27%
Inter-days	19.35	0.51	2.65%	35.87	0.58	1.61%
Inter-operators	21.20	0.43	2.04%	37.03	0.04	0.11%
Inter-laboratories	19.98	0.06	0.32%	35.97	0.26	0.73%
Inter-batches	19.60	0.10	0.49%	35.95	0.37	1.03%

Clinical performance

The clinical performance of the BIOSYNEX AMPLIQUICK Fecal Bacteriology kit was determined using 401 stool samples including positive and negative samples. The study was conducted at the French National Reference Center of Campylobacter (CNR des Campylobacters et Hélicobacters / CHU Bordeaux, France).

Samples used in this study were qualified using three PCR enteric bacterial panel kits that are CE-marked and/or bacterial cultures by the National Reference Center of Campylobacters.

Three different PCR kits from two different manufacturers were used as reference PCRs to cover all targets in the AMPLIQUICK panel.

The contingency table below shows the detection performance of the BIOSYNEX AMPLIQUICK Fecal Bacteriology test, in comparison with the reference PCR.

Sample qualification		Reference PCR	
		Positives	Negatives
BIOSYNEX AMPLIQUICK Fecal Bacteriology	Positives	344	1 (+1)*
	Negatives	2*	54

*One sample was found falsely negative for EHEC (stx1/stx2) and falsely positive for *Salmonella spp.*

The clinical sensitivity and specificity were calculated per bacterial target. The clinical sensitivity of five targets (Cholera toxin, *Vibrio spp.*, *C. difficile* hypervirulent strain, ETEC and *Plesiomonas shigelloides*) were not calculated as the number of positive tested samples is very low. Results are given in the table below.

Master Mix	Target	True Positive	True Negative	False Positive	False Negative	Sensitivity (95% CI)	Specificity (95% CI)
A	<i>Campylobacter spp</i>	196	205	0	0	100% (98.14-100%)	100% (98.22-100%)
	<i>Yersinia enterocolitica/Yersinia pseudotuberculosis</i>	22	105	0	0	100% (84.56-100%)	100% (96.55-100%)
	<i>Salmonella spp</i>	56	344	1	0	100% (93.62-100%)	99.71% (98.40-99.99%)
B	<i>C. difficile</i> Hypervirulent strain	0	102	0	0	-	100% (96.45-100%)
	<i>C. difficile</i> TcdB/CDT	11	102	0	0	100% (71.51-100%)	100% (96.45-100%)
	<i>Aeromonas spp</i>	35	63	0	0	100% (90.00-100%)	100% (94.31-100%)
C	Cholera toxins	0	110	0	0	-	100% (96.70-100%)
	<i>Vibrio spp</i>	2	110	0	0	-	100% (96.70-100%)
D	EIEC/ <i>Shigella spp</i>	39	362	0	0	100% (90.97-100%)	100% (98.99-100%)
	EPEC	46	37	0	0	100% (92.29-100%)	100% (90.51-100%)
	ETEC	8	74	1	0	-	98% (92.79-99.97%)
E	EAEC	36	47	0	0	100% (90.26-100%)	100% (92.45-100%)
	EHEC	32	319	0	2	94.12% (80.32-99.28%)	100% (98.85-100%)
	<i>Plesiomonas shigelloides</i>	2	64	0	0	-	100% (94.40-100%)


11. LIMITATIONS

- Failure to follow any instructions of the IFU may adversely affect the test performance and/or invalidate the test result.
- This test cannot differentiate between *Yersinia enterocolitica* and *Yersinia pseudotuberculosis* as they are detected in the same fluorescent channel.
- This test cannot differentiate between *C. difficile* strains carrying the sequence coding the Toxin B and/or coding the Binary toxin as they are detected in the same fluorescent channel.
- The possibility to obtain false positive or false negative results cannot be ruled out (e.g. in case of accidental contamination, in case of poor quality of the sample or in case the test procedure is not correctly followed). Hence, as with all diagnostic tests, results must be considered with other clinical information available to the physician.
- Liquid stool samples may induce false-negative results due to dilution. In case of liquid samples, we recommend performing the extraction of stool sample directly by adding 200µL of sample.

12. LITERATURE

1. Graves, N.S. (2013). Acute Gastroenteritis. Primary Care: Clinics in Office Practice 40, 727–741. 10.1016/j.pop.2013.05.006.
2. Bányai, K., Estes, M.K., Martella, V., and Parashar, U.D. (2018). Viral gastroenteritis. The Lancet 392, 175–186. 10.1016/S0140-6736(18)31128-0.
3. Centers for Disease Control and Prevention (CDC). “Outbreak of Salmonella Infections Linked to Backyard Poultry.” 2020.
4. European Food Safety Authority (EFSA). “The European Union summary report on trends and sources of zoonoses, zoonotic agents and food-borne outbreaks in 2017.” EFSA Journal, 2018; 16(12):5500.
5. Brenner FW, Villar RG, Angulo FJ, Tauxe R, Swaminathan B. 2000. Salmonella Nomenclature. J Clin Microbiol 38: 2011;2011:735308. doi: 10.1128/jcm.38.7.2465-2467.2000.
6. European Centre for Disease Prevention and Control. Healthcare-associated infections acquired in intensive care units. In: ECDC. Annual Epidemiological Report for 2020. Stockholm: ECDC; 2024.
7. Fukushima H, Shimizu S, Inatsu Y. Yersinia enterocolitica and Yersinia pseudotuberculosis Detection in Foods. J Pathog. 2011;2011:735308. doi: 10.4061/2011/735308. Epub 2011 Oct 5. PMID: 22567341; PMCID: PMC3335482.
8. Kaakoush NO, Castaño-Rodríguez N, Mitchell HM, Man SM. Global Epidemiology of Campylobacter Infection. Clin Microbiol Rev. 2015 Jul;28(3):687-720. doi: 10.1128/CMR.00006-15. PMID: 26062576; PMCID: PMC4462680.
9. World Health Organization (WHO). “Campylobacter.” Fact sheet, 2020.
10. Janda JM, Abbott SL. The genus Aeromonas: taxonomy, pathogenicity, and infection. Clin Microbiol Rev. 2010 Jan;23(1):35-73. doi: 10.1128/CMR.00039-09. PMID: 20065325; PMCID: PMC2806660.
11. Yuwono C, Wehrhahn MC, Liu F, Zhang L. 2023. Enteric Aeromonas Infection: a Common Enteric Bacterial Infection with a Novel Infection Pattern Detected in an Australian Population with Gastroenteritis. Microbiol Spectr 11:e00286-23. https://doi.org/10.1128/spectrum.00286-23.
12. Guh AY, Kuty PK. Ann Intern Med. 2018 October 02; 169(7)
13. Spigaglia P et al. Clostridium difficile causing pediatric infections: New findings from a hospital-based study in Italy. Anaerobe. 2017 Dec;48:262-268. doi: 10.1016/j.anaerobe.2017.10.008.
14. European Centre for Disease Prevention and Control. Antimicrobial resistance surveillance in Europe 2011. Annual Report of the European Antimicrobial Resistance Surveillance Network (EARS-Net). Stockholm: ECDC; 2012.
15. Janda JM, Abbott SL, McIver CJ. Plesiomonas shigelloides Revisited. Clin Microbiol Rev. 2016 Apr;29(2):349-74. doi: 10.1128/CMR.00103-15. PMID: 26960939; PMCID: PMC4786884.
16. World Health Organization (WHO). “Cholera.” Fact sheet, 2020.
17. Baker-Austin, C., Stockley, L., Rangdale, R. and Martinez-Urtaza, J. (2010), Environmental occurrence and clinical impact of Vibrio vulnificus and Vibrio parahaemolyticus: a European perspective. Environmental Microbiology Reports, 2: 7-18. https://doi.org/10.1111/j.1758-2229.2009.00096.x
18. Pakbin B, Brück WM, Rossen JWA. Virulence Factors of Enteric Pathogenic Escherichia coli: A Review. Int J Mol Sci. 2021 Sep 14;22(18):9922. doi: 10.3390/ijms22189922. PMID: 34576083; PMCID: PMC8468683.
19. Shigella species and characteristics: “Shigellose : symptômes, traitement, prévention.” Institut Pasteur.
20. Schmidt, H., Scheef, J., Huppertz, H.I., Frosch, M., and Karch, H. (1999). Escherichia coli O157:H7 and O157:Hx Strains That Do Not Produce Shiga Toxin: Phenotypic and Genetic Characterization of Isolates Associated with Diarrhea and Hemolytic-Uremic Syndrome. J. CLIN. MICROBIOL. 37.
21. Stephan, R., Zhang, W., Bielaszewska, M., Mellmann, A., and Karch, H. (2009). Phenotypic and Genotypic Traits of Shiga Toxin–Negative E. coli O157:H7/Hx Bovine and Porcine Strains. Foodborne Pathogens and Disease 6, 235–243. https://doi.org/10.1089/fpd.2008.0205.

13. INTERNATIONAL SYMBOLS

	Consult instructions for use or electronic instructions for use		Contains sufficient for <n> tests
	Catalogue number		<i>In vitro</i> diagnostic medical device
	Do not reuse		Bag of optical cap strips
	Temperature limits		Use-by-date
	Manufacturer		Batch code
	Keep away from sunlight		Importer
	Authorized representative in Switzerland		Do not use if package is damaged and consult the instructions for use
	Unique device identifier		Master Mix
	Positive control		Negative control
	Control in process		Master mix control
	Prefilled microplates format		Bulk tubes format
	Qualitative detection by PCR		

14. MANUFACTURER INFORMATION

BIOSYNEX S.A.

22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN – France
Tel : +33 3 88 78 78 87
www.biosynex.com

Contact France:


Tel.: +33 3 88 77 57 00
service.clients@biosynex.com

Contact other countries :

Tel.: +33 3 88 77 57 52
sales@biosynex.com

After-sales service:

Tel: +33 3 88 77 57 25
Tech.support@biosynex.com

BIOSYNEX SWISS S.A.

Route de Rossemaison 100
2800 DELEMONT - Switzerland

15. HISTORIC MODIFICATION

Record version	Amended paragraph	Change details
V1	All	Creation of the IFU

ANNEXE 1

Bacterial strains	Master mix A			Master mix B			Master mix C			Master mix D			Master mix E		
	<i>Campylobacter spp</i>	<i>Yersinia spp</i>	<i>Salmonella spp</i>	<i>C.difficile Hypervirulent</i>	<i>C.difficile TcdB/CDT</i>	<i>Aeromonas spp</i>	Cholera toxins	<i>Vibrio spp</i>	<i>E.coli O157</i>	<i>E.coli EIEC/Shigella spp</i>	<i>E.coli EPEC</i>	<i>E.coli ETEC</i>	<i>Plesiomonas shigelloides</i>	<i>E.coli EAEC</i>	<i>E.coli EHEC</i>
	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)
<i>Escherichia coli</i> (EAEC) MBC121 (VIRCELL)	-	-	-	-	-	-	-	-	-	-	-	-	-	Detected	-
<i>Escherichia coli</i> (EIEC) MBC122 (VIRCELL)	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Escherichia coli</i> (EPEC) MBC123 (VIRCELL)	-	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-
<i>Escherichia coli</i> (ETEC) MBC124 (VIRCELL)	-	-	-	-	-	-	-	-	-	-	-	Detected	-	-	-
<i>Escherichia coli</i> (VTEC) MBC022 (VIRCELL)	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-	Detected
<i>Escherichia coli</i> O111/NM (Zeptomatrix)	-	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-
<i>Shigella sonnei</i> DSMZ 5570	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Shigella boydii</i> DSMZ 7532	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Shigella flexneri</i> DSMZ 4782	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Shigella dysenteriae</i> DSMZ 103303	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Shigella flexneri</i> MBC089 (VIRCELL)	-	-	-	-	-	-	-	-	-	Detected	-	-	-	-	-
<i>Campylobacter upsaliensis</i> RM 3195 (ATCC)	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Campylobacter coli</i> RM 2228 (ATCC)	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Campylobacter fetus</i> NCTC 10842 (ATCC)	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Campylobacter jejuni</i> MBC088 (VIRCELL)	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Clostridium difficile</i> MBC043 (VIRCELL)	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 001	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 010	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 014	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 020	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 027	-	-	-	Detected	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 033	-	-	-	Detected	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> Ribotype 078	-	-	-	Detected	Detected	-	-	-	-	-	-	-	-	-	-
<i>Clostridioides difficile</i> hypervirulent Strain 4118 (ATCC)	-	-	-	Detected	Detected	-	-	-	-	-	-	-	-	-	-
<i>Aeromonas caviae</i> DSMZ 7323	-	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-
<i>Aeromonas dhakensis</i> DSMZ 18362	-	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-
<i>Aeromonas veronii</i> DSMZ 7386	-	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-
<i>Aeromonas hydrophila</i> DSMZ 30187	-	-	-	-	-	Detected	-	-	-	-	-	-	-	-	-
<i>Yersinia pseudotuberculosis</i> NCTC 2476 (ATCC)	-	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-
<i>Yersinia enterocolitica</i> MBC027 (VIRCELL)	-	Detected	-	-	-	-	-	-	-	-	-	-	-	-	-

Bacterial strains	Master mix A			Master mix B			Master mix C			Master mix D			Master mix E		
	<i>Campylobacter spp</i>	<i>Yersinia spp</i>	<i>Salmonella spp</i>	<i>C.difficile Hypervirulent</i>	<i>C.difficile TcdB/CDT</i>	<i>Aeromonas spp</i>	Cholera toxins	<i>Vibrio spp</i>	<i>E.coli O157</i>	<i>E.coli EIEC/Shigella spp</i>	<i>E.coli EPEC</i>	<i>E.coli ETEC</i>	<i>Plesiomonas shigelloides</i>	<i>E.coli EAEC</i>	<i>E.coli EHEC</i>
	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)	(FAM)	(HEX)	(TxRed)
<i>Salmonella enteritidis</i> MBC003 (VIRCELL)	-	-	Detected	-	-	-	-	-	-	-	-	-	-	-	-
<i>Salmonella typhi</i> MBC044 (VIRCELL)	-	-	Detected	-	-	-	-	-	-	-	-	-	-	-	-
<i>Plesiomonas shigelloides</i> GNI14 (ATCC)	-	-	-	-	-	-	-	-	-	-	-	-	Detected	-	-
<i>Vibrio parahaemolyticus</i> EB101 (ATCC)	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio vulnificus</i> DSMZ 10143	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> MBC118 (VIRCELL)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O1 EL Tor ctxB1 (tcpA El Tor)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O1 El Tor ctxB3 (tcpA El Tor)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O1 El Tor ctxB7 (tcpA El Tor)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O1 classical (ctxB1, tcp1 classical)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O139	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O37 (tcpA classical)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> O141 (ctxB, tcpA classical)	-	-	-	-	-	-	Detected	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> non-O1/non-O139 (ctxB & tcpA neg)	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio cholerae</i> non-O1/non-O139 (ctxB & tcpA neg)	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio vulnificus</i> Strain 1	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio vulnificus</i> Strain 2	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio parahaemolyticus</i> Strain 1	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio parahaemolyticus</i> Strain 2	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio alginolyticus</i>	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-
<i>Vibrio fluvialis</i>	-	-	-	-	-	-	-	Detected	-	-	-	-	-	-	-