



Molzylm

Add-On 10

Add-on kit for the extraction of pathogen DNA for volumes up to 10 ml



Key features

- Larger volume protocol (up to 10 ml) compatible with Sepsitest™, UMD-Universal, Sepsitest™ SelectNA and UMD SelectNA
- Various body fluids (e.g. whole blood, ascites, urine)
- Removal of human DNA and free floating DNA from dead bacteria

Kit	Application	Content	Order No.
Add-On 10	Add-On kit for the extraction of pathogen DNA for volumes (up to 10ml)	24 extractions	U-120-024
		48 extractions	U-120-048