

FlashDx[®] Feline GI Panel 1.0

RUO



| | |
|----------------------------|---|
| Test Reagent Kit | Feline GI Panel 1.0 |
| Catalog Number | FEL-P2-10 |
| Instrument | FlashDx-1000-E |
| Technology | SPRT™ (semi-Solid Phase Real-Time PCR) technologies |
| Targets | FPV: VP2 gene FCoV: N gene |
| Batch or On-Demand | On-Demand |
| Minimum Batch Size | 1 |
| Sample Types | Anal/Rectal swab Transport Media: UTM/VTM* or saline (0.9%) |
| Sample Extraction | Fully automated/integrated |
| Precision Pipetting | Not required |
| Turnaround Time | <1 hour |
| Hands-on Time | ~1 minute |
| Control | Internal Process Control |
| Kit Storage | 2-8 °C |
| Commercial Controls | See FlashDx Package Insert or contact FlashDx Customer Support |

* Sample transport media with guanidine buffer is not compatible.