

# FlashDx<sup>®</sup> Group B Strep (GBS) Test



<b>Test Reagent Kit</b>	GBS Test
<b>Catalog Number</b>	GBS-P1-10
<b>Instrument</b>	FlashDx-1000-E, SPRT™ 4T/8T/12T
<b>Technology</b>	SPRT™ (semi-Solid Phase Real-Time PCR) technologies
<b>Targets</b>	<i>cfb/sip</i> gene
<b>Batch or On-Demand</b>	On-Demand
<b>Minimum Batch Size</b>	1
<b>Sample Types</b>	<b>Specimen:</b> Vaginal swabs <b>Transport Media:</b> UTM* or saline (0.9%)
<b>Sample Extraction</b>	Fully automated/integrated
<b>Precision Pipetting</b>	Not required
<b>Turnaround Time</b>	<1 hour
<b>Hands-on Time</b>	~1 minute
<b>Control</b>	Internal Process Control
<b>Clinical Evaluation</b>	100% sensitivity   100% specificity
<b>Kit Storage</b>	2-8 °C
<b>Commercial Controls</b>	See FlashDx GBS Package Insert or contact FlashDx Customer Support

\*Sample transport media with guanidine buffer is not compatible.

