



Variant List – Chimeric ID Panel

SNP RS#	Chromosome	Chromosome Position	Average Heterozygosity
rs1035271	chr2	143903354	0.497959
rs1044910	chr12	96953146	0.499234
rs10454068	chr16	61073636	0.497558
rs10744034	chr12	7727955	0.491425
rs10771010	chr12	23616515	0.49975
rs10898438	chr11	86147211	0.493201
rs10956179	chr8	124418519	0.498871
rs1105359	chr3	59423798	0.493154
rs11136962	chr8	6058982	0.499908
rs11670999	chr19	29239107	0.495252
rs11671916	chr19	4304501	0.49917
rs11989001	chr8	61918958	0.499017
rs11989433	chr8	40685533	0.499053
rs1206523	chr6	25543493	0.496746
rs12086717	chr1	230120112	0.492775
rs12313812	chr12	60514576	0.497984
rs12318959	chr12	28629032	0.499971
rs12872438	chr13	83208497	0.498964
rs12919933	chr16	3087260	0.499839
rs13140054	chr4	146288252	0.499732
rs1387566	chr2	10688524	0.499551
rs1446940	chr8	141842656	0.499295
rs1451509	chr17	59365772	0.499053
rs1481847	chr8	71869260	0.49655
rs153749	chr5	171754251	0.49889
rs1854853	chr6	167119574	0.497558
rs1893023	chr11	121133442	0.5
rs1893673	chr18	57660068	0.49753
rs1958693	chr14	37875170	0.499368
rs1971391	chr5	10512218	0.499928
rs2026263	chr21	37338729	0.495781
rs217297	chr6	83671204	0.499999
rs2202828	chr9	103527145	0.499154
rs2236067	chr14	65781294	0.499824
rs2255054	chr1	218838908	0.492871
rs2256763	chr2	193476588	0.498982

Variant List – Chimeric ID Panel (PAGE 2 OF 3)

SNP RS#	Chromosome	Chromosome Position	Average Heterozygosity
rs2358996	chr1	113973095	0.485767
rs246718	chr5	141142844	0.499477
rs2548467	chr5	65905662	0.49917
rs2548468	chr5	65905996	0.496714
rs2561243	chr5	65907119	0.498964
rs2615519	chr6	113792854	0.493477
rs284661	chr19	41426215	0.491565
rs3745290	chr19	41384098	0.499958
rs4241782	chr4	183978608	0.498945
rs4322101	chr9	126666398	0.497803
rs4387937	chr3	147649264	0.499918
rs441460	chr6	25548060	0.49681
rs4586881	chr4	161922650	0.499846
rs4588273	chr21	14266424	0.499653
rs4722897	chr7	29292617	0.498964
rs4751376	chr10	131393565	0.499017
rs4810570	chr20	46891166	0.499527
rs643502	chr5	141205942	0.499663
rs643864	chr5	141196668	0.499673
rs6442180	chr3	10546328	0.499563
rs6455790	chr6	161997415	0.499809
rs652539	chr18	44883475	0.499897
rs6657751	chr1	102003830	0.498652
rs6684496	chr1	15026324	0.499539
rs6691155	chr1	3994007	0.499853
rs6693853	chr1	15026462	0.499587
rs6782792	chr3	184456742	0.499539
rs6811404	chr4	64683631	0.499928
rs685553	chr18	44888700	0.499902
rs6927989	chr6	165696936	0.495671
rs6938076	chr6	41205390	0.499121
rs6994603	chr8	124418421	0.498852
rs7070511	chr10	131393337	0.49931
rs7102303	chr11	95708643	0.499846
rs738940	chr16	19343810	0.499265
rs754666	chr1	230114435	0.493523
rs7652997	chr3	178658471	0.498794
rs7774941	chr6	1443156	0.499954

Variant List – Chimeric ID Panel (PAGE 3 OF 3)

SNP RS#	Chromosome	Chromosome Position	Average Heterozygosity
rs7779384	chr7	83219014	0.499928
rs7842255	chr8	61925502	0.498964
rs7848571	chr9	125378927	0.497776
rs7903919	chr10	66011188	0.499809
rs7967526	chr12	93471056	0.499502
rs8037351	chr15	94818558	0.497776
rs838640	chr3	143481847	0.499587
rs924758	chr4	111498403	0.499464
rs9291849	chr5	65907845	0.49931
rs9372042	chr6	157091621	0.499999
rs945246	chr9	120111878	0.489837
rs9458428	chr6	162003547	0.499653
rs9552210	chr13	20493802	0.499053
rs962197	chr19	55088658	0.491525
rs9729034	chr1	84824065	0.495818
rs9862216	chr3	147650629	0.499879
rs9901757	chr17	83005455	0.499996
rs9927268	chr16	19345047	0.499265

For Research Use Only. Not for use in diagnostic procedures.

3

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: agenabio.com

US +1.877.4.GENOME
EU +49.40.899676.0
AP +61.7.3088.1600
CN +86.21.6427.0566

