




EZER™ Human Metapneumovirus Antigen Rapid Test



INTENDED USE

The EZER™ Human Metapneumovirus Antigen Rapid Test is intended to be used to detect Human Metapneumovirus antigen in nasopharyngeal swab, or pharyngeal swab to support the diagnosis of Human Metapneumovirus infection.



Hangzhou Genesis Biodetection and Biocontrol Co., Ltd.

+86-571 878 18163

www.genesis-ivd.com

No.139, St.10th (East), Hangzhou Economic & Technological Development Zone, China

PERFORMANCE CHARACTERISTICS

Human Metapneumovirus (hMPV)	Relative Sensitivity	Relative Specificity	Accuracy
NP Swabbing	92.0%	100%	96.3%
Pharyngeal Swabbing	91.5%	100%	94.2%

AVAILABILITY

Product	EZER™ Human Metapneumovirus Antigen Rapid Test
Cat. No.	P213115
Test Device	20
Sterilized Swab	20
Extraction Tube	20
Nozzle	20
Tube Stand	1
Package Insert	1

TEST PROCEDURE

Sample Preparation




Nasopharyngeal swab

OR



Pharyngeal swab

Sample Preparation



RESULTS

